

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
HLA-A allele^{1,2,3}																																
*30:01:01-30:01:07, 30:38-30:39, 30:43, 30:48-30:49, 30:53-30:54, 30:58, 30:62-30:63, 30:65, 30:72, 30:75, 30:78N ⁴	1	2							9	11	13						17															
*30:02:01-30:02:04, 30:02:06-30:02:12, 30:50, 30:57, 30:64, 30:66-30:69	1		3						9	11	13																					
*30:02:05			3						9	11	13																					
*30:03	1		3	4					9	11	13																					
*30:04:01-30:04:02					5				9	11	13																					
*30:06					5	6			9	11	13																					
*30:07	1		3			6			9	11	13																					
*30:08	1	2					7																									
*30:09	1								9	11	12	13																				
*30:10	1		3					8	9	11	13																					
*30:11:01-30:11:02	1	2		4							11							17														
*30:12	1		3						9	10	11																					
*30:13	1		3						9	11			13	14																		
*30:14L	1	2							9	11	13			15			17															
*30:15	1	2							9	11	13					16	17															
*30:16	1	2							9	11	13	14					17															
*30:17					5				9	11	13						17															
*30:18	1	2							9	10	11						17															
*30:19	1	2			5				9	11	13						17															
*30:20	1	2							9	11	13						17	18														
*30:22	1		3																19													
*30:23	1	2							9	11	13						17		20													
*30:24	1	2							9	11	13						17			21												
*30:25	1		3						9	11	13													24								
*30:26		2							9	11	13						17					23										
*30:27N	1		3						9	11	13											22										
*30:28	1						7		9	11	13																					
*30:29					5				9	11	13		15																			
*30:30		2							9	11	13						17			21												
*30:31	1	2							9		13						17		19													
*30:32	1		3						9	11	13							18														
*30:33	1		3						9	11	13					16																
*30:34	1		3						9	11	13									20												
*30:35	1	2							9	11	12	13					17															
*30:36	1	2							9	11	13						17				22											
*30:37	1	2							9	11	13						17													30		
*30:40	1	2							9	11	13						17								25							
*30:41	1	2							9	11							17															
*30:42	1	2							9	11	13						17									26						
*30:44	1	2							9	11	13	14																				
*30:45	1		3						9	11	13																27					
*30:46					5				9	11	13	14																				
*30:47	1								9	11	13																					
*30:51	1		3						9	11	13																		28			
*30:52	1								9	11	13						17															
*30:55	1	2								10																						
*30:56	1	2							9	11	13						17										27					
*30:59N, 30:74	1	2							9	11	13						17												29			
*30:60	1	2							9		13						17															
*30:61	1		3						9	11	13																		29			
*30:70N	1		3						9	11	13								19													
*30:71	1	2		4						11								18														
*30:73N	1	2							9	11	13						17				22								29			
*30:76N	1		3						9	11	13																				31	
*30:77					5				9	11	13								19													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*01:02, 01:20				4							w																					
*01:07, 26:19, 29:14, 31:01:02-31:01:06, 31:01:08-31:02, 31:05-31:07, 31:09-31:43, 31:45-31:50, 31:52-																																

