

## HLA-A31 2020 8K7

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43		
<b>HLA-A allele<sup>1,2,3</sup></b>																																													
<b>*31:01:02:01-31:01:02:02, 31:01:02:04-31:01:06, 31:01:08-31:01:29, 31:01:31-31:01:33, 31:01:35-31:01:38, 31:42, 31:45, 31:47, 31:49-31:53, 31:57-31:58, 31:61, 31:63-31:66, 31:70, 31:73-31:78, 31:80, 31:82-31:86, 31:88, 31:90, 31:92-31:94, 31:99-31:101, 31:103-31:104, 31:106-31:108, 31:110, 31:112-31:114, 31:116-31:118, 31:121-31:122, 31:124, 31:127-31:130, 31:132-31:141N, 31:143-31:145, 31:147-31:163, 31:165-31:167</b>	1								9			12																																42	
*31:01:02:03N, 31:35 <sup>4</sup>	1								9			12									21																						42		
*31:01:07, 31:01:34	1								9			12																																42	
*31:01:30, 31:69, 31:120	1											12																																42	
*31:02:01-31:02:02, 31:91	1	2							9			12																																42	
*31:03	1		3	4	5							12																																42	
*31:04:01:01-31:04:01:02	1		3	4								12																																42	
*31:05	1					6			9			12																																42	
*31:06	1			4								12																																42	
*31:07	1	2					7		9																																			42	
*31:08	1	2					7		9																																			42	
*31:09	1							8	9			12																																	42
*31:10	1						7		9																																			42	
*31:11	1								9	10		12																																42	
*31:12	1								9		11	12																																42	
*31:13	1								9			12	13																															42	
*31:14N	1								9			12		14																														42	
*31:15	1								9			12			15																													42	
*31:16, 31:46 <sup>4</sup>	1								9			12				16																												42	
*31:17	1								9			12					17																											42	
*31:18	1								9			12						18																										42	
*31:19	1								9			12							19																									42	
*31:20	1								9			12								20																								42	
*31:21	1								9												21																							42	
*31:22	1								9			12										22																						42	
*31:23	1								9			12											23																					42	

Negative Control











