

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
HLA-A allele ^{1,2,3}																																	
*33:01:01:01 -33:01:10, 33:49-33:50, 33:67, 33:89, 33:91, 33:121-33:122, 33:127	1	2	3																														
*33:03:01 -33:03:02, 33:03:04-33:03:32, 33:35, 33:37, 33:40-33:41, 33:43, 33:45- 33:48, 33:55-33:56, 33:59-33:63, 33:66, 33:71-33:72, 33:75-33:76, 33:78-33:79, 33:81, 33:87-33:88, 33:93-33:95, 33:97- 33:105, 33:107-33:108, 33:110, 33:112- 33:118, 33:120, 33:124, 33:126, 33:128, 33:130	1			4																													
*33:03:03Q, 33:86 ⁵	1			4																												30	
*33:04	1	2			5																												
*33:05	1	2	3			6																											
*33:06	1			4			7																										
*33:07	1	2	3					8																									
*33:08, 31:99 ⁵				4					9	10																							
*33:09, 02:309, 26:22, 66:09 ⁴				4					9																								
*33:10	1			4							11																						
*33:11, 33:80N ⁵	1			4								12																					
*33:12	1			4									13																				
*33:13	1			4						10																							
*33:14	1			4					9																								
*33:15	1			4										14																			
*33:16, 33:64 ⁵	1	2	3												15																		
*33:17	1			4												16																	
*33:18:01-33:18:02	1																17																
*33:19	1													14																			
*33:20	1			4		6																											
*33:21				4			7									16																	
*33:22	1			4					9	10																							
*33:23	1			4											15																		
*33:24	1			4				8																									

rol

*26:60N, 31:66, 31:89, 31:115				4								12																				
*26:68				4	5																											
*26:99, 31:29, 31:48				4											14																	
*26:107N, 32:98				4																								31				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*29:03											11						17															
*29:07, 29:49														14			17															
*29:09, 29:59																	17							24								
*29:12, 29:92						6			9								17															
*29:19								8									17															
*29:22									9	10							17															
*29:48								8	9	10							17															
*29:61																	17														31	
*29:66, 31:03										10							17															
*31:02, 31:07-31:08, 31:91, 31:109				4						10																						
*31:05, 32:13				4							11																					
*31:24, 32:33:01				4																				24								
*31:54				4				8																								
*31:94				4		6																										
*31:110, 32:10, 32:63, 74:28				4															19													
*32:51				4																											30	
*32:53				4																20												
*68:29	1											12																				
<i>B*51:136</i>																								23								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control