

Well No.	1	2	3	4	5	6	7	8	9	10	11
<b>HLA-A allele</b> <sup>1,2,3</sup>											
*34:01:01-34:01:02	1	2									
*34:02:01-34:02:03	1		3								
*34:03	1			4							
*34:04	1		3		5						
*34:05	1	2				6					
*34:06	1			4			7				
*34:07			3				7				
*34:08	1		3					8			
*34:09			3					8			
*34:10N	1		3							10	
*34:11	1	2							9		
*01:13, 01:17, 03:63, 03:88, 11:01:01-11:11, 11:13-11:16, 11:20-11:27, 11:29-11:39, 11:41-11:52Q, 11:54-11:95, 11:97, 11:99N-11:105, 11:107-11:120, 11:122-11:158, 11:160-11:176, 25:02, 26:13, 26:19, 26:33, 66:01:01-66:01:02, 66:04, 66:06-66:11, 66:13-66:14, 66:17-66:19, 69:02, C*07:335	1										
*01:51, 02:55, 02:135, 02:309, 02:454, 03:24, 25:03, 25:13, 26:20, 26:30, 26:65, 32:15, 66:02-66:03, 66:12, 66:16, 68:71								8			
*02:91, 02:322, 03:94, C*07:81, C*07:243						6					
*03:01:19				4				8			
*11:96							7				
*25:09, 26:14, 26:28, 26:73				4							
*26:18, 31:03				4			7				
*26:48, 26:69		2									
*31:01:07, B*15:82, B*15:260, B*40:186, C*03:186					5						
*31:04				4			7	8			
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>

Negative Control