

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
HLA-A allele ^{1,2,3}												
*34:01:01-34:01:02, 34:12, 34:14, 34:16	1	2										11
*34:02:01-34:02:04, 34:13, 34:15	1		3									11
*34:03	1			4								11
*34:04	1		3		5							11
*34:05	1	2				6						11
*34:06	1			4			7					11
*34:07			3				7					11
*34:08	1		3					8				11
*34:09			3					8				11
*34:10N	1		3							10		
*34:11	1	2							9			
*34:17	1	2		4								11
*01:13, 01:17, 01:176, 03:63, 03:88, 11:01:01:01-11:01:67, 11:01:69-11:11, 11:13-11:16, 11:20-11:27, 11:29-11:39, 11:41-11:52Q, 11:54-11:95, 11:97, 11:99N-11:105, 11:107-11:120, 11:122-11:158, 11:160-11:177, 11:179-11:190, 11:192-11:246, 11:249, 11:251N, 25:02, 26:13, 26:19, 26:33, 29:66, 66:01:01:01-66:01:03, 66:04, 66:06-66:11, 66:13-66:14, 66:17-66:20, 66:23-66:24, 66:27N, 69:02, C*07:335	1											
*01:51, 02:55, 02:135, 02:644, 03:24, 25:03, 25:13, 25:30, 26:20, 26:30, 26:65, 32:15, 66:02-66:03, 66:12, 66:16, 66:21, 66:25-66:26Q, 68:71								8				
*02:91, 02:322, 03:94, C*07:81, C*07:243, C*16:76						6						
*02:309, 02:454								8				11
*03:01:19				4				8				11
*03:248			3									
*11:96							7					
*11:191	1		3									
*25:06, 26:09, 26:91												11
*25:09, 26:14, 26:28, 26:73				4								
*26:18				4			7					
*26:48, 26:69		2										
*26:112				4			w					
*31:01:07, B*15:82, B*15:260, B*40:186:01, C*03:186:01, C*12:57:02					5							
*31:03				4			7					11
*31:04, 74:01:03				4			7	8				11
*66:22	1			4								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Negative Control