

INTERPRETATION TABLE						
HLA-A36 SSP subtyping						
Amplification patterns of the A*3601 to A*3604 alleles						
	Tube					
	1	2	3	4	5	6
Length of spec.	75	225	235	230	235	85
PCR product						
Length of int.	800	1070	800	1070	1070	1070
pos. control¹						
5'-primer²	527	363	363	355	363	527
	5'-TgC ^{3'}	5'-ATA ^{3'}	5'-ATA ^{3'}	5'-CCC ^{3'}	5'-ATA ^{3'}	5'-TgC ^{3'}
3'-primer³	559	545	555	545	559	570
	5'-CgT ^{3'}	5'-AgA ^{3'}	5'-CCA ^{3'}	5'-AgA ^{3'}	5'-CCg ^{3'}	5'-CAC ^{3'}
Tube No.	1	2	3	4	5	6
HLA-A allele						
*3601	+	+				+
*3602	+		+			+
*3603	+			+		+
*3604		+			+	+
*01010101-0102, 0104N, 0107-0111N, 0113, 0116N-0118N, 0121-0130, 0318		+				
*0114	W	+				
*03010101-0304, 0307-0309, 0311N-0317, 0319-0338, 3204			+			
*110101-110202, 1105-1107, 1109-1113, 1116-1119, 1121N-1124, 1129-1131					+	W
*1103, 1108, 1114, 1120, 1125, 1132					+	
*1104, 111501-111502, 1126						W
*2422				+		
HLA-A allele						
Tube No.	1	2	3	4	5	6