

Well No.	1	2	3	4	5	6	7	8
HLA-A allele^{1,2}								
*36:01, 03:187 ³	1	2				6		
*36:02	1		3			6		
*36:03	1			4		6		
*36:04, 01:72 ³		2			5	6		
*36:05	1	2				6	7	
*01:01:01:01-01:02, 01:04N, 01:07-01:11N, 01:13, 01:16N-01:18N, 01:21-01:33, 01:35-01:52N, 01:54-01:62, 01:64, 01:67-01:01:71, 01:74-01:88, 01:90-01:97, 01:99-01:126, 01:128-01:129, 01:131-01:135, 01:137-01:170, 03:18, 03:135, 11:94		2						
*01:14	w	2						
*01:73, 24:22, 24:160, B*15:90, B*38:45				4				
*01:89		2		4				
*03:01:01:01-03:01:29, 03:01:31-03:01:34, 03:01:36-03:01:47, 03:01:51-03:04:03, 03:07-03:09, 03:11N-03:17:02, 03:19-03:39, 03:41, 03:43-03:74, 03:76-03:94, 03:96-03:97, 03:99-03:104, 03:106-03:121, 03:123-01-03:134, 03:136-03:166, 03:168N-03:176, 03:178N-03:186, 03:188-03:203, 03:205-03:207, 03:209, 11:130, 30:89, 32:04			3					
*11:01:01:01-11:01:20, 11:01:22-11:01:43, 11:01:45-11:01:61, 11:01:63-11:03, 11:05-11:14, 11:16-11:25, 11:29-11:33:02, 11:36-11:37, 11:40-11:49, 11:51-11:52Q, 11:54-11:89, 11:91-01-11:93, 11:95-11:100, 11:102-11:111, 11:113-11:117, 11:119-01-11:129, 11:131-11:138, 11:140-11:142, 11:144-11:154, 11:156, 11:158-11:169, 11:171-11:181, 11:183-11:208N, 11:210N					5			
*11:112, 11:211		2			5			
*11:155					5	6		
*31:62	1					6		
*68:41, B*57:65, C*04:31, C*06:137						6		
Well No.	1	2	3	4	5	6	7	8

Negative Control