

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-A allele^{1,2,3}																
*66:01:01:01-66:01:07, 66:20, 66:23-66:24, 66:29-66:32, 66:35-66:38, 66:41-66:42, 66:44-66:46	1	2	3			6										
*66:02:01:01-66:02:01:03, 66:21, 66:25, 66:34			3	4										14		
*66:03:01:01-66:03:01:03, 66:39N				4			7							14		
*66:04	1	2	3			6				11						15
*66:05	1					6	8									
*66:06, 66:33	1	2	3			6			9							
*66:07	1	2	3			6	8			11						15
*66:08, 66:17 ⁴	1	2	3			6					12					
*66:09	1	2	3			6			10							
*66:10	1	2	3			6				11						
*66:11	1	2	3											13		
*66:12	1		3												14	
*66:13, 66:19 ⁴	1	2	3			6							13			
*66:14	1	2	3			6			9	11						
*66:15	1		3					8								
*66:16			3	4										14	15	
*66:18	1	2	3			6										15
*66:22	1	2	3			6							13	15		
*66:26Q			3										13	14		
*66:27N	1	2	3		5	6										
*66:28N				4			7			10						14
*66:40		2	3	4		6										
*66:43			3	4					9				13	14	15	
*01:01:01:01-01:01:27, 01:01:29-01:01:55, 01:01:58-01:01:61, 01:01:63-01:01:106, 01:01:108-01:01:123, 01:01:125-01:04:01:02N, 01:06, 01:08-01:12, 01:14-01:16N, 01:18N-01:27N, 01:29-01:33, 01:35-01:70, 01:72-01:99, 01:101-01:104, 01:106-01:133, 01:135-01:143, 01:145-01:175, 01:177-01:189, 01:191-01:193, 01:195-01:203, 01:205-01:228Q, 01:230-01:249, 01:251-01:288, 01:290N-01:297, 01:300-01:301Q, 01:303-01:325, 01:327-01:333, 01:335-01:419, 02:644, 02:749, 03:41, 11:17, 11:40, 11:98, 25:01:09, 25:01:17, 25:06, 25:43, 25:50, 25:52, 25:59, 26:01:59, 26:04, 26:09, 26:75, 26:91, 26:204:01:01-26:204:01:02, 31:164, 32:62, 33:13, 36:01:01:01-36:12, 74:10, 80:01:01:01-80:01:02, 80:03-80:09N								8								
*01:01:56	1							8		11			13			
*01:13, 01:176, 03:63, 03:88, 11:01:01:01-11:01:50, 11:01:52-11:01:67, 11:01:69-11:01:72, 11:01:74-11:01:75, 11:01:77-11:05, 11:07:01-11:11, 11:13-11:16, 11:20-11:27, 11:29-11:39, 11:41-11:44, 11:46-11:52Q, 11:54:01-11:95, 11:97, 11:99N-11:105, 11:107-11:110, 11:112-11:120, 11:122-11:158, 11:160-11:177, 11:179-11:189, 11:191-11:214, 11:216-11:222, 11:224-11:249, 11:251N-11:252, 11:254-11:270, 11:272Q-11:286, 11:288-11:290, 11:292-11:300, 11:302:01-11:304, 11:306-11:307, 11:309-11:311, 11:313Q-11:324, 11:326-11:395, 11:397N-11:398, 11:400N-11:428, 30:125, 30:140		2				6										
*01:17, 01:302, 03:382, 11:01:51, 11:45, 11:111, 11:190, 11:223, 11:253, 11:271, 11:287N, 11:305, 11:396, 29:66, C*07:335		2														
*01:28, 01:299, 02:346, 02:952, 11:19, 11:121, 11:308, 74:30, 80:02						6	8									
*01:134					5		8									
*01:229, 11:312, 26:78, 33:184					5	6	8									
*02:01:09, 02:05:05, 02:06:07, 02:07:12, 02:50, 02:76:02-02:76:03, 02:122, 02:243:02-02:243:03, 02:528:02, 02:591:01, 03:09, 03:89:02, 03:91N, 03:108, 03:172, 03:198, 03:330N, 11:18, 24:28, 24:89, 24:379, 29:19, 29:48, 29:159, 32:112N, 33:24, 33:154N, 68:05, 68:15, 68:20, 68:255, 74:06, 74:21, 74:39					5											
*02:03:01:01-02:03:10, 02:25, 02:38, 02:117, 02:148, 02:171:01-02:171:02, 02:230, 02:253, 02:258, 02:264:01:01-02:264:01:02, 02:267, 02:280-02:281, 02:315, 02:345, 02:355, 02:370, 02:412, 02:431, 02:447, 02:463, 02:466, 02:480, 02:505, 02:529, 02:544, 02:568, 02:595, 02:612, 02:633-02:634, 02:641, 02:666, 02:684, 02:711, 02:774-02:775N, 02:809, 02:833N, 02:849-02:850, 02:854, 02:860, 02:914, 02:949, 02:971, 02:976-02:977, 02:979, 02:988, 02:997Q, 02:1000, 02:1002										10						
*02:03:11					5					10						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*02:11:01:01-02:11:12, 02:34-02:35:03, 02:56:01-02:56:02, 02:62, 02:69, 02:78, 02:103, 02:128, 02:297-02:298, 02:308, 02:457, 02:490N, 02:494, 02:503, 02:514N, 02:579-02:580, 02:584, 02:589, 02:701, 02:705, 02:714, 02:738, 02:751, 02:754, 02:757, 02:773N, 02:808, 02:815, 02:821, 02:893, 02:906-02:907, 02:923, 02:926, 02:932, 02:959, 02:964, 02:986, 02:991, 02:1031N, 03:242, 03:312, 23:13, 24:07:01:01-24:07:04, 24:24, 24:112, 24:131, 24:288, 24:290, 24:294Q, 24:339, 24:347:01-24:347:02, 24:387, 24:406, 24:453, 24:457, 24:477, 24:489, 24:510, 24:528, 24:541, 24:544, 24:566, 24:569, 24:574, 29:155, 68:01:01:01-68:02:23, 68:04:01-68:04:02, 68:06-68:14, 68:16-68:19, 68:21:01-68:30, 68:32-68:56:02, 68:58:01-68:89, 68:91-68:108, 68:110-68:179, 68:181N-68:205N, 68:207-68:226, 68:228-68:248, 68:250-68:254, 68:256-68:276Q, 68:278-68:293Q, 69:01:01:01-69:01:04, 69:03-69:09			3													
*02:16, 02:131:01:01-02:131:03, 02:487, 02:560, 02:693							7									
*02:19, 02:294, 02:453, 02:690, 32:54, 32:142, 33:151, 33:204, 68:206												12				
*02:135, 25:13, 25:82, 26:65	1														14	

Negative Control

*02:309, 02:454								9						14		
*02:427					6		8		10							
*02:557									10		12					
*02:741, 34:10N-34:12, 34:16, 69:02		2	3		6											
*02:940			3								12					
*03:01:19, 03:103:02, 31:04:01:01-31:04:02, 31:123, 74:01:03								9				13	14	15		
*11:01:68, 11:01:76, 11:96, 11:106, 11:301, 11:325					6											
*11:06, 11:215N		2			5	6										
*23:122, 24:109, 24:260							w									
*24:19			3			6										
*24:44						6	w									
*24:315, 24:374							w									
*24:535													13			
*25:01:01:01-25:01:08, 25:01:10-25:01:16, 25:01:18-25:01:21, 25:03-25:05, 25:07, 25:10-25:12N, 25:14, 25:16, 25:18, 25:20-25:26, 25:27:02-25:29, 25:31-25:35, 25:37-25:42N, 25:44-25:49N, 25:51, 25:53-25:58, 25:60-25:63, 25:66-25:70, 25:72-25:74, 25:76-25:80, 26:01:01:01-26:01:21, 26:01:23-26:01:27, 26:01:29-26:01:58, 26:01:60-26:02:02, 26:05, 26:07:01-26:08:03, 26:10-26:12:01:02, 26:15-26:17, 26:20:01-26:20:02, 26:23-26:27, 26:31-26:32, 26:34-26:42, 26:45-26:46, 26:48, 26:50-26:63, 26:66-26:72, 26:74, 26:76-26:77, 26:79-26:82, 26:84-26:90, 26:93-26:110, 26:113-26:129, 26:131, 26:133-26:145N, 26:147, 26:149-26:150, 26:152-26:176, 26:178-26:192, 26:196-26:203, 26:206:01N-26:216, 26:218-26:223, 43:01-43:03	1						8									
*25:02, 25:75, 26:33	1	2														
*25:08, 26:47	1						8							15		
*25:09, 26:14, 26:18, 26:28, 26:73, 26:146	1						8					13		15		
*25:15, 25:17, 25:64, 25:71, 25:81, 26:01:22, 26:01:28, 26:64, 26:83, 26:132, 26:148, 26:151, 26:194-26:195	1															
*25:19:01-25:19:02, 25:30, 26:43:01-26:43:02	1						8	9								
*25:27:01, 26:130	1						8	9		11						
*25:36, 26:205							8			11						
*25:65, 26:29, 26:49	1						8			11						
*26:03:01:01-26:03:01:02, 26:06, 26:21, 26:92, 26:111, 26:177, 26:193	1				5	6	8									
*26:13, 26:19	1	2	3													
*26:22	1						8		10							
*26:30	1				5								14			
*26:112:01-26:112:02	1						8	9				13		15		
*26:217	1						8					13				
*31:03							8	9				13		15		
*34:01:01:01-34:01:07, 34:05, 34:18-34:19, 34:23, 34:27		2	3		6						?					
*34:02:01:01-34:02:07, 34:04, 34:08, 34:13, 34:15, 34:20-34:21, 34:24-34:26N, 34:28-34:30N		2	3		6			9								
*34:03, 34:06		2	3		6			9				13		15		
*34:07			3		6			9								
*34:09			3					9					14			
*34:14		2			6					?						
*34:17		2	3		6							13		15		
*34:22		2	3		6			9						15		
*68:227		2	3													
*68:294			3										14			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control