

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-A allele																								
*680101-680107	1	2																				22	23	
*68020101-68020103, 680202	1		3														17							
*680301- 680302		2		4																		22	23	
*6804		2			5																	22	23	
*6805		2		4		6																22	23	
*6806	1					7																		
*6807	1							8														22	23	
*6808	1	2							9														23	
*6809	1	2								10													23	
*6810	1										11											22	23	
*6811N	1	2										12										22	23	
*6812	1	2											13		?							22	23	
*6813	1													14								22	23	
*6814	1										11				15							22	23	
*6815			3			6										16	17							
*6816	1	2														16						22	23	
*6817	1	2									11											22	23	
*6818N	1		3						9								17							
*6819	1	2												15								22	23	
*6820		2		4		6											17					22	23	
*6821	1	2																18				22	23	
*6822	1	2																19				22	23	
*6823	1	2																	20			22	23	
*6824	1	2																			21	22	23	
*6825	1	2			5					10												22	23	
*6826	1	2								10								18						
*6827	1																17							
*6828	1		3														17				21			
*6829	1	2																				22	23	
*6830	1			4																		22	23	
*6831			3												16	17								
*6832	1	2				7																22	23	
*6833	1	2																				22	23	
*6834	1		3														17							
*6835	1	2																				22	23	
*6836		2																				22	23	
*6837	1	2																				22	23	
*6838	1	2																				22	23	
*6839	1	2	3																			22	23	
*6840	1		3												16	17								
*6841	1	2												14								22	23	
*6842	1	2																				22	23	
*01010101-0104N, 0107-0111N, 0114-0118N, 0121-0124, 0126-0129																		18				w		
*0106, 0125, 2903																		18						
*0112, 0119										10								18						
*0113, 0318	1																	18				w		
*0120																		18		20		w		
*0130, 0222, 0250, 9204, 92362302, 2406, 2611N, 2614, 2618, 2628, 3210																	17							
*0202, 0205, 0208, 0247, 0263, 9202, 9215																	w							

020301-020302										17					
*020605	4			6											
*0213, 0226, 0238, 0299															
*0234-023503, 0256, 0262, 0278, 9203, 0319, 1122, 1123, 1131, 2419, 3018, 3406	1														
*0249										17			21		
*0252				7											
*0255													24		
*027402															
*9214										20			23		
*9217										17					
*9221, 2418							15								
*9222				6			11				16	17			
*9229							11								
*9235, 250101-2503, 2505, 2506, 260101-260104, 2602, 2603, 2605-260702, 2609, 2610, 2612, 2613, 2615-2617, 2619-2621, 2623-2626, 2629-2635, 4301, 6603, 6605												17			
*03010101-03010103, 030103-0303N, 0305-0307, 0309, 0311N-0316, 0320-0323, 0325-0331, 0333-0335, 0337-0339	1												23		
*030102	1						9							23	
*0304	1									18			23		
*0308, 0332, 0336N, 29010101-290201, 290203, 2904, 2906-2917, 3103, 3104, 3106, 3204														23	
*0310, 1104	1						10						23		
*0317	1								15					23	
*0324	1	2											23		
*110101-110109, 1105-1107, 1111-1113, 111501-111502, 1121N, 1126, 1129, 1130, 1132	1							10							
*110201-110203, 1114, 1116-1119							10								
*1103, 1120	1						10								
*1108	1												22		
*1109	1						w								
*1110	1	2					10								
*1124	1						9								
*1125	1														
*1127	1						10					18			
*2417, 2441							10						18		
*2482															
2504, 260105-260106, 2608, 2627															
*2604												17			
*2622													17		
*2905							10							23	

[illegible]

[illegible]

[illegible]