

Well No.	1	2	3	4	5	6
HLA-A allele^{1,2}						
*69:01	1	2	3			
*69:02	1			4		
*69:03	1	2	3		5	
*01:20, 01:66, 02:346, 02:427, 26:22, 66:09				4		
*02:01:01:01-02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:70, 02:01:72-02:01:81, 02:01:83-02:01:86, 02:01:88-02:01:109, 02:01:111-02:01:115, 02:03:01-02:03:08, 02:06:01-02:06:18, 02:09, 02:11:01-02:11:07, 02:16, 02:20:01-02:22:02, 02:24:01-02:26, 02:28-02:32N, 02:36, 02:40:01-02:43N, 02:45-02:46, 02:48-02:49, 02:51-02:52, 02:59-02:61, 02:64, 02:66-02:72, 02:74:01-02:77, 02:79:01-02:83N, 02:85-02:93, 02:95-02:97:02, 02:99, 02:101:01-02:101:02, 02:104-02:107, 02:109, 02:111, 02:113:01N-02:114, 02:116-02:121, 02:123-02:129, 02:131-02:134, 02:136-02:141, 02:143-02:147, 02:149-02:151, 02:153:01-02:153:02, 02:158-02:171:02, 02:173-02:178, 02:181-02:185, 02:187-02:190, 02:192-02:193, 02:195-02:198, 02:200-02:206, 02:208, 02:210-02:211:02, 02:214-02:218, 02:220-02:223N, 02:225N, 02:227N-02:228, 02:230, 02:233-02:241, 02:243:01-02:243:02, 02:246-02:254, 02:256-02:260, 02:262-02:264, 02:266-02:267, 02:269-02:270, 02:272-02:279, 02:281, 02:283-02:285, 02:287-02:299, 02:301N, 02:304-02:318, 02:321N-02:323, 02:325-02:333, 02:336, 02:339-02:343, 02:347-02:349, 02:351-02:358, 02:360-02:368, 02:370-02:372, 02:374-02:375, 02:377-02:397, 02:399-02:402, 02:404-02:407, 02:409-02:412, 02:414-02:416, 02:418-02:420, 02:422-02:425, 02:428, 02:430-02:431, 02:434-02:436, 02:438-02:448, 02:454-02:471, 02:475-02:476N, 02:479-02:483, 02:485-02:488, 02:490N-02:491, 02:493-02:494, 02:497-02:498, 02:500Q-02:501N, 02:503-02:506N, 02:508-02:510, 02:512, 02:514N-02:516N, 02:518-02:525N, 02:528-02:530, 02:533-02:539, 02:542, 33:22			3			
*02:12-02:13, 02:19, 02:27, 02:37-02:38, 02:44, 02:54, 02:142, 02:226N, 02:280, 02:417, 02:543, 24:14			3		5	
*02:34-02:35:03, 02:56:01-02:56:02, 02:62, 02:78, 03:95, 29:22, 68:105		2	3			
*02:55	1		3			
*02:103, 03:01:01:01-03:07, 03:09-03:28, 03:31, 03:33, 03:35, 03:37-03:40, 03:42-03:61, 03:64-03:87, 03:90-03:94, 03:96-03:106, 03:109-03:151, 03:153-03:171, 03:173-03:175, 03:177-03:197N, 03:199-03:209, 11:12, 24:24, 24:290, 29:01:01:01-29:18, 29:21, 29:23-29:29, 29:31-29:33, 29:35-29:47, 29:49-29:65, 29:67-29:73, 30:01:01-30:01:10, 30:08, 30:11:01-30:11:02, 30:14L-30:20, 30:23-30:26, 30:30, 30:35-30:43, 30:48-30:49, 30:52-30:56, 30:58-30:60, 30:62-30:63, 30:65, 30:71-30:75, 30:78N-30:79, 30:81-30:83, 30:86-30:89, 30:91-30:92, 31:89, 32:17, 34:09, 66:02, 66:12, 66:16, 68:01:01:01-68:01:17, 68:01:19-68:02:08, 68:06-68:14, 68:16-68:19, 68:21:01-68:30, 68:32-68:39, 68:41-68:71, 68:73-68:84, 68:86-68:89, 68:91-68:92, 68:94N-68:96, 68:98-68:104:02, 68:106-68:108, 68:110-68:124, 74:13		2				
*24:93					5	
Well No.	1	2	3	4	5	6

Negative Control