

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-A allele ^{1,2,3}																
*74:01:01-74:01:03	1	2														
*74:02:01:01-74:02:01:02	1		3													
*74:03	1	2							9							
*74:04	1	?	?	4												
*74:05	1	?	?		5											
*74:06	1	?	?			6										
*74:07	1		3				7									
*74:08	1	?	?					8								
*74:09	1	?	?							10						
*74:10		?	?								11					
*74:11	1	2										12				
*74:12N	1	?	?										13			
*74:13	1	2												14		
*74:14N	1	?	?						9		11					
*74:15	1	?	?													15
*74:16:01-74:16:02	1	?	?					8								
*74:17	1	?	?							10				13		
*74:18-74:20, 74:22-74:27	1	?	?													
*74:21	1		3			6										
*01:01:68, 01:69:02, 03:01:03				w												
*01:83:02, 02:01:09, 02:05:05, 02:06:07, 02:11:02, 02:35:03, 02:76:02, 02:243:02, 03:23:01, 03:198, 11:70:02, 11:121, 24:02:76, 24:310:02, 26:07:02, 33:01:02, 33:08-33:09, 68:30				4												
*03:01:39	1		3												14	
*03:27, 33:81										10			13			
*29:01:01:01-29:01:06, 29:01:08-29:03, 29:11, 29:40, 29:46, 29:56- 29:58, 29:70, 32:01:01-32:01:17, 32:01:19-32:01:21, 32:01:23-32:02, 32:05-32:06, 32:17, 32:19N, 32:61, 32:65, 32:68-32:69			3													
*29:04-29:10, 29:12-29:15, 29:17-29:18, 29:20-29:38, 29:42-29:45, 29:47, 29:49-29:55, 29:59-29:65, 29:67-29:69, 29:71-29:83, 32:03, 32:07-32:16, 32:18, 32:20-32:60, 32:63-32:64, 32:66-32:67, 32:70- 32:74, 32:76-32:84			?													
*29:16, 32:75			?							10			13			
*29:19	1		?			6										
*29:39	1		3				7									15
*29:48	1		?	4		6										
*29:66, 32:62			?								11					
*68:01:28	1		?												14	
C*04:01:51, C*06:02:15											11					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control