

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-A allele^{1,2,3}																
*74:01:01-74:01:08, 74:25, 74:28-74:29, 74:31	1	2														
*74:02:01:01-74:02:01:02	1		3													
*74:03	1	2							9							
*74:04	1	?	?	4												
*74:05	1	?	?		5											
*74:06, 74:21	1		3			6										
*74:07	1		3				7									
*74:08, 74:16:01 ⁴	1	?	?					8								
*74:09	1	2								10						
*74:10		?	?								11					
*74:11	1	2										12				
*74:12N	1	?	?										13			
*74:13	1	2														14
*74:14N	1	?	?						9		11					
*74:15	1	?	?													15
*74:16:02	1	2						8								
*74:17	1	?	?							10				13		
*74:18-74:20, 74:22-74:24, 74:26-74:27	1	?	?													
*74:30	1		3								11					
*74:32N	1	2			5				9							
*01:01:68, 01:69:02, 03:01:03				w												
*01:83:02, 02:01:09, 02:05:05, 02:06:07, 02:11:02, 02:35:03, 02:76:02, 02:243:02-02:243:03, 02:528:02, 03:23:01, 03:198, 11:70:02, 11:121, 23:01:21, 24:02:76, 24:310:02, 26:07:02, 33:01:02, 33:08-33:09, 68:01:41, 68:30				4												
*03:01:39, 68:01:36, 68:130:02	1		3													14
*03:27, 33:81, 33:122										10			13			
*29:01:01:01-29:01:06, 29:01:08-29:02:22, 29:02:24, 29:02:26-29:04, 29:09-29:10:01, 29:11-29:12, 29:23, 29:30, 29:32, 29:37, 29:40, 29:46, 29:56-29:58, 29:70-29:71, 29:84, 29:87, 29:90, 29:92-29:95, 29:100, 29:109-29:110, 29:114-29:121, 32:01:01:01-32:01:17, 32:01:19-32:01:21, 32:01:23-32:03, 32:05-32:08, 32:11Q, 32:15, 32:17, 32:19N, 32:24, 32:26:02-32:27N, 32:30:02-32:31, 32:44, 32:54-32:55:01, 32:60-32:61, 32:65, 32:68-32:69, 32:86, 32:92N-32:94, 32:98, 32:101Q-32:105, 32:107, 32:109-32:112N, 32:114, 32:116-32:118			3													
*29:02:23			3	4												
*29:05-29:08N, 29:10:02, 29:13-29:15, 29:17-29:18, 29:20-29:22, 29:24-29:29, 29:31, 29:33-29:36, 29:38, 29:42-29:45, 29:47, 29:49-29:55, 29:59-29:65, 29:67-29:69, 29:72-29:83, 29:85-29:86, 29:88-29:89, 29:91, 29:96-29:99, 29:101-29:108, 29:111-29:113, 32:09-32:10, 32:12-32:14, 32:16, 32:18, 32:20-32:23, 32:25-32:26:01, 32:28-32:30:01, 32:32-32:43, 32:45N-32:53, 32:55:02-32:59, 32:63-32:64, 32:66-32:67, 32:70-32:74, 32:76-32:85, 32:87-32:91, 32:95-32:97, 32:99-32:100, 32:108, 32:113			?													
*29:16, 32:75			?							10			13			
*29:19	1	?				6										
*29:39	1		3				7									15
*29:48	1		?	4		6										
*29:66			?								11					
*31:138								8								
*32:62			3								11					
*32:106		2														
*33:80N					5											
*68:01:28	1		?													14
C*04:01:51, C*06:02:15, C*07:02:76											11					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control