

| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 |
|---|----------|----------|----------|----------|----------|----------|----------|----------|----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| HLA-A allele^{1,2,3} | | | | | | | | | | | | | | | | |
| *74:01:01:01-74:01:08, 74:25, 74:28-74:29, 74:31 | 1 | 2 | | | | | | | | | | | | | | |
| *74:02:01:01-74:02:01:02, 74:34 | 1 | | 3 | | | | | | | | | | | | | |
| *74:03 | 1 | 2 | | | | | | | 9 | | | | | | | |
| *74:04 | 1 | ? | ? | 4 | | | | | | | | | | | | |
| *74:05 | 1 | ? | ? | | 5 | | | | | | | | | | | |
| *74:06, 74:21 | 1 | | 3 | | | 6 | | | | | | | | | | |
| *74:07 | 1 | | 3 | | | | 7 | | | | | | | | | |
| *74:08, 74:16:01 ⁴ | 1 | ? | ? | | | | | 8 | | | | | | | | |
| *74:09 | 1 | 2 | | | | | | | | 10 | | | | | | |
| *74:10 | | ? | ? | | | | | | | | 11 | | | | | |
| *74:11 | 1 | 2 | | | | | | | | | | 12 | | | | |
| *74:12N | 1 | ? | ? | | | | | | | | | | | 13 | | |
| *74:13 | 1 | 2 | | | | | | | | | | | | | | 14 |
| *74:14N | 1 | ? | ? | | | | | | 9 | | 11 | | | | | |
| *74:15 | 1 | ? | ? | | | | | | | | | | | | | 15 |
| *74:16:02 | 1 | 2 | | | | | | 8 | | | | | | | | |
| *74:17 | 1 | ? | ? | | | | | | | 10 | | | | 13 | | |
| *74:18-74:20, 74:22-74:24, 74:26-74:27 | 1 | ? | ? | | | | | | | | | | | | | |
| *74:30 | 1 | | 3 | | | | | | | | 11 | | | | | |
| *74:32N | 1 | 2 | | | 5 | | | | 9 | | | | | | | |
| *74:33 | 1 | 2 | | | | | 7 | | 9 | | | | | | | |
| *74:35 | 1 | | 3 | | 5 | | | | | | | | | | | |
| *01:01:68, 01:69:02, 03:01:03 | | | | w | | | | | | | | | | | | |
| *01:83:02, 02:01:09, 02:05:05, 02:06:07, 02:07:12, 02:11:02, 02:35:03, 02:76:02, 02:243:02-02:243:03, 02:528:02, 03:23:01, 03:198, 11:70:02, 11:121, 23:01:21, 24:02:76, 24:310:02, 26:07:02, 33:01:02, 33:08-33:09, 68:01:41, 68:30 | | | | 4 | | | | | | | | | | | | |
| *03:01:39, 68:01:36, 68:130:02 | 1 | | 3 | | | | | | | | | | | | | 14 |
| *03:27, 33:81, 33:122 | | | | | | | | | | 10 | | | 13 | | | |
| *29:01:01:01-29:01:06, 29:01:08-29:02:22, 29:02:24, 29:02:26-29:04, 29:09-29:10:01, 29:11-29:12, 29:23, 29:30, 29:32, 29:37, 29:40, 29:46, 29:56-29:58, 29:70-29:71, 29:84, 29:87, 29:90, 29:92-29:95, 29:100, 29:109-29:110, 29:114-29:132, 32:01:01:01-32:01:17, 32:01:19-32:01:21, 32:01:23-32:03, 32:05-32:08, 32:11Q, 32:15, 32:17, 32:19N, 32:24, 32:26:02-32:27N, 32:30:02-32:31, 32:44, 32:54-32:55:01, 32:60-32:61, 32:65, 32:68-32:69, 32:86, 32:92N-32:94, 32:98, 32:101Q-32:105, 32:107, 32:109-32:112N, 32:114, 32:116-32:126N | | | 3 | | | | | | | | | | | | | |
| *29:02:23 | | | 3 | 4 | | | | | | | | | | | | |
| *29:05-29:08N, 29:10:02, 29:13-29:15, 29:17-29:18, 29:20-29:22, 29:24-29:29, 29:31, 29:33-29:36, 29:38, 29:42-29:45, 29:47, 29:49-29:55, 29:59-29:65, 29:67-29:69, 29:72-29:83, 29:85-29:86, 29:88-29:89, 29:91, 29:96-29:99, 29:101-29:108, 29:111-29:113, 32:09-32:10, 32:12-32:14, 32:16, 32:18, 32:20-32:23, 32:25-32:26:01, 32:28-32:30:01, 32:32-32:43, 32:45N-32:53, 32:55:02-32:59, 32:63-32:64, 32:66-32:67, 32:70-32:74, 32:76-32:85, 32:87-32:91, 32:95-32:97, 32:99-32:100, 32:108, 32:113 | | | ? | | | | | | | | | | | | | |
| *29:16, 32:75 | | | ? | | | | | | | 10 | | | 13 | | | |
| *29:19 | 1 | | ? | | | 6 | | | | | | | | | | |
| *29:39 | 1 | | 3 | | | | 7 | | | | | | | | | 15 |
| *29:48 | 1 | | ? | 4 | | 6 | | | | | | | | | | |
| *29:66 | | | ? | | | | | | | | 11 | | | | | |
| *31:138 | | | | | | | | 8 | | | | | | | | |
| *32:62 | | | 3 | | | | | | | | 11 | | | | | |
| *32:106 | | | 2 | | | | | | | | | | | | | |
| *33:80N | | | | | 5 | | | | | | | | | | | |
| *68:01:28 | 1 | | ? | | | | | | | | | | | | | 14 |
| C*04:01:51, C*06:02:15, C*07:02:76 | | | | | | | | | | | 11 | | | | | |
| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 |

Negative Control