

INTERPRETATION TABLE			
HLA-A*80 SSP typing			
	Well		
	1	2	3
Length of spec.	135	155	165
PCR product			
Length of int.	800	1070	1070
pos. control¹			
5'-primer²	176	176	176
	5' -gCA 3'	5' -gCA 3'	5' -gCA 3'
3'-primer³	270	292	299
	5' -ACA 3'	5' -gTT 3'	5' -CCA 3'
Well No.	1	2	3
HLA-A allele			
*80:01	1	2	
*80:02	1		3
HLA-A allele			
Well No.	1	2	3