

## INTERPRETATION TABLE

### HLA-A\*80 SSP typing

	Well		
	1	2	3
Length of spec.	135	155	165
PCR product			
Length of int.	800	1070	1070
pos. control <sup>1</sup>			
5'-primer <sup>2</sup>	176	176	176
	5' -gCA 3'	5' -gCA 3'	5' -gCA 3'
3'-primer <sup>3</sup>	270	292	299
	5' -ACA 3'	5' -gTT 3'	5' -CCA 3'
Well No.	1	2	3
HLA-A allele <sup>4</sup>			
*80:01:01:01-80:01:01:02	1	2	
*80:02, 30:56 <sup>5</sup>	1		3
HLA-A allele			
Well No.	1	2	3