

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
HLA-B allele^{1,2,3}																																
*18:01:01:01-18:01:18, 18:55, 18:59-18:60, 18:62-18:66, 18:70-18:71, 18:76, 18:80, 18:82-18:91	1	2																			21											
*18:02	1		3																		21											
*18:03		2		4																	21											
*18:04	1	2			5									14							21					26						
*18:05	1	2				6															21											
*18:06	1	2																		19	21											
*18:07:01-18:07:02	1	2					7														21											
*18:08	1	2							9												21											
*18:09		2								11											21											
*18:10	1							8																								
*18:11	1																															
*18:12	1	2										12									21											
*18:13	1	2											13								21											
*18:14	1													14							21											
*18:15	1	2													15				18		21											
*18:17N	1	2														16					21											
*18:18	1	2																	17		21											
*18:19	1	2											13		15	16					21											
*18:20	1	2																			20											
*18:21	1							8											w						23							
*18:22	1																				21	22										
*18:23N, 18:28	1	2							10												21											
*18:24	1	2			5																21											
*18:25	1	2												14							21					26						
*18:26	1	2																			21			24								
*18:27	1								9												21											
*18:29		2																			18	21										
*18:30	1	2											13		15				w		21											
*18:31	1	2													15				17		21											
*18:32	1	2																		18	21											
*18:33	1	2																			21	22										
*18:34, 18:52	1	2																			21									30		
*18:35	1							8																	25							
*18:36, 18:67	1	2								11											21											
*18:37, 18:46	1	2																			21							27				
*18:38, 18:45	1	2																			21								28			
*18:39	1																				19	21										
*18:40, 18:51	1	2																			21			23								
*18:41	1	2					7					12									21											
*18:42-18:43	1	2																			21					26						
*18:44:01-18:44:02	1												13								21											
*18:47	1	2																			21									29		
*18:48	1	2																			21				25							
*18:49	1	2																			21											31
*18:50	1	2						8													21									29		
*18:53	1	2																			19	21		23								
*18:54		2								11											21				25							
*18:56, 18:58	1									11											21											
*18:57	1	2								11					15				18		21											
*18:61	1	2																			21			23	24							
*18:68	1	2												14							21											
*18:69	1																				21	w										
*18:72		2			5																21											
*18:73	1	2										12									21								27			
*18:74N	1	2					7														21							27	28			
*18:75	1	2							9												21	22										
*18:77	1	2										12									21	22										
*18:78	1	2																			21							27	28			
*18:79	1	2					7														21											31
*18:81	1	2						8													21											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Well No.	1	2	3</																													

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
----------	---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----