

Well No.

HLA-B allele^{1,2,3}

1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6
---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---

*35:55
*35:56
*35:57
*35:58, 35:519
*35:59:01-35:59:02
*35:60
*35:61:01
*35:61:02, 35:556Q
*35:62
*35:63
*35:64:01, 35:331
*35:65Q
*35:66
*35:67
*35:68:01
*35:68:02
*35:69
*35:70
*35:71
*35:72
*35:74
*35:75, 35:205:03, 35:521
*35:76
*35:77
*35:78
*35:79
*35:80
*35:81
*35:82
*35:83
*35:84
*35:85
*35:86
*35:87
*35:88
*35:89
*35:90, 35:441
*35:91
*35:92
*35:93

*35:94
*35:95
*35:96
*35:97
*35:98
*35:99
*35:100
*35:101:01-35:101:02
*35:102
*35:103
*35:104
*35:105
*35:106
*35:107
*35:108:01-35:108:02
*35:109
*35:110, 35:558
*35:111
*35:112:01
*35:113
*35:114

Negative Control

Well No.	HLA-B allele ^{1,2,3}	W														+														Negative Control							
		1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6
*40:330																																					
*40:339																																					
*40:343																																					
*40:368																																					
*40:377																																					
*40:422																																					
Well No.																																					
HLA-B allele ^{1,2,3}	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	
*40:429																																					
*40:458, 44:09, 44:46:01-44:46:02, 44:75, 44:129, 44:131, 44:254, 44:48N	+																																				
*40:496																																					
*40:522																																					
*40:523																																					
*40:537																																					
*41:01:01-41:02:03, 41:02:05-41:03:02, 41:06-41:07, 41:09, 41:11, 41:13-41:17, 41:20, 41:22-41:23, 41:25, 41:28-41:29, 41:31-41:42, 41:44-41:45N, 41:53-41:59, 41:61, 41:64-41:81																																					
*41:04	+																																				
*41:10, 41:52																																					
*41:18	+																																				
*41:21, 41:63																																					
*41:24					+																																
*41:27, 41:47																																					
*41:30																																					
*41:43																																					
*41:46																																					
*41:48																																					
*41:49																																					
*41:50																																					
*41:51																																					
*41:62																																					
*42:01:01-42:02:02, 42:04-42:05:02, 42:07-42:08, 42:12, 42:14- 42:15, 42:18, 42:21-42:24, 42:27-42:31, 42:33-42:34																																					
*42:06																																					
*42:09					+																																
*42:10, 42:17																																					
*42:13																																					
*42:16:01																																					
*42:19																																					
*42:25																																					
*42:26																																					
*42:32																																					
Well No.																																					
HLA-B allele ^{1,2,3}	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	

The figure displays a grid of genotype data for 100 individuals (columns) across approximately 100 markers (rows). Shaded cells represent positive genotypes, while white cells represent negative genotypes. A vertical column on the far right is labeled "Negative Control".

Marker Labels:

- *57:25, *57:26, *57:31:02, *57:33, 57:40, 57:48, 57:77, 57:85, 57:110, 57:134
- *57:58, *57:63, *57:67:01, *57:67:02, *57:76, *57:80, *57:82, *57:83, *57:93, *57:107, *57:113, *57:118, *57:119, *57:140, *57:145, *57:147, *57:161
- *58:01:01-58:01:01:20, 58:01:04, 58:01:06-58:01:15, 58:01:17-58:01:19, 58:01:21, 58:01:23-58:01:35, 58:01:37-58:01:44, 58:04, 58:13-58:15, 58:21, 58:23-58:24, 58:29, 58:31N-58:33, 58:36, 58:39N-58:41, 58:45:01-58:45:02, 58:47-58:58, 58:62-58:63, 58:66-58:72N, 58:74-58:75, 58:79-58:81, 58:83, 58:85-58:87, 58:89-58:90, 58:92-58:95, 58:97-58:101, 58:103, 58:107, 58:109-58:113, 58:115, 58:117-58:120, 58:123, 58:125-58:126, 58:128N-58:131, 58:133N, 58:135-58:136, 58:139-58:143
- *58:01:02, *58:01:03, 58:01:16, 58:02:01:01-58:02:01:03, 58:16:02, 58:25-58:26, 58:38, 58:43, 58:60-58:61, 58:104, 58:127, 58:132
- *58:01:05, 58:82, 58:102, 58:116, 58:134
- *58:01:20, *58:01:22, 58:01:45, 58:105-58:106
- Well No.**
- HLA B allele^{1,2,3}**

Marker Legend:

- White box: Negative genotype
- Shaded box: Positive genotype
- Text label: Specific marker or well identifier
- W: Denotes a wild-type allele











