

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	
HLA-B allele^{1,2,3}																															
*38:01:01-38:01:02, 38:01:04-38:01:07, 38:01:09-38:01:12, 38:01:14, 38:38-38:40, 38:42, 38:51, 38:53, 38:56-38:60, 38:66, 38:69, 38:71, 38:73, 38:77-38:78, 38:80N, 38:83N	1	2	3	4									13	15														27			
*38:01:03, 38:01:13, 38:67	1	2	3	4									13	15																	
*38:01:08			2	3									13	15														27			
*38:01:15	1	2	3										13	15														27			
*38:02:01-38:02:02, 38:02:04-38:02:08, 38:43-38:44, 38:47-38:49, 38:62, 38:72, 38:75-38:76, 38:79, 38:82	1	2	3		5								13	15																	
*38:02:03			2	3	5								13	15																	
*38:03	1				5	6	7						13																		
*38:04			2	3	5	7							13	15																	
*38:05			2	3	4			8					13																27		
*38:06			2						9				13	15													26	27			
*38:07	1	2							9	10			13	15													26	27			
*38:08	1	2	3		5						11			15																	
*38:09	1	2	3	4								12	13	15																	
*38:10	1	2	3	4										14	15				19				23					27			
*38:11, 38:19 ⁴	1	2	3	4									13	15	16														27		
*38:12	1	2	3	4							10		13	15															27		
*38:13	1		3	4							11		13																27		
*38:14	1		3	4								12	13																27		
*38:15	1	2	3		5			8					13	15																	
*38:16	1			4							10		13	15															27		
*38:17	1	2	3										13	15	16									24							
*38:18	1	2	3		5							11		15																	
*38:20:01-38:20:02	1	2	3	4				8					13	15															27		
*38:21	1	2	3	4			7		9				13	15															27		
*38:22	1	2	3	4				8						15															27		
*38:23	1	2	3		5								13	15			17												27		
*38:24	1	2	3	4									13	15				18											27		
*38:25			2	3			7						13	15															27		
*38:26	1			4								12	13																27		
*38:27	1	2	3	4								11		13	15														27		
*38:28	1			4									13	15						20									27		
*38:29	1	2	3		5									15							22								27		
*38:30	1	2	3	4										15					19		21								27		
*38:31	1	2	3	4									13	15												25			27		
*38:32	1	2	3	4									13	14	15									23					27		
*38:33			2	3	4			8		10			13																27		
*38:34N	1	2	3	4		6	7						13	15															27		
*38:35	1	2	3		5	6							13	15															27		
*38:36	1	2	3	4										15					19										27		
*38:37	1	2	3	4										15				18											27		
*38:41	1		3	4								12	13																		
*38:45	1	2	3		5									15																	
*38:46	1								10				13	15																	
*38:50	1	2	3		5									15					19												
*38:52, 38:65, 38:70	1	2	3	4										15															27		
*38:54	1	2		4									13	15															27		
*38:55Q	1	2	3	4									13	15								22							27		
*38:61	1	2	3	4									13	15														27	28		
*38:63	1	2	3	4										15				19													
*38:64	1	2	3		5								13	15															28		
*38:68Q	1	2	3	4									13	15														27	29		
*38:74	1	2			5								13	15																	
*38:81	1	2	3			6							13	15														27			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	

Negative Control

