

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
<b>HLA-B allele<sup>1,2,3</sup></b>																
*42:01:01:01-42:01:06, 42:22-42:24, 42:28-42:33	1	2				6			9							
*42:02:01:01-42:02:02		2	3			6			9							
*42:04	1			4		6			9							
*42:05:01-42:05:02, 07:349 <sup>5</sup>	1	2			5	6			9							
*42:06	1	2				6	7		9							
*42:07, 42:26, 08:94, 08:226 <sup>5</sup>		2				6		8								
*42:08	1	2				6			9	10						
*42:09		2	3				7		9							
*42:10	1	2				6			9		11					
*42:11		2				6		8				12				
*42:12	1	2				6			9				13			
*42:13	1					6			9					14		
*42:14, 42:25 <sup>4</sup>	1	2				6			9						15	
*42:15	1	2				6			9						14	
*42:16:01-42:16:02, 42:20	1	2					7		9							
*42:17		2	3			6			9	11						
*42:18		2	3			6			9					14		
*42:19, 07:388 <sup>5</sup>	1								9	11						
*42:21	1	2				6										
*42:27		2	3			6										
*07:02:01:01-07:02:05, 07:02:07-07:02:27, 07:02:29-07:02:97, 07:05:01:01-07:07:03, 07:09:01-07:09:02, 07:12, 07:14:01-07:15, 07:17-07:18:02, 07:20-07:24, 07:26, 07:28, 07:30-07:31:02, 07:34-07:36, 07:39-07:49N, 07:51-07:52, 07:54-07:55, 07:58-07:59, 07:61-07:64, 07:66-07:68:03, 07:70, 07:73-07:74, 07:76-07:82, 07:84, 07:86-07:99, 07:101-07:115, 07:117-07:124, 07:126-07:142, 07:144-07:145, 07:147-07:159, 07:161N-07:176, 07:178-07:179, 07:181N-07:182N, 07:184-07:185, 07:187-07:205, 07:208-07:218, 07:220-07:227, 07:229:01-07:239, 07:241-07:252, 07:254-07:258, 07:260-07:262:02, 07:264-07:291, 07:293-07:305, 07:307-07:312, 07:314-07:339, 07:341-07:348, 07:350-07:367, 07:369-07:382, 07:384-07:387, 07:389-07:399, 07:401-07:427, 07:429-07:441, 07:443-07:452, 07:454-07:462, 07:464-07:474, 14:21, 15:255, 15:509, 38:26, 54:21, 55:10, 56:16, 67:01:01-67:01:03, 67:03-67:07, 81:01:01-81:04N, 81:06-81:10, 83:01, C*01:240	1								9							
*07:03, 07:08:01-07:08:02, 07:16, 07:32, 07:38, 07:69, 07:85:01-07:85:02, 07:180, 07:207, 07:228:01-07:228:02, 07:306, 07:400, 13:14, 13:83, 14:01:01:01-14:01:05, 14:01:07-14:01:13, 14:07N-14:08:02, 14:10, 14:12, 14:14, 14:19, 14:26, 14:32, 14:40, 14:46-14:47, 14:49, 14:54-14:55, 14:65, 14:71, 14:77-14:78, 14:80-14:83, 14:86-14:87, 14:95, 14:97, 14:102-14:105, 14:108Q, 14:112-14:113N, 15:71, 15:175, 15:381, 15:556, 18:01:01:01-18:01:12, 18:01:14-18:01:50, 18:01:52-18:03:02, 18:05:01:01-18:13, 18:15, 18:17N-18:24, 18:26-18:36, 18:38-18:48, 18:50-18:55, 18:57:01-18:67, 18:69, 18:71-18:78, 18:80-18:106Q, 18:109-18:112, 18:114-18:124, 18:126-18:131:01:02, 18:133-18:141, 18:143-18:154N, 18:156:01:01-18:181, 18:183-18:196, 18:198-18:217, 18:219-18:226, 27:12:01:01-27:12:01:03, 27:16, 27:23, 27:92, 27:119, 27:153, 27:242, 35:50, 35:84, 35:162, 35:197, 35:231, 35:267, 35:269, 35:411, 35:423, 35:494, 37:02, 37:07, 37:25, 37:39-37:40, 37:66, 37:83, 37:94, 37:97, 37:102, 38:01:01:01-38:02:11, 38:04, 38:06-38:18, 38:20:01-38:25, 38:27-38:29, 38:31-38:32, 38:35-38:53, 38:55Q-38:83N, 38:85-38:88, 38:90-38:104, 38:160-38:173Q, 39:01:01:01-39:01:01:29, 39:01:03:01-39:01:37, 39:05:01:01-39:07, 39:09:01:01-39:11, 39:15-39:19:02, 39:22, 39:25N-39:28, 39:30-39:36, 39:38Q, 39:40:01N-39:48, 39:50-39:51, 39:53-39:62, 39:64-39:71, 39:73-39:75, 39:77-39:87N, 39:89-39:98, 39:100, 39:102-39:104, 39:106-39:109, 39:111-39:119, 39:121-39:123, 39:125-39:128, 39:130-39:143, 39:145-39:151, 39:153, 39:155, 39:157N-39:161N, 39:163, 39:165-39:167, 39:169-39:177, 39:179-39:180, 39:182-39:183, 39:186-39:194Q, 40:04:01:01-40:04:02, 40:06:01:01-40:06:14, 40:06:16-40:06:22, 40:06:24-40:06:31, 40:11:01:01-40:11:02, 40:20:01:01-40:20:01:02, 40:26, 40:27:03-40:28, 40:44, 40:64:01:01-40:64:02, 40:68, 40:70:01-40:70:02, 40:74-40:75, 40:83, 40:86, 40:95-40:96, 40:98-40:99, 40:109, 40:120, 40:127, 40:131, 40:148-40:149, 40:159, 40:161-40:162, 40:165, 40:167, 40:177, 40:184, 40:190, 40:210, 40:225:01-40:225:02, 40:230, 40:243-40:244, 40:256N, 40:268-40:270, 40:275, 40:294-40:295, 40:306-40:307, 40:311, 40:313, 40:317, 40:330, 40:340-40:341, 40:343, 40:348-40:351, 40:354, 40:361N-40:364, 40:366, 40:373-40:374, 40:376, 40:387, 40:389, 40:408, 40:411, 40:423, 40:425, 40:427, 40:429, 40:434, 40:445, 40:451, 40:453-40:454, 40:458-40:460, 40:464, 40:468, 40:471-40:473, 40:476, 40:478, 40:482, 40:484, 40:487N-40:488N, 40:494, 40:519-40:520, 40:526, 44:103, 44:131, 44:212, 44:224, 48:06, 48:31, 48:41, 49:31, 50:43, 51:57, 51:103, 51:116, 51:218, 52:105, 53:21								8								
*07:04:01-07:04:02, 07:25, 07:146	1	2			5		7		9							
*07:11, 07:57, 07:75:01:01-07:75:01:02, 07:160, 07:240, 07:253, 07:340, 07:463	1								9					14		
*07:19	1	2					7		9	10						
<b>Well No.</b>	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*07:27, 07:50, 07:453, 15:571, 15:582, 18:37:01-18:37:02, 27:18, 27:29, 38:03, 39:02:01:01-39:02:03, 39:08, 39:13:01-39:13:02, 39:23, 39:39:01-39:39:02, 39:49, 39:88, 39:101, 39:105, 39:110, 39:162, 39:164, 40:93, 44:130, 48:01:01:01-48:03:02, 48:05, 48:07-48:12, 48:14-48:21, 48:23, 48:25-48:30, 48:32-48:40, 48:42-48:50, 48:52, 48:54-48:56								8				12				

Negative Control

*07:29, 07:71, 07:116, 07:186, 07:313, 07:383, C*07:335										9							
*07:33:01-07:33:03, 07:53, 07:60, 07:100, 07:442	1									9	10						
*07:56:01-07:56:02, 07:65, 07:72, 07:177, 07:206, 07:263, 15:76, 15:101, 35:76, 35:561, 44:90, 44:356, 54:01:01:01-54:01:11, 54:03-54:06, 54:08N-54:20, 54:22-54:23, 54:25-54:29, 54:31-54:46, 55:01:01:01-55:03, 55:05, 55:07-55:09, 55:11-55:15, 55:17, 55:19, 55:21-55:24, 55:26-55:33, 55:35-55:46, 55:48, 55:50, 55:52-55:55N, 55:57-55:63, 55:65-55:77, 55:79-55:85, 55:87-55:109, 55:111-55:117N, 55:119-55:120, 55:122-55:134, 56:01:01:01-56:07, 56:09-56:12, 56:15, 56:18-56:22, 56:24-56:32, 56:34, 56:36-56:41, 56:43, 56:45-56:52, 56:55:01:02-56:66, 56:68-56:88, 56:90, C*03:542, C*14:62	1																
*07:143			3							9							
*08:01:01:01-08:01:22, 08:01:24-08:01:62, 08:01:64-08:03, 08:05, 08:08N, 08:10-08:11, 08:15, 08:18-08:19N, 08:21-08:24, 08:26:01-08:27, 08:29-08:36, 08:38:01-08:39, 08:41-08:46, 08:48, 08:50-08:52, 08:56:01-08:59:02, 08:61-08:68, 08:72N-08:73, 08:75, 08:78, 08:80-08:83, 08:85-08:86N, 08:88, 08:90-08:93, 08:95-08:102, 08:104-08:106, 08:109, 08:112-08:115, 08:117-08:122, 08:124-08:128, 08:130-08:131, 08:133, 08:135-08:142, 08:144, 08:146-08:155, 08:157-08:163, 08:165-08:170, 08:172-08:174, 08:176-08:180, 08:182-08:183, 08:185-08:186, 08:188-08:189:02, 08:191-08:199, 08:201-08:204, 08:206-08:225, 08:228-08:231, 08:233-08:244, 08:246-08:250, 08:252N-08:253, 08:256-08:263, 08:265-08:273, 08:275-08:280, 08:282-08:284Q, 08:286-08:289, 08:292-08:294, 08:297, 08:299-08:306, 35:87, 41:02:01:01-41:02:11, 41:11, 41:13, 41:19, 41:23, 41:27, 41:31, 41:36, 41:38-41:43, 41:45N-41:47, 41:55, 41:58, 41:65, 41:68-41:69, 41:71-41:72, 41:74, 44:166	2					6											
*08:01:23, 08:01:63, 08:14, 08:28, 08:74, 08:77, 08:87, 08:108, 08:111, 08:164, 08:227, 08:245, 08:251, 08:264, 08:290, 41:10, 41:49, 41:52	2																
*08:04:01-08:04:02, 08:54:01-08:54:02, 08:110, 08:285	2					6					12						
*08:07, 08:37, 08:200, 41:04, 41:18, 41:51	2						7										
*08:09, 08:84, 08:254, 15:83:01:01-15:83:01:02, 41:01:01:01-41:01:08, 41:05, 41:07, 41:09, 41:12, 41:14, 41:16-41:17, 41:20, 41:22, 41:25-41:26, 41:28-41:30, 41:32-41:35, 41:37, 41:44, 41:53-41:54, 41:56Q-41:57, 41:59-41:61, 41:64, 41:66-41:67, 41:70, 41:73, 41:75-41:79, 51:08:01:01-51:08:03, 51:08:05-51:08:06, 51:20, 51:36, 51:44N, 51:97, 51:141, 51:153:01-51:153:02, 51:263, 51:275, 51:342, 51:362, 51:371, 52:19				4			6										
*08:12:01-08:12:03, 08:184, 08:232, 41:03:01-41:03:02, 44:14, 44:62, 57:09							6										
*08:16, 08:60, 08:129, 08:181, 41:48							6									14	
*08:17	2					6					12					14	
*08:20:01-08:20:02, 08:281, 15:516, 15:576, 27:83											11						
*08:47, 08:76, 08:145, 08:295	2					6										14	
*08:49							7				11					14	
*08:53:01-08:53:02, 08:71, 08:296, 41:24	2					6					11						
*08:69	2					6	7										
*08:89, 08:107, 15:51, 15:179:01-15:179:02, 15:199, 15:218Q, 15:518, 35:38, 35:115, 35:169, 35:260, 35:545, 37:55, 37:60, 37:79N, 37:85, 37:99, 41:63, 44:17, 44:51, 44:123, 44:140, 44:230, 53:22, 57:24, 57:119, C*05:62, C*05:146, C*05:163, C*07:626, C*07:892, C*08:142, C*15:25								7									
*08:103, 08:134, 08:175, 08:274, 44:213	2					6					10						
*08:116, 08:298, 41:15	2					6											15
*08:123	2					6				9							14
*08:132	2					6				9							
*08:187	2				5	6											
*08:190	2					6		w									
*08:291							7				11						
*13:46, 41:21, 44:15:01:01-44:15:01:02, 44:18:01:01-44:18:01:02, 44:355, 45:01:01:01-45:02, 45:04-45:05, 45:07-45:13, 45:15-45:25, 45:27-45:28N, 45:30				4			7										
*13:61, 51:373, 52:92																	15
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
*13:71, 13:108, 14:05, 14:13, 15:07:01:01-15:07:03, 15:45, 15:55, 15:126, 15:197:02, 15:207, 15:324, 15:329, 15:331, 15:360, 15:383, 15:405, 15:421, 15:431, 15:450, 15:524, 27:07:01-27:07:06, 27:11, 27:24, 27:32-27:34, 27:43, 27:138, 27:180, 27:219, 27:233, 27:251, 35:05:01:01-35:05:06, 35:22:01:01-35:22:01:02, 35:31, 35:51:01-35:51:02, 35:58, 35:72, 35:89, 35:97, 35:114, 35:199, 35:232, 35:368, 35:389, 35:416, 35:427N, 35:440, 35:459N, 35:472, 35:481, 35:498-35:499, 35:509, 35:519, 35:541, 35:572, 39:76, 40:27:02, 40:57-40:58, 40:164, 40:200, 40:219, 40:226, 40:420, 40:457, 40:470, 40:496, 40:504, 41:62, 44:54, 44:135, 44:184, 46:12, 51:64, 51:81, 51:148, 53:14, 53:67, 58:18, 58:27, C*02:60, C*03:05:01:01-03:05:01:03, C*03:25, C*03:27, C*03:167, C*03:296:01-03:296:02, C*03:335, C*03:386, C*03:407, C*03:482, C*03:494, C*03:555, C*04:443, C*06:264, C*07:294, C*07:526:01, C*07:530:01, C*07:681, C*07:850, C*08:08:01-08:08:02										10							
*13:97, 14:73, 15:03:01:01-15:03:09, 15:03:11, 15:03:13, 15:47:01-15:47:02, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131-15:132, 15:151, 15:156, 15:158, 15:173, 15:210, 15:220:01:01-15:220:04, 15:235, 15:242:01-15:242:02, 15:253, 15:259, 15:266, 15:274, 15:281-15:282, 15:369, 15:379, 15:389, 15:397, 15:444, 15:448, 15:475, 15:502, 15:511, 15:517, 15:527, 15:539, 15:560, 15:588, 15:610, 15:614:01:01-15:614:01:02, 15:616Q-15:617, 15:628, 15:650, 18:125, 35:361, 40:12:01-40:12:02, 40:46, 40:430, 44:40, 44:44, 44:156, 44:210:01-44:210:02, 48:22, 48:57, 49:73, 50:75, 52:16, 55:18, 56:44																	12

Negative Control

*13:154, 41:50, 45:14				4		6	7										
*14:53, 18:14, 38:19, 39:03:01:01-39:03:02, 39:14:01:01-39:14:01:03, 39:24:01-39:24:03, 39:29, 39:37:01:01-39:37:01:02, 39:120, 39:144, 39:152, 39:185, 40:02:01:01-40:03:02, 40:05:01:01-40:05:02, 40:08:01:01-40:09, 40:13, 40:18-40:19, 40:24, 40:27:01, 40:29, 40:35:01-40:35:03, 40:37, 40:39-40:40, 40:50:01:01-40:50:01:02, 40:56, 40:71, 40:78, 40:82, 40:85, 40:89-40:91, 40:94, 40:97, 40:104-40:105, 40:107, 40:111, 40:115, 40:119, 40:122, 40:133Q, 40:142N-40:145, 40:157, 40:169, 40:173-40:174, 40:176, 40:180-40:181, 40:189, 40:201-40:203, 40:205-40:206, 40:209, 40:211, 40:214, 40:220, 40:224, 40:229, 40:232, 40:246, 40:248, 40:254-40:255, 40:266-40:267, 40:271, 40:274, 40:276, 40:283-40:284, 40:287, 40:289-40:293, 40:296-40:297, 40:302-40:305, 40:309, 40:314, 40:320, 40:322, 40:331, 40:334:01-40:335, 40:337N, 40:345N, 40:347, 40:356, 40:359, 40:369, 40:371, 40:375, 40:378, 40:384, 40:390, 40:399N, 40:405, 40:409-40:410, 40:412-40:415, 40:432, 40:435, 40:437-40:438N, 40:440-40:442, 40:444, 40:455, 40:465, 40:467, 40:481N, 40:489, 40:493, 40:500, 40:505, 40:511N, 40:514N-40:515, 40:518, 40:521-40:522, 40:527							8		10								
*15:68, 39:184, 48:04:01:01-48:04:02, 48:13, 48:24, 48:51, 48:53										8		10	12				
*15:530, C*01:59, C*01:118, C*01:157, C*06:82, C*06:210, C*07:49, C*07:124, C*07:155, C*07:210, C*07:238, C*07:247																	14
*18:49, 18:79, 18:108, 18:132, 18:142, 39:20, 39:124, 39:168										8							14
*18:56, 18:113, 37:01:01:01-37:01:29, 37:03N-37:04:02, 37:06:01-37:06:02, 37:08, 37:10-37:11, 37:13-37:24, 37:26-37:27, 37:29-37:38, 37:41-37:54, 37:56-37:59, 37:61, 37:63-37:65, 37:67-37:78, 37:80-37:82N, 37:86N-37:93, 37:95-37:96, 37:98, 37:100-37:101, 37:103-37:107								7	8								
*27:157, 40:166, 54:02, 55:16, 56:35				3													
*35:530, 44:342, 53:15												11					14
*37:05								7	8								14
*37:09	2							7	8		10						
*37:12, 37:84								6	8								
*37:28								7	8				12				
*37:62								7	w								
*38:30, 39:156								6	7	8							
*38:84, 39:99, 39:178										w							
*40:136, 40:231, 40:279, 41:08								5	6								
*41:06				4				6									15
*44:20, 44:100, 45:03, 45:29, 51:08:04, 51:274				4													
*44:106, 44:158	2										10						
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	
*44:117, 44:226, C*07:04:01:01-07:04:16, C*07:04:19-07:04:28, C*07:11-07:12, C*07:45, C*07:63, C*07:101, C*07:139, C*07:142, C*07:181, C*07:199:01-07:199:02, C*07:272, C*07:323-07:324, C*07:328-07:329N, C*07:336, C*07:338, C*07:354-07:355, C*07:357-07:358, C*07:361, C*07:364-07:365, C*07:378, C*07:394-07:395, C*07:406, C*07:420, C*07:426, C*07:428, C*07:459, C*07:466-07:467, C*07:480, C*07:487, C*07:501, C*07:523, C*07:534-07:535, C*07:562-07:563, C*07:569, C*07:585-07:586, C*07:600:01N-07:600:02N, C*07:622, C*07:651, C*07:655-07:656, C*07:664, C*07:672N, C*07:674, C*07:693, C*07:698, C*07:751N, C*07:780, C*07:797N, C*07:831, C*07:838-07:839N, C*07:852, C*07:858, C*07:868, C*07:876, C*07:881N, C*07:895-07:897, C*07:901, C*07:917, C*07:926, C*07:941, C*07:948, C*07:951, C*07:974Q, C*07:1010, C*07:1012, C*07:1033																	
*45:06			3	4				7									
*51:354				4			6	8									
*55:04, 55:49, 55:51, 55:86, 55:110	1									10							
*55:20	1			4			6										
*55:56, 56:13	1			4			6	7									
*82:01:01:01-82:04	1						6	7		9							
C*07:403								w									14
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	

Negative Control

Negative Control