

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-B allele^{1,2}																
*45:01:01:01-45:01:12, 45:15-45:21, 45:23, 45:25, 45:27-45:28N, 45:30	1	2	3													
*45:02	1	2			5											
*45:03	1	2		4												
*45:04		2	3			6										
*45:05	1	2	3				7									
*45:06	1	2	3					8								
*45:07	1	2	3	4												
*45:08		2	3			6			9							
*45:09	1	2								10						15
*45:10		2	3			6					11					
*45:11	1	2	3									12				
*45:12	1	2	3										13			
*45:13	1	2	3											14		
*45:14	1	2	3													15
*45:22, 44:15:01:01-44:15:01:02, 44:18:01:01-44:18:01:02, 44:355 ³	1		3													
*45:24		2	3													
*45:29, 50:02:01:01-50:02:01:02, 50:21, 50:56 ³	1	2														
*07:20, 07:140, 07:153, 15:14, 15:91, 15:131, 15:161, 15:551, 15:634, 18:169, 35:45, 35:71, 35:169, 35:452, 44:10, 44:14, 44:31, 44:43:01-44:43:02, 44:45, 44:51, 44:62, 44:77, 44:82, 44:107, 44:123, 44:135, 44:184, 44:213, 46:17, 57:07, 57:93, 58:07																15
*07:78, 07:84, 07:429, 13:02:01:01-13:02:36, 13:08-13:09, 13:14-13:16, 13:18-13:19, 13:27, 13:30-13:35, 13:37-13:38, 13:40-13:42, 13:44-13:45, 13:47, 13:49N, 13:54-13:56:02N, 13:58-13:59, 13:65-13:70, 13:72, 13:74-13:75, 13:81-13:82, 13:84-13:85, 13:88-13:89, 13:91, 13:93-13:94, 13:96-13:105, 13:111-13:119, 13:121, 13:123Q-13:128, 13:130, 13:132-13:134, 13:136-13:137N, 13:139N-13:140, 13:143, 13:145, 13:147, 13:149-13:153, 13:155-13:160, 13:165-13:167, 13:169, 13:172, 27:14, 27:81, 27:130, 35:534, 39:50, 40:93, 40:96, 40:109-40:110, 40:190, 40:340, 40:350, 49:07, 51:10, 51:16, 51:31, 51:34, 51:82, 51:93, 52:08, 52:67, 54:26, 55:09, 55:22, 55:24, 59:09, 73:01:01:01-73:03, 78:08						6					11					
*07:133, 40:166, 44:90, 44:356		2						8								
*07:225, 13:86, 35:02:01:01-35:02:07, 35:02:09-35:02:19, 35:81, 35:83, 35:95, 35:129N, 35:146, 35:157, 35:162, 35:172, 35:182-35:184, 35:211, 35:220, 35:258, 35:270, 35:285, 35:311, 35:323:01-35:323:02, 35:335, 35:339, 35:357, 35:361, 35:366-35:367, 35:372, 35:374, 35:377-35:379, 35:387, 35:391, 35:394, 35:396, 35:403, 35:407, 35:410, 35:419, 35:443, 35:468, 35:483, 35:487, 35:492, 35:501, 35:505, 35:508N, 35:522, 35:546, 35:550, 35:553, 35:557, 35:562, 53:19, 53:36, 58:65													13			
*08:09, 08:84, 08:254, 13:46, 13:53, 13:131, 14:37, 15:42, 15:83:01:01-15:83:01:02, 15:485, 35:60, 35:276, 39:06:01:01-39:06:07, 39:34, 39:57, 39:62, 39:64, 39:83:01-39:83:03, 39:90, 39:109, 39:127, 39:129, 39:132, 39:140-39:141, 39:154, 39:178, 39:186-39:187, 41:14, 41:60, 42:04, 46:18, 49:76, 51:21, 51:36, 51:101, 51:143, 51:153:01-51:153:02, 51:157, 51:197, 51:280, 51:295, 51:300, 51:340, 52:37, 52:71, 52:84, 54:01:01:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10-54:25, 54:27-54:32, 54:34-54:37, 54:39, 54:41-54:46, 55:01:01:01-55:01:06, 55:01:08-55:03, 55:05, 55:07, 55:10-55:13, 55:15-55:21, 55:23, 55:25-55:48, 55:50, 55:52-55:55N, 55:57-55:62, 55:64-55:85, 55:87-55:109, 55:111-55:119, 55:121, 55:123-55:129, 56:23, 56:54, 59:01:01:01-59:03, 59:05-59:08, 59:10N-59:11, 78:07, C*05:152, C*07:463, C*07:758, C*12:183, C*14:92						6			9							
*08:28, 08:35, 08:37, 08:89, 08:107, 08:291, 15:51, 15:179:01-15:179:02, 15:199, 15:218Q, 18:113, 35:115, 35:260, 37:35, 37:56, 42:09, 42:32, 44:308, 57:24, C*07:68, C*07:140:01-07:140:02, C*07:151, C*07:302, C*07:929										10						
*08:155, 15:46, 15:53, 15:106, 15:143, 15:212, 15:525, 18:48, 18:107, 35:19, 35:47, 35:63, 35:217, 35:274, 35:317, 35:450, 39:69, 40:01:01-40:05:02, 40:07-40:11:02, 40:14:01-40:16:01:02, 40:18, 40:20:01:01-40:20:01:02, 40:22N-40:40, 40:42-40:43, 40:45, 40:48-40:52, 40:54-40:69, 40:71-40:74, 40:76-40:78, 40:80-40:82, 40:84-40:85, 40:87:01-40:92, 40:94, 40:97-40:102, 40:104-40:108, 40:111-40:116, 40:118N-40:126, 40:128-40:130:02, 40:133Q-40:136, 40:138-40:147, 40:149-40:156, 40:158, 40:160:01-40:160:02, 40:163-40:164, 40:168-40:176, 40:179-40:184, 40:186:01-40:187, 40:189, 40:191-40:200, 40:202-40:206, 40:210-40:229, 40:231-40:243, 40:245:01-40:245:02, 40:247-40:255, 40:257-40:264, 40:266-40:267, 40:271-40:274, 40:276-40:280, 40:282-40:288, 40:290-40:291N, 40:293, 40:296-40:297, 40:299-40:305, 40:309-40:310, 40:312, 40:314-40:316, 40:319-40:329, 40:331-40:339, 40:343-40:347, 40:351-40:353:02, 40:355-40:360, 40:367-40:372N, 40:375, 40:377-40:386, 40:390-40:393, 40:395-40:407, 40:409-40:410, 40:412-40:422, 40:424, 40:428N-40:429, 40:431-40:433, 40:435-40:444, 40:446-40:448, 40:450, 40:452, 40:455-40:458, 40:461-40:463, 40:465-40:467, 40:469-40:470, 40:472, 40:474-40:475:02, 40:477, 40:479-40:481N, 40:483N, 40:485-40:486, 40:489-40:491, 40:493, 40:495-40:518, 40:521, 41:02:01:01-41:02:10, 41:03:01-41:04, 41:08, 41:10-41:11, 41:13, 41:15, 41:18-41:19, 41:23, 41:27, 41:31, 41:36, 41:38-41:43, 41:45N, 41:47-41:49, 41:52, 41:55, 41:58, 41:62-41:63, 41:65, 41:68-41:69, 41:71-41:72, 41:74, 44:09, 44:46:01-44:46:02, 44:75, 44:129, 44:131, 44:254, 44:277, 47:02, 48:41, 50:12, 50:17, 50:51		2														
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control

