

Well No.	1	2	3	4	5	6	7
HLA-B allele							
*47010101-47010102	1	2				6	
*4702			3		5	6	
*4703		2			5	6	
*4704		2		4			7
*4705	1	2					7
*1506, 152701-152703, 1584, 9509, 1827, 8201, 8202						6	
*1546, 1553, 9506, 9543, 3519, 3547, 3563, 400101-4011, 401401-4016, 4018, 4020, 4022N-4036, 4038-4040, 4042-4045, 4048-4078, 4080-4092, 4094, 4095, 4101-4108, 4409, 4446, 4501-4509, 5001, 5002, 5004			3		5		
*2701-270402, 270502-2706, 2708-2710, 2712, 2713, 2716-2718, 2720, 2723, 2726, 2727, 2729, 2731, 2735-2737, 2739-2742, 2744-2746, 2750, 3702							7
*3706	1	2					
*4013, 4019, 4411, 4418, 4425, 4450, 490101-490102, 4904, 4905		2					
*4037			3				
*4047, 4096, 44020101-4405, 4407, 4410, 4412-4417, 4419N-4424, 4426-4439, 4441-4443, 4445, 4447-4449, 4451-4456N, 4458N, 4459, 4461N-4467, 4469, 4902		2		4			
*4079					5		
Length of spec. PCR product	135	150	150	145	135	135	235
Well No.	1	2	3	4	5	6	7
Specific Reaction							
1st match: B*							
2nd match: B*							
Control Reaction							
Failed Controls							
False Positive							
False Negative							