

Well No.	1	2	3	4	5	6	7	8	9	10
HLA-B allele^{1,2,3}										
*47:01:01:01-47:01:02	1	2				6				
*47:02			3		5	6				
*47:03		2			5	6				
*47:04		2		4			7			
*47:05	1	2					7			
*47:06	1	2				6		8		
*47:07	1	2				6	7			
*47:08	1	2		4		6				
*47:09	1	2		4		6			9	
*07:133, 15:46, 15:53, 15:106, 15:143, 15:212, 18:48, 35:19, 35:47, 35:63, 35:154, 35:217, 39:69, 40:01:01-40:11:02, 40:14:01-40:16, 40:18, 40:20, 40:22N-40:40, 40:42-40:45, 40:48-40:78, 40:80-40:92, 40:94-40:95, 40:97-40:99, 40:101-40:108, 40:111-40:116, 40:118N-40:136, 40:138-40:156, 40:158-40:184, 40:186-40:187, 40:189, 40:191-40:200, 40:202-40:206, 40:210-40:245, 40:247-40:264, 40:266-40:280, 40:282-40:288, 40:290-40:291N, 40:293-40:295, 41:01:01-41:13, 41:15-41:33, 44:09, 44:46, 44:75, 44:90, 44:129, 44:131, 45:01:01-45:14, 50:01:01-50:02, 50:04-50:05, 50:07-50:19, 50:31-50:39			3							
*07:197, 27:01-27:05:15, 27:05:17-27:05:19, 27:05:21-27:06, 27:08-27:10, 27:12-27:13, 27:16-27:18, 27:20, 27:23, 27:26-27:27, 27:29, 27:31, 27:35-27:37, 27:39-27:42, 27:44-27:46, 27:48-27:61, 27:64N-27:69, 27:72-27:75, 27:77-27:80, 27:82-27:124, 27:126, 27:128-27:129, 27:131, 37:02, 48:21, 48:26							7			
*08:14, 27:15, 27:28, 27:62, 27:71, 44:127, A*02:298, A*02:406				4						
*15:06, 15:27:01, 15:27:03, 15:84, 15:109, 15:195, 15:327, 18:27					5	6				
*15:27:02, 82:01-82:02:02, A*23:23						6				
*15:162				4					9	
*35:23, 35:258, A*23:36, A*24:32					5					
*37:06:01-37:06:02, 40:188	1	2								
*40:13, 40:19, 40:109, 40:117, 40:292, 44:02:19, 44:05:03, 44:11, 44:18, 44:25, 44:50:01-44:50:02, 44:95, 44:152, 44:170, 44:196, 49:01:01-49:01:06, 49:04-49:20, 49:22-49:25, 49:27-49:31, 51:112, 53:30, 57:45, 57:51, 57:69		2								
*40:47, 40:96, 40:110, 40:157, 40:201, 44:02:01:01-44:02:12, 44:02:14-44:02:18, 44:02:20-44:05:02, 44:05:04, 44:07, 44:10, 44:12-44:17, 44:19N-44:24, 44:26-44:39, 44:41:01-44:43:02, 44:45, 44:47-44:49, 44:51-44:56N, 44:58N-44:59:02, 44:61N-44:74, 44:76-44:82, 44:84:01-44:89, 44:91-44:94, 44:96-44:126:02, 44:128, 44:132-44:133, 44:135-44:151, 44:153-44:155, 44:157-44:169, 44:171N-44:195N, 44:197-44:200, 44:202-44:204, 49:02		2		4						
*40:100, 50:20			3		5					
Well No.	1	2	3	4	5	6	7	8	9	10

Negative Control