

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
<b>HLA-B allele<sup>1,2,3</sup></b>																				
*48:01:01-48:01:06, 48:32-48:33, 48:35-48:38, 48:41	1	2	3								11									
*48:02:01-48:02:03	1			4							11									
*48:03:01-48:03:02, 48:40	1	2			5						11									
*48:04:01-48:04:02	1		3							10	11									
*48:05	1					6			9		11		13			16				
*48:06		2	3				7											18		
*48:07	1	2	3					8			11									
*48:08	1								9		11					16				
*48:09, 48:11 <sup>5</sup>	1	2	3								11			14						
*48:10, 39:13:02, 42:11 <sup>4</sup>	1										11		13							
*48:12	1										11		13		15					
*48:13	1									10	11		13							
*48:14	1										11		w	14	15					
*48:15	1	2	3								11					16				
*48:16		2	3					8			11									
*48:17	1	2			5			8			11									
*48:18, 48:29 <sup>5</sup>	1	2	3								11				15					
*48:19, 48:27 <sup>5</sup>	1	2	3						9		11									
*48:20	1	2	3			6					11									
*48:21	1	2			6						11									
*48:22, 48:30	1	2	3								11	12								
*48:23	1	2			5						11	12								
*48:24	1		3						9	10	11									
*48:25	1			4							11									19
*48:26	1	2			5	6					11									
*48:28	1	2	3													16				
*48:31	1	2	3																18	
*48:34	1	2	3								11								17	
*48:39		2	3								11									
*07:02:01-07:02:25, 07:02:27-07:02:38, 07:02:40-07:03, 07:08, 07:10, 07:15-07:16, 07:20-07:24, 07:28-07:31, 07:35, 07:37, 07:39, 07:41-07:42, 07:44N, 07:46-07:47, 07:49N, 07:55-07:59, 07:61-07:64, 07:67N, 07:70-07:76, 07:79, 07:82-07:83, 07:85-07:96:02, 07:98-07:99, 07:101-07:104, 07:106-07:111N, 07:113-07:114, 07:116-07:121, 07:124-07:130, 07:133-07:134, 07:136:01-07:136:02, 07:141-07:142, 07:144-07:145, 07:148-07:154, 07:156-07:161N, 07:166-07:168, 07:170-07:173, 07:175, 07:177, 07:179, 07:184-07:192, 07:194-07:195, 07:199-07:201N, 07:203-07:205, 07:208, 07:211-07:212, 07:215-07:218, 07:221, 07:224, 07:229-07:234, 07:236, 07:238, 07:240-07:247, 07:250-07:254, 07:256-07:257, 07:260, 07:263, 07:265-07:268, 07:272N-07:277, 07:280-07:281, 08:156							7		9							16				
*07:02:26, 07:02:39, 07:04, 07:26, 07:43, 07:45, 07:52, 07:66, 07:84, 07:131, 07:135N, 07:146-07:147, 07:178, 07:196, 07:198, 07:227, 07:239							7									16				
*07:05:01-07:06:02, 07:32, 07:69, 07:97, 07:105, 07:112, 07:138, 07:140, 07:176, 07:182N, 07:206, 07:210, 07:213, 07:249, 07:258, 07:264, 07:269-07:270, 07:278, 07:283, 08:79							7		9				13			16				
*07:07, 07:12, 07:14, 07:18:01-07:18:02, 07:77, 07:163, 07:193, 07:248					5		7									16				
*07:09, 07:11							7		9				14		16			w		
*07:13, 07:54, 07:65, 07:132, 07:181N, 07:259							7		9						16					
*07:17							7		9				14		16					
*07:19, 07:25, 07:48, 07:164, 07:223, 07:226, 08:01:12, 08:01:23, 08:07, 08:12:02, 08:14, 08:64, 08:74, 08:83, 08:87, 08:89, 08:107, 08:111, 14:01:01-14:03, 14:06:01-14:09, 14:11-14:12, 14:14-14:48, 14:50-14:52, 14:54-14:55, 38:02:06, 38:29, 39:01:19, 39:07, 39:15, 39:17, 39:25N, 39:43-39:44, 39:55, 39:91, 39:96, 42:16, 42:20, 50:15							7													
*07:27, 07:50							7		9		11					16				
*07:33:01, 07:60, 07:139							7		9	10										
*07:33:02-07:33:03, 07:51, 07:122, 07:155, 07:165, 07:220, 07:261							7		9											
*07:34, 07:80, 07:225							7						13		16					
*07:36, 07:38, 07:68:01-07:68:03, 07:81, 07:180							7		9						15	16				
*07:40, 08:20, 08:40, 08:70, 08:77							7		9				13							
*07:53							7		9	10			13							
*07:78							7	8					13		16					
*07:100							7		9	10					15					
*07:115							7	8	9							16				
*07:123									9				13		16					
*07:137					5		7						13		16					
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>
*07:143							7		9			12				16				

Negative Control

*07:162, 07:255					5	7									16			w		
*07:169						7		9							16	17				
*07:174, 07:222, 07:279	2	3				7		9							16					
*07:183						w		9							16					
*07:197					5	6	7								16					
*07:202	2	3				7		9												
*07:207						7		9				13			16		18			
*07:209, 08:01:01-08:01:11, 08:01:13-08:01:22, 08:01:24-08:02, 08:04:01-08:05, 08:08N-08:09, 08:11-08:12:01, 08:12:03-08:13, 08:15-08:16, 08:18-08:19N, 08:21-08:22, 08:24-08:27, 08:29-08:34, 08:36, 08:39, 08:41-08:48, 08:50-08:51, 08:53:01-08:53:02, 08:55-08:56, 08:58-08:59, 08:72N-08:73, 08:75-08:76, 08:80-08:82N, 08:85-08:86N, 08:88, 08:90-08:100, 08:102-08:106, 08:108-08:110, 08:112-08:115, 08:117-08:132, 08:135-08:146, 08:149-08:155, 08:157-08:166, 38:02:01-38:02:05, 38:02:07, 38:08, 38:15, 38:17-38:18, 38:23, 38:35, 38:43-38:48, 38:50, 38:62, 38:64, 39:01:01-01-39:01:01:04, 39:01:03-39:01:18, 39:01:20-39:01:23, 39:04-39:06:05, 39:09:01-39:10:01, 39:11-39:12, 39:16, 39:18-39:20, 39:26-39:28, 39:30-39:36, 39:38Q-39:41, 39:45-39:48, 39:50-39:54, 39:57-39:75, 39:77, 39:79-39:90, 39:92-39:95N, 39:97N, 39:99-39:100, 39:102-39:104, 39:106-39:108, 39:111-39:119, 39:121, 42:01:01-42:01:03, 42:04-42:07, 42:12-42:13, 42:15, 42:19, 42:21-42:24, 67:01:01-67:01:03, 67:03-67:07						7						13								
*07:214					5	7	8					13			16					
*07:219								9							15	16				
*07:228						7		9							16		18			
*07:235, 08:49, 14:10			4			7														
*07:237						7		9					14	15	16			w		
*07:262						7		9				13			15	16				
*08:03, 08:52, 08:71, 08:78, 38:01:01, 38:01:03-38:01:12, 38:05-38:07, 38:09-38:14, 38:16, 38:21, 38:24, 38:26-38:28, 38:30-38:34N, 38:36-38:37, 38:39-38:42, 38:51, 38:53-38:55Q, 38:57-38:61, 38:65, 42:10						7						13			15					
*08:10, 08:147-08:148N, 15:09, 15:229, 15:243, 15:252, 15:303, 15:361, 27:20, 27:27, 27:130, 27:153, 39:10:02, 39:78, 40:04, 40:06:01:01-40:06:09, 40:06:11-40:06:15, 40:70:01-40:70:02, 40:99, 40:244, 40:268-40:270, 44:166, 46:11, 46:18, 46:61, 51:01:01:01-51:01:18, 51:01:20-51:01:30, 51:01:32-51:01:49, 51:01:51-51:02:02, 51:02:04-51:03, 51:05-51:06:03, 51:08:01-51:16, 51:19-51:24:05, 51:26, 51:28-51:36, 51:38-51:41N, 51:43-51:44N, 51:48-51:55, 51:57-51:63:02, 51:65-51:74, 51:76-51:80, 51:82-51:83, 51:86-51:89, 51:91, 51:95-51:115, 51:117-51:121, 51:123-51:129, 51:131-51:134, 51:136-51:138, 51:140-51:145, 51:147, 51:149N-51:150, 51:152-51:162, 51:164-51:171, 51:173Q-51:176, 51:178N-51:184N, 51:186-51:194, 51:196, 51:198-51:205, 51:207-51:213, 52:01:13, 52:06:01-52:06:03, 52:22-52:23, 54:01:01-54:05N, 54:07-54:08N, 54:10, 54:12-54:13, 54:17-54:21, 54:23-54:26, 54:28-54:37, 55:01:01-55:03, 55:05, 55:07-55:13, 55:15-55:22, 55:24-55:26, 55:28-55:40, 55:42-55:45, 55:47-55:48, 55:52-55:57, 55:59-55:67, 55:69-55:78, 55:80, 56:01:01:01-56:02:01, 56:04-56:08, 56:10, 56:13-56:17, 56:19N-56:30, 56:33-56:49, 56:51, 57:82, 58:08:01-58:08:02, 58:20, 59:01:01:01-59:09, 67:02, 78:01:01-78:01:02, 78:03, 78:07-78:08, 82:01-82:03, C*01:90, C*01:116, C*02:68, C*02:125, C*03:317, C*06:21, C*15:24, C*15:91, C*15:128												13								
*08:17, 08:38, 08:54:01-08:54:02, 08:101, 38:03-38:04, 39:02:01-39:02:02, 39:08, 39:13:01, 39:22, 39:101, 39:105, 39:110, 40:317, 41:01:01-41:03:02, 41:05-41:07, 41:09, 41:11-41:14, 41:16-41:17, 41:19-41:20, 41:22-41:23, 41:26-41:36, 41:38-41:45N, 41:47-41:48, 45:12, 45:14, 50:36	1														13					
*08:23, 08:133						7						13					18			
*08:28, 08:37						7						14								
*08:35						7						13	14							
*08:84, 38:49, 39:56, 39:98, 39:109						7						w								
*08:116, 08:134, 42:08, 42:14						7		10				13								
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>
*13:01:01-13:01:12, 13:06-13:07N, 13:12:01-13:13:02, 13:17, 13:20, 13:25-13:26:02, 13:28-13:29, 13:36, 13:39, 13:43, 13:50-13:52, 13:57, 13:60-13:61, 13:63N, 13:73, 13:76N-13:80, 13:83, 13:87, 13:90, 13:92, 13:95, 15:02:01-15:02:11, 15:13:01-15:13:02, 15:20, 15:25:01-15:25:03, 15:36, 15:77, 15:85, 15:88-15:89, 15:106, 15:112, 15:121, 15:144, 15:154, 15:165, 15:170, 15:213-15:214, 15:240, 15:250, 15:265, 15:271, 15:283, 15:289, 15:291, 15:301-15:302N, 15:308, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:357-15:358, 15:366, 15:374, 15:378, 15:393, 15:399, 35:20:01-35:20:02, 35:28, 35:226, 35:296, 37:55, 44:60, 44:83-44:84:01, 44:194, 44:205:01-44:205:02, 44:225, 53:09, 53:17:01-53:17:02, 53:38, 53:40				4							11									

Negative Control







*40:19, 40:292	?								10		13	?								
*40:21, 40:137	?	2			5					11										
*40:22N, 40:33, 40:38, 40:42-40:43, 40:46-40:47, 40:49, 40:54-40:55, 40:61, 40:63, 40:67, 40:74, 40:76, 40:81, 40:84, 40:87:01-40:88, 40:100-40:102, 40:112-40:114:02, 40:116, 40:118N, 40:123, 40:124:02-40:126, 40:128, 40:132, 40:134-40:135, 40:141, 40:149, 40:151, 40:154, 40:156, 40:160:01-40:160:02, 40:168, 40:170-40:171, 40:175, 40:179, 40:182, 40:186:02-40:188, 40:191-40:197, 40:199, 40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:218, 40:221-40:223, 40:227-40:228, 40:231, 40:234-40:236, 40:239-40:242, 40:245, 40:247, 40:249, 40:257, 40:273, 40:279, 40:281-40:282, 40:288, 40:299-40:300, 40:308, 40:310, 40:315, 40:319, 40:324-40:325, 40:327-40:329, 40:333, 40:338N-40:339, 40:344	?	2		5																
*40:23	1	2			5													16		
*40:25, 40:166		2			5	?														
*40:27:01-40:27:02, 40:29, 40:37, 40:50, 40:56-40:58, 40:78, 40:82, 40:90, 40:104, 40:107, 40:111, 40:115, 40:119, 40:133Q, 40:145, 40:157, 40:164, 40:169, 40:173-40:174, 40:176, 40:181, 40:189, 40:200-40:203, 40:205-40:206, 40:211, 40:219, 40:224, 40:226, 40:246, 40:274, 40:276, 40:283, 40:287, 40:289, 40:291N, 40:293, 40:296-40:297, 40:320, 40:322, 40:331, 40:334-40:335, 40:337N	?								10		13									
*40:28	?			4								13								
*40:30, 40:34, 40:233, 40:298	?	2																		
*40:31, 40:45, 40:80	?	2	3																	
*40:32, 40:98	?								9									16		
*40:35:01-40:35:02, 40:71	?								10		13	15								
*40:36, 40:48, 40:280	?	2			5		8													
*40:44, 40:92, 40:120, 40:130:01-40:131, 40:159, 40:163, 40:178, 40:307, 40:330, 40:336, 40:341	?																			
*40:52	?	2			5						12									
*40:59	1	2			5						12									
*40:62, 40:65, 40:108, 40:238, 40:312, 40:326	?				5															
*40:68						?						13						18		
*40:72:01-40:72:02	1				5							w								
*40:77	1				5		8					w								
*40:83	?									11	12	13								
*40:105	?								10			13	14							
*40:106		2			5	?												18		
*40:109, 40:340	?											13	?							
*40:117	?	2			5								?							
*40:158	?				5							13		16						
*40:162	?											13						18		
*40:183, 40:198	?	2			5									16						
*40:185	?	2			5					11	12									
*40:209, 40:248	?								10	11	12	13								
*40:220	?								10			13						18		
*40:232						?			10			13						18		
*40:290	?								10			w								
*40:305									10			w								
*40:313						?							13							
*41:08	1				5							13								
*41:15	1								10			13								
*41:24	1											13	14	15						
*41:25, 41:49	1											w								
*42:02:01:01-42:02:01:02, 42:18						7					12	13								
*42:09						7					12	14								
*42:17						7					12	13	15							
*44:18, 49:01:01-49:01:08, 49:05-49:14, 49:16-49:25, 49:27-49:46	1												15							
*44:110, 44:211				4					10											
*44:150				4	5													16		
*46:06		2			5						12									
*46:12									9	10			14	15						
*49:03	1										11	12		15						
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>
*49:04:01-49:04:02	1														15				19	
*50:17	1			4																
*51:01:19, 51:75, 51:206, 56:02:02, A*26:01:19													w							
*53:14, 58:18										10				14	15					
*57:39				4								12	13							
*57:42				4								w								
*73:01-73:02						w						12								

Negative Control

Control

