

HLA-B allele ^{1,2,3}	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	
*52:01:01:01-52:01:28, 52:01:30-52:01:32, 52:01:34-52:01:50, 52:22-52:24, 52:33-52:35, 52:37-52:38, 52:40-52:42, 52:44-52:48, 52:50-52:52, 52:54-52:61, 52:64-52:66, 52:69-52:70, 52:72, 52:76, 52:78-52:87, 52:89N-52:96N, 52:98-52:107, 52:109-52:110N	1		3									12	13									21		
*52:01:29, 52:01:52, 52:39, 52:62, 52:68, 52:97	1		3									12	13											
*52:01:33, 52:28-52:29			3									12	13									21		
*52:01:51, 52:108	1		3									12	13									21		
*52:02:01-52:02:02	1			4								12	13									21		
*52:03			3		5	6						12	13									21		
*52:04, 52:09:01-52:09:02 ⁴	1	2	3									12	13									21		
*52:05	1		3		5	6						12	13									21		
*52:06:01	1		3						9		11		13									21		
*52:06:02-52:06:03	1		3								11		13									21		
*52:07:01-52:07:02	1		3									12		14								21		
*52:08	1		3							10		12	13									21		
*52:10:01-52:10:03, 52:71, 52:74			3			6						12	13									21		
*52:11:01-52:11:02			3		5							12	13									21		
*52:12, 52:26, 52:63	1		3									12	13		15									
*52:13	1		3				7					12	13									21		
*52:14:01-52:14:02	1		3						9			12	13									21		
*52:15	1		3		6							12	13									21		
*52:16	1											12	13								20	21		
*52:17	1		3									11	12	13								21		
*52:18	1								10			12	13									21		
*52:19			3				7					12	13									21		
*52:20	1		3					8				12	13									21		
*52:21:01-52:21:02	1		3		5							12	13									21		
*52:25:01-52:25:02	1		3									12	13										22	
*52:27, 52:77	1		3									12	13								20	21		
*52:30	1		3									12	13									21		
*52:31:01-52:31:02, 51:310 ⁵	1		3									12	13		16							21		
*52:32	1											12	13									21		
*52:36	1		3	4								12	13									21		
*52:43	1		3									12	13									21		
*52:49N	1		3						9			12	13									21		
*52:53	1		3									12	13									21		
*52:67			3		5					10		12	13									21		
*52:73, 52:111	1		3									13										21		
*52:75	1		3						9			12	13		15									
*52:88, 53:28 ⁵	1		3									12	13		16									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	
*07:02:01:01-07:02:05, 07:02:07-07:02:10, 07:02:12-07:05:02, 07:05:04-07:26, 07:28-07:37:02, 07:39-07:49N, 07:51-07:57, 07:59-07:77, 07:79-07:83, 07:85:01-07:217, 07:219-07:399, 07:401-07:404, 07:406-07:428, 07:430-07:452, 07:454-07:474, 08:01:01:01-08:01:68, 08:05, 08:07-08:08N, 08:10-08:16, 08:18-08:51, 08:53:01-08:53:02, 08:55-08:77, 08:79:01-08:83, 08:85-08:109, 08:111-08:116, 08:118-08:125, 08:127-08:190, 08:192-08:254, 08:256-08:270Q, 08:272-08:273, 08:275-08:284Q, 08:286-08:297, 08:299-08:306, 13:01:15, 13:07N, 13:36, 13:39, 14:41N, 14:43, 14:46, 14:56, 14:59, 14:62, 14:90, 14:95, 15:01:01:01-15:01:01:65Q, 15:01:03-15:01:04, 15:01:06-15:02:12, 15:03:12, 15:05:01:01-15:08:02, 15:10:01:01-15:12:03, 15:14-15:15:01:02, 15:18:01:01-15:19, 15:21:01:01-15:21:01:02, 15:23, 15:25:01:01-15:35:02, 15:39:01:01-15:40:02, 15:42, 15:44-15:46, 15:48-15:53, 15:55-15:58, 15:60, 15:63, 15:65-15:66, 15:70-15:82, 15:84-15:86, 15:88-15:90, 15:92-15:94N, 15:96-15:97, 15:99, 15:101, 15:104-15:122, 15:124-15:126, 15:128-15:129, 15:133-15:136, 15:138-15:150, 15:152-15:154, 15:159-15:161, 15:163-15:167, 15:169-15:172, 15:174-15:176, 15:178-15:195, 15:197:01-15:201, 15:203-15:207, 15:209N, 15:211-15:215, 15:217-15:219, 15:221, 15:223-15:227, 15:229, 15:231-15:234, 15:236-15:238, 15:240-15:241, 15:244-15:250, 15:255, 15:257-15:265, 15:267, 15:269-15:272N, 15:275:01-15:276, 15:279-15:280, 15:283, 15:285-15:292, 15:294N-15:304N, 15:306-15:309, 15:311-15:322, 15:324-15:334, 15:336-15:338, 15:340-15:355, 15:357-15:360, 15:363:01-15:363:02, 15:365-15:367, 15:370-15:375N, 15:377Q-15:378, 15:380N-15:388, 15:390-15:395, 15:398-15:402, 15:404-15:407, 15:409-15:410, 15:412-15:417, 15:419-15:422, 15:425-15:428, 15:430-15:434, 15:436-15:441, 15:443, 15:445, 15:447, 15:449-15:461, 15:463N-15:474, 15:476-15:498, 15:501, 15:503-15:515, 15:518-15:522, 15:524-15:526, 15:528N-15:530, 15:533:01-15:538, 15:540N-15:544N, 15:547-15:549N, 15:551-15:554, 15:557-15:559, 15:561-15:567, 15:569-15:570, 15:572-15:574, 15:576-15:581, 15:583-15:585, 15:587, 15:589-15:598, 15:600-15:601, 15:604N-15:608, 15:611-15:612, 15:615, 15:618, 15:620-15:625, 15:627, 15:629-15:641, 15:643, 15:646-15:649, 15:651-15:652, 18:04:01-18:04:02, 18:07:01-18:07:03:02, 18:29, 18:33, 18:79-18:80, 18:98:01-18:98:02, 18:102, 18:107, 18:125, 18:155, 18:182N, 18:218, 27:01-27:04:02, 27:04:04-27:05:30, 27:05:32-27:13:02, 27:15-27:21:02, 27:24-27:28, 27:30-27:80, 27:82-27:101, 27:103-27:129, 27:131-27:152, 27:154-27:196, 27:198-27:207, 27:209-27:214, 27:216-27:250, 27:252-27:265, 35:31, 35:32:03, 35:35, 35:43:01, 35:43:03-35:44, 35:46, 35:57, 35:60, 35:79, 35:87, 35:135, 35:185, 35:208:02, 35:279, 35:475-35:476, 37:60, 37:67, 38:01:01:01-38:02:11, 38:04-38:05, 38:08-38:72, 38:75-38:104, 38:160-38:165N, 38:167-38:173Q, 39:01:01:01-39:01:01:29, 39:01:03:01-39:01:38, 39:03:01:01-39:05:05, 39:07, 39:09:01:01-39:12:01:02, 39:14:01:01-39:20, 39:22, 39:24:01-39:32, 39:35-39:38Q, 39:40:01N-39:45, 39:47-39:48, 39:51, 39:56, 39:58, 39:61, 39:62, 39:65, 39:82, 39:84, 39:87N, 39:89, 39:91																								

Negative Control

*51:40:01-51:40:02, 51:108					5	6					11	13									21		
*51:54					5	w	8				11	13									21		
*51:61:01-51:61:02	1				5						11	13									21		
*51:62, 51:106:01-51:106:02, 51:223, 51:225	1										11	13									22		
*51:65	1										11	13			17						21		
*51:69	1										11	13			18						21		
*51:78:01-51:78:02, 51:341	1						8				11	13									21		
*51:92:01-51:92:02	1							9			11	13	15								21		
*51:93					5	6			10		11	13									21		
*51:97							7				11	13									21		
*51:104:01-51:104:02	1		4								11	13									21		
*51:139						6					11	13		16							21		
*51:185	1	2									11	13									21		
*51:197											11	13	15								21		
*51:217						w					11	13									21		
*51:227, 51:281	1										11	13		16							21		
*51:237:01-51:237:02	1										11	14									21		
*51:274							7				11	13	15								21		
*51:340						6						13									21		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
*51:359, 53:01:28, 53:59, A*24:328											11												
*51:362	1						7				11		13									21	
*51:371							7	8			11		13										
*53:05:01-53:05:02, 53:16, 53:33		2									11		13										
*53:09, 58:137								8					13										
*53:11-53:13, 53:31, 53:36, 53:67								8			11		13										
*53:17:01-53:17:02			3									12	13										
*53:38			3					8				12	13										
*56:06, A*33:231, C*01:30, C*02:142, C*03:384, C*05:216:01-01-05:216:01:02, C*06:207, C*08:51, C*08:114, C*08:243, C*08:246, C*12:87, C*14:76, C*15:39, C*15:130, C*15:195	1																						
*56:34														14									22
*57:01:11, 57:91				4																			
*57:01:32				4	5	6								14	16								
*57:02:01:01-57:03:05, 57:05, 57:07, 57:09, 57:17, 57:39, 57:42, 57:46, 57:57, 57:63, 57:66, 57:70, 57:80, 57:82, 57:84, 57:94, 57:96, 57:101, 57:124, 57:137, 57:155				4										14								21	
*57:60		2		4										14									
*58:27-58:28:02, 58:65													13									21	
*58:36				4									13										
*58:92							7						13										
*59:03									9		11			14									
*59:09										10	11			14									
*73:01:01:01-73:01:01:02, 73:03									9	10			w										
*73:02									9	10													
*78:05-78:06	1		3											13									
*78:08							6				10			13	14								
*78:10							6							13									
A*23:129, A*24:146																						21	
A*68:01:52																16							
C*02:02:23, C*12:03:31, C*14:02:37, C*14:92						5																	
C*03:464, C*07:463, C*07:758, C*12:183															15								
C*05:152									9														
C*15:02:53						5										16							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23

Negative Control