

| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 |
|---|---|---|---|---|---|---|---|---|---|----|----|----|----|----|----|----|
| HLA-B allele | | | | | | | | | | | | | | | | |
| *530101-530104 | 1 | | | | | | | | 9 | 10 | 11 | | | | | 16 |
| *5302 | 1 | 2 | | | | 6 | | | 9 | 10 | 11 | | | | | 16 |
| *5303 | | | 3 | | | | | | 9 | 10 | 11 | | | | | 16 |
| *5304 | 1 | | | 4 | | | | | 9 | | | | | | | 16 |
| *5305 | 1 | | | | 5 | 6 | | | 9 | 10 | 11 | | | | | 16 |
| *5306 | 1 | 2 | | | | 6 | | 8 | | | 11 | | | | | 16 |
| *5307 | 1 | | | 4 | | | | 7 | 9 | | | | | | | 16 |
| *5308 | 1 | | | | | | | 8 | | | 11 | | | | | 16 |
| *5309 | | | | | 5 | | | | 9 | 10 | 11 | | | | 15 | 16 |
| *5310 | 1 | | | | | | 7 | | 9 | 10 | 11 | | | 14 | | 16 |
| *5311 | | | | | 5 | | | | 9 | w | 11 | | | | | |
| *5312 | | | | | 5 | 6 | | | 9 | 10 | 11 | | | | | 16 |
| *5313 | | | | | 5 | | | | 9 | 10 | 11 | | | | | 16 |
| *5314 | 1 | | | | | | | | 9 | 10 | 11 | | | 14 | | 16 |
| *5315 | 1 | | 3 | | | | | | 9 | 10 | 11 | | | | | 16 |
| *5316 | 1 | | | | 5 | | | | 9 | 10 | 11 | | | | | |
| *070201-0708, 0710, 0713, 0715, 0716, 0719-0724, 0727-0733, 0735-0742, 0744-0747, 0749N, 0750, 0752-0766, 0806, 0820, 1403, 1530, 1533, 1545, 1563, 1577, 1578, 1599, 9516, 9541, 2725, 3544, 3918, 3936, 4015, 4016, 4023, 4032, 4805, 4808, 5709, 5713 | | | | | | | | 8 | | | | | | | | |
| *0709, 0711, 0717 | | | | | | | | 8 | | | | | | 14 | | |
| *0712, 0714, 071801-071802 | | | 3 | | | | | 8 | | | | | | | | |
| *0726, 0734, 0743, 1320, 1401-140202, 1404, 140601-1409, 1509, 1510, 15170101-151702, 1537, 1572, 1590, 270401-270402, 2706, 2710, 2715, 2718, 2720, 2721, 2724, 2740, 3810, 3935, 4005, 4026, 4051, 4815, 5150, 5202, 5517, 5528, 560501-5606, 5621, 7801-7803 | | | | | | | | 8 | | | | | | | | 16 |
| *0748, 0813, 1301, 1306, 1307N, 1311-1313, 1317, 1321, 1322, 1586, 2701-2703, 270502-270510, 2707-2709, 2711-2714, 2716, 2717, 2719, 2723, 2726-2737, 2739, 2742, 2743, 3531, 3533, 3565Q, 3588, 380101-3805, 3808, 3809, 3811-3817, 39010101-39010102L, 390103-390202, 3904-390602, 3910, 3912-3917, 391901-3920, 3922, 3923, 3925N-3931, 3934, 3938Q-3942, 3944, 3945, 400101-400203, 4004, 40060101-401403, 4018, 4019, 4021, 4022N, 4024, 4025, 4027, 4029, 4031, 4033, 4035-4037, 4039, 4040, 4042-4050, 4053-4058, 4061, 4062, 4064-4070, 407201-4091, 4431, 446401-446402, 47010101-4705, 480101-480102, 480301-4804, 4806, 4809-4813, 4816, 4818, 5404, 5411, 5414, 5504, 5513, 5415, 5527, 5610, 5612, 5902, 5903, 670101-6702, 7301, 8101N-8104 | | | | | | | | | | | | | | | | 16 |
| *0821, 3518, 44020101-440205, 4405, 4408, 4409, 4411, 4412, 4414, 4419N, 4420, 4422-4425, 4427-442802, 4433-4435, 4441, 4442, 4444, 4448, 4449, 4451-4453, 4455, 4458N, 4459, 4462, 4463, 5702, 5704, 5705, 5712, 8301 | | | | | | | | | 9 | | | | | | | |
| *0825, 1548, 2738, 350201-350202, 350401-350403, 350901-350902, 3512, 3522, 3534, 3539, 3583, 4063, 4092, 440301-440303, 4407, 4413, 4426, 4429, 4430, 4432, 4436-4440, 4445-4447, 4450, 4454, 4457, 4460, 4461N, 5508, 5615, 570301-570302, 5707, 7804 | | | | | | | | | 9 | | | | | | | 16 |

| | | | | |
|---|---|-----|---------|-------|
| *0828, 0835, 4206, 4209 | | | | 14 |
| *130201-130204, 1308Q, 1309, 1314, 1315, 1319, 3560, 5401, 5402, 5405N, 5407, 5408N, 5412, 5413, 550201-550202, 5507, 5510, 5512, 5516, 5518, 5519, 5522, 5526, 5530, 5619N, 5623, 5901 | | | 12 | 16 |
| *1303, 4611, 5403, 5601, 5607, 5608, 5614, 5616, 5617, 5620 | 4 | | 9 12 | 16 |
| *1304, 1505, 1520, 1591, 9523, 9536, 350101-350107, 3507, 3510, 3516, 3517, 3519, 3523, 3525-3527, 3529, 3530, 3532, 3537, 3540N-3542, 3546-3550, 3552, 3554, 3557, 3563, 356801-3569, 3571, 3576-3578, 3582, 3590-3592, 3594, 4443, 4613, 4802, 5609, 570101-570102, 5706, 5708, 5710, 5711, 5715, 5716 | | | 9 10 11 | 16 |
| *1310 | | | w w | 16 |
| *1316, 1542, 550101-550104, 5503, 5505, 5511, 5515, 5524, 5525, 5529, 5531 | | 8 | 12 | 16 |
| *1318, 4508, 5410 | | | 12 | |
| *1405 | | 8 | | 14 16 |
| *15010101-150104, 150106-150108, 1504, 1508, 151101-1512, 1514, 1515, 1519, 1524, 1526N-1528, 1534, 1535, 1543, 1546, 1547, 1549, 1550, 1553, 1554, 1556, 1557, 1560, 1565, 1566, 1570, 1571, 1576, 1579N, 1581, 1582, 1585, 1592, 1594N, 1596, 1597, 9501, 9502, 9504, 9509, 9510, 9513, 9517, 9518, 9520, 9522, 9525, 9535, 9540, 9542-9544, 9546-9549N, 351401-351402, 3543, 3567, 3579, 3586, 460101-4605, 4607N, 4610, 4614-4616, 5406, 5603 | | 8 | 11 | |
| *1502, 1506, 1513, 1521, 1525, 1536, 1544, 1589, 9506, 9512, 9521, 4905 | | 8 | w | 16 |
| *1503, 1516, 1518, 1523, 1529, 1539, 1540, 1561, 1562, 1564, 1567, 1569, 1574, 1580, 1593, 1595, 1598, 9503, 9508, 9519, 9527, 9531, 9532, 9534, 9538, 9539, 3511, 4608, 4904, 5514 | | 8 | 11 | 16 |
| *1507, 1568, 9526, 4612 | | 8 | 11 | 14 |
| *1531 | | | 9 w w | 16 |
| *1532 | | 7 8 | 11 | |
| *1538, 1587 | | 6 8 | 11 | |
| *1551, 3562, 4417, 4509, 4609 | | | 9 11 | |
| *1552, 9514, 9524, 2741, 3515, 3535, 3593, 3907, 5409, 5414, 5523, 5532, 5618, 570103 | | | 10 11 | 16 |
| *1555 | | | 9 w w | 14 16 |
| *1558, 1573, 9537 | 4 | 8 | | |
| *1575, 9505, 9511N, 9545 | | | 11 | |
| *1584, 350801-350803, 3545, 3561 | | | 9 w 11 | |
| *1588 | | | w | 16 |
| * 9507, 3553N, 3564 | | | 9 10 | 16 |
| *9515 | | 6 8 | 11 | 16 |
| *9528, 9529 | | 8 | w | |
| *9533 | 4 | 8 | | 16 |
| *180101-1806, 1808, 1812, 1817N, 1823N-1825, 1827, 1828 | | 6 | 10 11 | 15 16 |
| *1807, 1929, 3943 | | 6 | 10 11 | 16 |
| *1809 | | 5 6 | 10 11 | 15 16 |
| *1810, 1811 | | | 10 11 | 15 16 |
| *1813 | | 6 | w 11 | 15 |

| | | | | |
|---|-----|-----|---------|----------|
| *1814 | | 6 | 10 11 | 14 15 16 |
| *1815 | | 6 8 | 11 | 15 16 |
| *1818 | | 6 7 | 10 11 | 15 16 |
| *1819 | | 6 8 | 11 | 15 |
| *1820 | | w | w w | 15 16 |
| *1821 | | 8 | | 15 |
| *1822 | 2 | 6 | 10 11 | 15 16 |
| *1826 | | | | 15 16 |
| *3503, 3506, 3513, 3536, 3555, 3556, 3559, 3570, 3574, 3575, 3584, 3585, 3932, 4071, 4410, 5602, 5604, 5611, 5622 | 4 | | 9 | 16 |
| *3505, 3551 | | | 9 10 11 | 14 16 |
| *352001-352002, 3528 | | | 9 10 11 | 15 16 |
| *3521 | 2 | 6 8 | 11 | 16 |
| *3524, 5714 | 2 | 6 | 9 10 11 | 16 |
| *3538, 3933, 4502, 8201, 8202 | 4 | | 9 | |
| *3558 | | 8 | 11 | 14 16 |
| *3566 | | | w 11 | 14 |
| *3572 | | | 9 10 11 | 14 15 16 |
| *3580 | | | 9 10 11 | |
| *3581 | 2 | | 9 | 16 |
| *370101-370104, 3703N, 3708 | 3 | 7 | | 15 |
| *3702 | 3 | | | 15 16 |
| *3704 | 2 3 | 7 | | 15 |
| *3705, 3711, 3714 | | 7 | | 15 |
| *3706 | | 7 | | |
| *3707 | 3 | 8 | | 15 |
| *3709, 3712 | 3 | | | 15 |
| *3710 | | 5 7 | | 15 |
| *3713 | 3 | 7 8 | | 15 |
| *3806, 3807, 5110, 5131, 5153 | 1 | | | 16 |
| *3903, 3924, 3937, 4807 | | | | 14 16 |
| *3909 | | 7 | | 16 |
| *4030, 4034, 4817 | 3 | | | 16 |
| *4059 | 3 | | 10 | 16 |
| *4060 | | | w | 16 |
| *4406, 5105, 5108, 5120, 5144N | 1 | | 9 | |
| *4003 | | | 10 | 14 16 |
| *4020, 4038, 4052, 4606 | | | 10 | 16 |
| *4028 | 2 | 8 | | 16 |
| *4415 | 4 | 6 | 9 12 | |
| *4418, 4501, 4504-4507, 5613 | 4 | | 9 12 | |
| *4503 | | | 9 12 | |
| *4814 | | | 10 11 | 14 16 |
| *4901, 4902, 5001, 5002, 5004 | 4 | 8 | 12 | 16 |
| *4903 | 1 4 | 8 | 12 | 15 16 |
| *510101-5103, 5106, 5111N, 5112, 5114, 5116-5118, 5121-512402, 5126-5128, 5130, 5132-5135, 5138, 5139, 5141N, 5143, 5148, 5149, 5151, 5152, 5157, 520601-520602 | 1 | 8 | | 16 |
| *5104, 5146, 5156 | 1 2 | 8 | | 16 |
| *5107, 520101-520104, 5204, 5205, 5207, 5208, 5211-5213 | 1 | 8 | | 15 16 |
| *510901-510902, 5119, 5140, 5142 | 1 | | 9 | 16 |

