

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	
HLA-B allele ^{1,2,3}																													
*57:01:01-57:01:04, 57:01:06-57:01:18, 57:38, 57:48, 57:50, 57:54, 57:56, 57:59-57:60, 57:64-57:65, 57:67-57:68	1	2	3										13																
*57:01:05, 57:47		2	3										13																
*57:02:01-57:02:02, 57:42	1			4	5							12	13																
*57:03:01-57:03:02, 57:46, 57:66, 57:70	1		3	4								12	13																
*57:04	1				5	6							13					19											
*57:05	1			4									13																
*57:06	1	2	3			6							13																
*57:07	1		3	4				8				12	13																
*57:08	1	2	3						10				13																
*57:09	1			4					9		11	12	13																
*57:10	1	2	3										13	14															
*57:11													13		15														
*57:12	1			4	5							12	13			16													
*57:13		2			5				9				13				17	19											
*57:14	1	2	3								11		13																
*57:15, 57:37, 57:55 ⁵	1	2	3				7						13																
*57:16	1	2	3					8																					
*57:17	1		3	4								12	13	14															
*57:18	1	2	3										13			16													
*57:19	1	2			5								13																
*57:20	1	2	3																20										
*57:21	1	2	3										13							21									
*57:22	1	2	3						9				13				17												
*57:23	1	2	3										13				18												
*57:24	1	2									11		13																
*57:25, 57:41 ⁵	1	2	3										13					19											
*57:26	1	2	3										13						20										
*57:27	1	2	3										13				17												
*57:28N	1				5								13																
*57:29	1	2	3										13	15															
*57:30	1	2			5	6							13																
*57:31		2	3										13				17												
*57:32	1		3			6							13																
*57:33	1	2	3										13	15						21									
*57:34	1	2	3										13								22								
*57:35	1	2	3										13											23					
*57:36	1	2	3																				23						
*57:39	1		3	4								12	13				18												
*57:40	1	w	3										13								22		24						
*57:43	1	2	3								11		13					19											
*57:44	1	2	3										13													25			
*57:45, 57:62, 57:69, 35:208, 58:14 ⁴	1	2	3																										
*57:49	1	2	3					8					13																
*57:51	1	2	3														18												
*57:52	1	2											13																
*57:53	1		3										13																
*57:57	1		3	4					9			12	13				17												
*57:58	1	2	3										13														26		
*57:61	1	2											13															27	
*57:63	1			4	5				9			12	13				17												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	
*07:02:01-07:24, 07:26-07:47, 07:49N-07:50, 07:52-07:119, 07:121-07:136, 07:138-07:154, 07:156-07:163, 07:165-07:201N, 07:203-07:208, 07:210-07:218, 08:20, 08:53:01-08:53:02, 08:79, 13:16, 13:31, 13:48, 13:62, 14:02:09, 14:05-14:06:02, 14:08:01-14:08:02, 14:13, 14:37, 15:01:01:01-15:01:04, 15:01:06-15:01:35, 15:03:01-15:04:02, 15:06-15:12, 15:14-15:19, 15:23-15:24:02, 15:26N-15:30, 15:32:01-15:35, 15:37-15:40, 15:42-15:43, 15:45-15:47:02, 15:49-15:50, 15:53-15:54, 15:56-15:58, 15:60-15:61, 15:63-15:74, 15:76, 15:78:01-15:79N, 15:81-15:82, 15:87, 15:90, 15:92-15:99, 15:101-15:104, 15:108-15:110, 15:113, 15:115-15:120, 15:122, 15:125-15:129, 15:131-15:135, 15:137-15:138, 15:140-15:143, 15:146-15:150, 15:152-15:153, 15:156-15:161, 15:163-15:164, 15:166-15:169, 15:171-15:175, 15:177-15:178, 15:180-15:184, 15:187, 15:189-15:193, 15:195-15:199, 15:201-15:203, 15:205-15:212, 15:215-15:217, 15:219-15:221, 15:225-15:234, 15:236, 15:238-15:239, 15:241-15:249, 15:251-15:264																													

Negative Control

