

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
HLA-B allele ^{1,2,3}																																	
*57:01:01:01-57:01:04, 57:01:06-57:01:44, 57:38, 57:50, 57:54, 57:56, 57:59-57:60, 57:64-57:65, 57:67:02-57:68, 57:71-57:72, 57:75, 57:81, 57:86-57:88, 57:90-57:92, 57:95, 57:97, 57:100, 57:102-57:108, 57:111, 57:113-57:117, 57:121-57:123, 57:125-57:127, 57:130N, 57:132-57:133, 57:135-57:136, 57:138-57:139N, 57:141-57:144	1	2	3										13																			31	
*57:01:05, 57:47, 57:74		2	3										13																			31	
*57:02:01:01-57:02:02, 57:42	1			4	5							12	13																			31	
*57:03:01:01-57:03:05, 57:46, 57:66, 57:70, 57:84, 57:96, 57:101, 57:124, 57:137	1		3	4									12	13																		31	
*57:04:01, 57:04:03	1				5	6							13						19													31	
*57:04:02	1				5	6							13						19													31	
*57:05	1			4									13									23										31	
*57:06	1	2	3			6							13																			31	
*57:07	1		3	4				8					12	13																		31	
*57:08	1	2	3							10			13																			31	
*57:09	1			4					9	11	12		13																			31	
*57:10	1	2	3										13	14																		31	
*57:11													13		15																	31	
*57:12	1			4	5								12	13		16																31	
*57:13		2			5				9				13				17	19														31	
*57:14:01-57:14:02	1	2	3							11			13																			31	
*57:15, 57:37, 57:55 ⁴	1	2	3				7						13																			31	
*57:16	1	2	3					8																								31	
*57:17	1		3	4									12	13	14																	31	
*57:18	1	2	3										13				16															31	
*57:19	1	2			5								13																			31	
*57:20	1	2	3										13								20											31	
*57:21	1	2	3										13									21										31	
*57:22	1	2	3						9				13					17														31	
*57:23	1	2	3										13						18													31	
*57:24	1	2								11			13																			31	
*57:25, 57:41 ⁴	1	2	3										13						19													31	
*57:26	1	2	3										13								20											31	
*57:27	1	2	3										13				17															31	
*57:28N	1				5								13																			31	
*57:29	1	2	3										13	15																		31	
*57:30	1	2			5	6							13																			31	
*57:31:01-57:31:02	1	2	3										13				17															31	
*57:32	1		3				6						13																			31	
*57:33	1	2	3										13	15							21											31	
*57:34	1	2	3										13									22										31	
*57:35	1	2	3										13										23									31	
*57:36	1	2	3										13										23									31	
*57:39	1		3	4									12	13															28		31		
*57:40	1	w	3										13									22	24								31		
*57:43	1	2	3								11		13					19														31	
*57:44	1	2	3										13																25			31	
*57:45, 57:62, 57:69, 57:89, 57:109, 57:120, 57:128	1	2	3										13																			31	
*57:48, 57:85, 57:110	1	2	3										13																			31	
*57:49, 57:93	1	2	3					8					13																			31	
*57:51	1	2	3										13																28			31	
*57:52, 57:140	1	2											13																			31	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*57:53, 57:76	1		3										13																			31	
*57:57	1		3	4					9				12	13					17													31	
*57:58	1	2	3										13														26					31	
*57:61	1	2											13															27				31	
*57:63	1			4	5				9				12	13					17													31	
*57:67:01	1	2	3										13										23									31	
*57:73	1	2	3										13																28			31	
*57:77	1	2	3										13																			31	
*57:78	1	2	3										13										24									31	
*57:79N	1	2	3										13																			31	
*57:80	1		3	4		</																											

<p>13:87, 13:90, 13:92, 13:95, 13:106-13:107, 13:109-13:110, 13:122, 13:129, 13:135, 13:138, 13:141-13:142, 13:144, 13:146, 13:148, 13:161N-13:163, 15:88, 15:250, 15:283, 18:22, 18:69, 18:105, 18:149-18:150, 18:204, 27:19:01:01-27:19:01:02, 27:30, 27:127, 27:150, 35:01:01:01-35:01:54, 35:01:56-35:01:61, 35:01:63-35:01:67, 35:03:01:01-35:03:33, 35:07-35:08:14, 35:10, 35:13, 35:15:01-35:17:04, 35:19-35:20:02, 35:23-35:30:01:02, 35:33-35:36, 35:38-35:42:02, 35:45, 35:47-35:50, 35:52, 35:54-35:57, 35:59:02, 35:61:01-35:63, 35:65Q, 35:69-35:71, 35:74, 35:76-35:78, 35:80, 35:82, 35:84-35:85, 35:90-35:94, 35:98, 35:100-35:101:02, 35:103-35:113, 35:115-35:116, 35:120-35:126, 35:128, 35:130N-35:134N, 35:136-35:148, 35:150:01-35:150:02, 35:152-35:153, 35:155-35:156, 35:158-35:161, 35:163-35:166, 35:168-35:171, 35:173:01N-35:173:02N, 35:175-35:181, 35:186-35:187, 35:189-35:198, 35:200, 35:202-35:204, 35:206, 35:209-35:210, 35:212, 35:214-35:219, 35:221-35:225, 35:227-35:229, 35:231, 35:236-35:249, 35:252-35:257, 35:259, 35:261-35:264, 35:267-35:269, 35:271-35:272, 35:274-35:275, 35:278-35:284, 35:286-35:291, 35:293-35:299, 35:301-35:308, 35:310, 35:312-35:315, 35:318-35:320, 35:322, 35:324-35:330, 35:332-35:334, 35:336, 35:338, 35:340-35:342, 35:344-35:356, 35:358-35:360, 35:362-35:365, 35:370-35:371, 35:373, 35:375-35:376, 35:380:01:01-35:383, 35:385-35:386, 35:388, 35:392-35:395, 35:397-35:402, 35:404-35:408, 35:411-35:412, 35:414-35:415, 35:417-35:418, 35:420-35:424, 35:426, 35:428Q-35:432, 35:434-35:438, 35:441-35:442, 35:444-35:458, 35:460Q-35:462, 35:464-35:467, 35:469-35:470, 35:473, 35:478-35:480, 35:482, 35:484-35:486, 35:488-35:491, 35:493-35:497, 35:500, 35:502-35:504, 35:506-35:507, 35:510-35:511, 35:513N-35:518, 35:520, 35:524-35:529, 37:25, 37:35, 38:20:02, 40:458, 44:02:01:01-44:02:33, 44:02:35-44:02:57, 44:02:59-44:14, 44:16-44:17, 44:19N, 44:21-44:30, 44:32-44:35, 44:37:01-44:40, 44:42-44:46:02, 44:48-44:52N, 44:55-44:61N, 44:63-44:64:02, 44:66-44:75, 44:78, 44:80-44:81, 44:83-44:98, 44:101-44:105, 44:107-44:134, 44:136-44:137, 44:139-44:145, 44:147-44:149N, 44:151-44:152, 44:154-44:157, 44:159-44:165, 44:167-44:181, 44:183, 44:185-44:188, 44:190-44:191, 44:194:01-44:196, 44:198N-44:212, 44:214-44:229, 44:231-44:235, 44:237N-44:245, 44:247-44:248, 44:250-44:304, 44:306N, 44:308-44:310N, 44:312-44:316, 44:318-44:338, 44:340-44:342, 44:344-44:346, 44:437-44:439, 44:441-44:449N, 44:451, 44:453-44:458, 44:460-44:490, 44:492-44:505, 44:507-44:516, 44:519-44:524, 48:02:01-48:02:03, 51:42, 53:01:01:01-53:05:02, 53:07, 53:09-53:13, 53:15-53:18, 53:20-53:27, 53:29-53:35, 53:37-53:40, 53:42-53:43:02, 53:45-53:46, 53:48N-53:66, 56:09, 56:11, 56:60, 58:01:01:01-58:01:02, 58:01:04-58:01:15, 58:01:17-58:01:39, 58:04-58:05, 58:09-58:11, 58:13, 58:15, 58:17N, 58:21-58:24, 58:29, 58:31N-58:35, 58:37, 58:39N-58:40, 58:42, 58:44-58:45:02, 58:47-58:53, 58:55-58:59:01, 58:62-58:63, 58:66-58:72N, 58:74-58:77, 58:79-58:82, 58:84-58:90, 58:92-58:103, 58:105-58:113, 58:115-58:120, 58:123-58:126, 58:128N-58:131, 58:133N, 83:01, C*03:102, C*03:263:01-03:263:03, C*03:514, C*15:126</p>	<p>1</p>								
<p>*08:60, 08:129, 13:21, 35:09:01-35:09:03, 35:18, 35:88, 35:230, 35:321, 35:372, 37:12, 40:04:01:01-40:04:02, 40:64:01:01-40:64:02, 40:68, 40:129, 40:137, 40:160:01-40:160:02, 40:343, 41:48, 42:13, C*03:317, C*03:382, C*15:24</p>	<p>4</p>								
<p>*08:76, 37:01:08, 37:01:10, 37:81, 40:99, 48:17, C*05:10</p>	<p>w</p>			<p>16</p>					
<p>*08:119, 39:92</p>				<p>17</p>	<p>20</p>				
<p>*08:181</p>	<p>4</p>		<p>11 12</p>				<p>24</p>		
<p>*08:262, 40:284, A*03:174, A*11:65, A*24:247, C*01:170, C*04:83, C*07:511, C*15:121, C*16:162</p>									
<p>*13:01:18, 13:26:01, 15:20, 15:393, 15:407, 15:587, 15:599, 35:46, 35:207, 46:33</p>	<p>1</p>					<p>23</p>			
<p>*13:02:06, 13:72, 13:127, C*01:73</p>						<p>23</p>			
<p>*13:13:01-13:13:02, 13:86, 35:02:01:01-35:02:19, 35:04:01:01-35:04:03, 35:12:01:01-35:12:04, 35:81, 35:83, 35:95, 35:129N, 35:149, 35:154, 35:157, 35:162, 35:172, 35:182-35:184, 35:201, 35:211, 35:220, 35:233, 35:258, 35:266, 35:270, 35:285, 35:309, 35:311, 35:316, 35:323, 35:335, 35:339, 35:357, 35:361, 35:366, 35:374, 35:377-35:379, 35:387, 35:391, 35:396, 35:403, 35:410, 35:419, 35:425, 35:443, 35:463, 35:468, 35:474, 35:477, 35:483, 35:487, 35:492, 35:501, 35:505, 35:508N, 35:522, 37:19:01-37:19:02, 37:84, 44:62, 44:77, 51:139, 53:19, 53:36, 56:12, 58:28:02, 58:65</p>	<p>1</p>	<p>4</p>							
<p>*13:20, 14:10, 14:102, 15:02:01:01-15:02:12, 15:13:01-15:13:03, 15:21:01:01-15:21:01:02, 15:25:01:01-15:25:04, 15:36, 15:44, 15:62, 15:77, 15:80, 15:89:01-15:89:02, 15:106, 15:112, 15:121, 15:139, 15:144, 15:154, 15:165, 15:170, 15:204, 15:213, 15:223, 15:240, 15:265, 15:271, 15:289, 15:291, 15:297, 15:301-15:302N, 15:308, 15:319, 15:325, 15:328, 15:330, 15:341, 15:345, 15:357-15:358, 15:366, 15:374, 15:378, 15:394, 15:399, 15:402, 15:404, 15:418, 15:420, 15:425, 15:437-15:438, 15:454N, 15:458, 15:460, 15:463N, 15:468, 15:491, 15:495, 15:505, 15:537, 15:543-15:544N, 15:557, 15:564, 15:566-15:567, 15:569, 15:574, 15:584N, 15:589, 15:592, 15:594, 15:596N, 35:11:01-35:11:03, 35:14:01-35:14:03, 35:21, 35:96, 35:167, 35:188, 35:213, 35:277, 35:337:01-35:337:02, 35:343, 35:390N, 35:413, 35:439:01:01-35:439:01:02, 35:512, 44:76, 44:79:01-44:79:02, 44:146, 44:150, 44:452, 46:67, 48:25, 49:61, 50:17, 51:227, 51:281, 52:88, 53:06, 53:08:01-53:08:02, 53:28, 58:19, 58:78, 58:91</p>	<p>1</p>	<p>9</p>							
<p>*14:01:01:01-14:02:08, 14:02:10-14:04, 14:07N, 14:09, 14:11-14:12, 14:14-14:19, 14:21-14:23, 14:25-14:36, 14:38-14:40, 14:42, 14:44-14:52, 14:54, 14:57-14:58, 14:60-14:70Q, 14:72N-14:89, 14:91-14:95, 14:97-14:101N, 14:103, 15:510, 35:475, 46:84, 50:51, 51:165, 55:108, 58:06</p>			<p>9</p>	<p>15</p>					
<p>*14:20</p>			<p>9</p>	<p>15</p>	<p>22</p>				
<p>*14:24</p>			<p>9</p>	<p>15 16</p>					
<p>*14:41N, 14:43, 14:59, 14:90, 18:44:01-18:44:02, 18:81, 39:79, 39:145, 39:158, 40:427, 44:440, 54:38, 58:02:01:01-58:02:02, 58:07, 58:25, 58:38, 58:43, 58:60, 58:104, 58:132, C*01:32:01-01:32:02, C*06:20, C*14:82, C*16:98, C*16:102, C*16:110</p>				<p>15</p>					
<p>*14:104</p>		<p>7</p>	<p>9</p>	<p>15</p>					

Negative Control

