

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-B allele ^{1,2}																
*57:01:01-57:01:15	1	2									11					
*57:02:01-57:02:02, 57:12, 57:42	1					6										
*57:03:01-57:03:02, 57:17, 57:39	1	2				6										
*57:04	1		3							10						
*57:05, 57:11, 57:28N, 58:36	1															
*57:06, 57:18, 57:27	1	2	3								11					
*57:07, 57:46	1	2			5	6										
*57:08	1	2				6					11	12				
*57:09	1					6	7									
*57:10, 57:44	1	2						8			11					
*57:13	1						7			10	11					
*57:14, 57:25, 57:31, 57:50	1	2					7				11					
*57:15, 57:29	1	2		4							11					
*57:16		2			5						11					
*57:19, 57:30, 57:52, 57:61	1										11					
*57:20, 57:62		2		4							11					
*57:21, 57:33	1	2							9		11					
*57:22, 57:37, 57:41, 57:43	1	2								10	11					
*57:23, 57:26, 57:34, 57:64	1	2			5						11					
*57:24	1						7				11					
*57:32	1	2	3													
*57:35, 57:38, 57:60	1	2									11	12				
*57:36		2									11	12				
*57:40	1	2							9		w					
*57:45, 57:51		2									11		13			
*57:47, 57:54, 57:59	1	2									11			14		
*57:48, 57:55-57:56	1	2									11				15	
*57:49, 57:58	1	2									11		13			
*57:53	1	2					7									
*57:57	1	2				6				10						
*57:63	1					6				10						
*57:65	1	2				6					11					
*07:120, 15:214, 18:81, 40:150				4												
*14:20, 35:127									9							
*15:33, 15:116, 15:248, 40:63, 40:92, 44:169, 49:22, 51:126, C*03:87, C*03:129, C*05:27, C*05:39, C*17:07														14		
*15:87												12				
*35:208, 58:14		2					7				11					
*40:30, 40:34		2				6	7									
*44:153		2	3		5		7			10						
*55:14		2			5		7			10	11					
*58:12, A*02:42, A*02:285, A*02:310													13			
A*33:12															15	
C*06:72					5		7			10	11					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control