

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-B allele^{1,2,3}																								
*58:01:01-58:01:14, 58:32-58:35, 58:37, 58:40-58:42, 58:45, 58:47-58:48	1	2											13											
*58:02, 58:38, 58:43	1		3								11		13											
*58:04		2		4									13											
*58:05	1	2			5								13											
*58:06, 58:46	1		3			6	7						13											
*58:07	1		3					8			11		13											
*58:08:01	1					6	7		9				13											
*58:08:02	1					6	7						13											
*58:09, 58:24 ⁴	1	2								10			13											
*58:10N, 58:13 ⁴	1	2										12	13											
*58:11	1	2												14										
*58:12	1	2											13		15									
*58:14	1	2											13								21			
*58:15	1	2		4									13											
*58:16:01-58:16:02	1												13				17							
*58:17N	1	2											13					18						
*58:18	1												13				17	19						
*58:19	1	2				6	7						13											
*58:20	1													14										
*58:21, 58:29 ⁴	1	2											13								20			
*58:22	1	2											13				17							
*58:23	1	2							9				13											
*58:25	1		3								11		13					18						
*58:26	1												13											23
*58:27	1												13			16								
*58:28	1	2											13		16									
*58:31N	1	2											13											23
*58:36		2											13		15									
*58:39N		2											13								20			
*58:44		2							9				13											
*07:02:01-07:02:34, 07:02:36-07:04, 07:07-07:08, 07:10, 07:13, 07:15-07:16, 07:19-07:24, 07:26-07:31, 07:33:01-07:33:02, 07:35-07:39, 07:41-07:47, 07:49N-07:50, 07:52, 07:54-07:57, 07:59-07:68:01, 07:68:03, 07:70-07:79, 07:81-07:89, 07:91-07:96:02, 07:98-07:104, 07:106-07:111N, 07:113-07:122, 07:124-07:136, 07:139, 07:141-07:154, 07:156-07:161N, 07:165-07:175, 07:177-07:181N, 07:183-07:198, 07:200, 07:203-07:204, 13:31, 14:02:07, 15:01:01:01-15:01:04, 15:01:06-15:01:33, 15:03:01, 15:03:03, 15:04, 15:06, 15:08-15:12, 15:14-15:19, 15:23-15:24:02, 15:26N-15:30, 15:32-15:35, 15:37-15:39:01, 15:40, 15:43, 15:45-15:47:02, 15:49-15:50, 15:53-15:54, 15:56-15:58, 15:60-15:61, 15:63-15:67, 15:69-15:74, 15:76, 15:78:01-15:79N, 15:81-15:82, 15:87, 15:90, 15:92-15:99, 15:101, 15:103-15:104, 15:108-15:110, 15:113, 15:115-15:120, 15:122, 15:125, 15:127-15:129, 15:131-15:135, 15:137-15:138, 15:140-15:143, 15:146-15:150, 15:152-15:153, 15:156-15:161, 15:163-15:164, 15:166-15:169, 15:171-15:175, 15:177-15:178, 15:180-15:184, 15:187, 15:189-15:193, 15:195-15:199, 15:201-15:203, 15:205-15:206, 15:208-15:212, 15:215-15:217, 15:219-15:221, 15:225-15:227, 15:229-15:234, 15:236, 15:238-15:239, 15:241-15:249, 15:251-15:264, 15:266-15:270, 15:272N-15:274, 15:276-15:282, 15:285-15:288, 15:290, 18:19, 18:21, 18:30, 18:57, 27:04:01-27:04:03, 27:06, 27:10, 27:15, 27:18, 27:20-27:21, 27:25, 27:40, 27:54, 27:63, 27:66N, 27:68-27:69, 27:79, 27:86, 27:91-27:92, 27:100, 27:103, 27:105-27:107, 35:43:01, 35:79, 35:185, 37:07, 39:18, 39:36, 40:23, 40:26, 40:32, 40:51, 40:98, 40:148, 40:158, 40:161, 40:183, 40:198, 46:01:01-46:05, 46:07N-46:08, 46:10, 46:14-46:17, 46:20, 46:22-46:24, 46:26-46:32, 46:34-46:42, 48:08, 49:01:02, 49:22, 50:01:03, 54:06, 55:21, 55:56, 56:03						6							14											
*07:02:35, 07:68:02, 14:03, 14:12, 50:32, 54:20, 55:09, 55:52						6	w						14											
*07:05:01-07:06, 07:32, 07:34, 07:40, 07:53, 07:69, 07:80, 07:90, 07:97, 07:105, 07:112, 07:123, 07:138, 07:140, 07:176, 07:182N, 07:201N, 08:20, 08:53:01-08:53:02, 08:79, 27:24, 40:05, 40:15-40:16, 40:174, 48:05, 48:15						6							14		16									
*07:09, 07:11, 07:17, 15:07:01-15:07:02, 15:68, 15:126, 15:207, 46:12						6							14					19						
*07:12, 07:14, 07:18:01-07:18:02, 07:162-07:163, 07:199						6							14				17							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control

