

INTERPRETATION TABLE					
HLA-B*67 SSP subtyping					
Amplification patterns of the B*67 alleles					
	Well				
	1	2	3	4	
Length of spec.	210	325	165	160	Length of spec.
PCR product					PCR product
Length of int.	800	1070	800	1070	Length of int.
pos. control ¹					pos. control ¹
5'-primer ²	103	44	419	165	5'-primer ²
	5'-CCT ^{3'}	5'-ggC ^{3'}	5'-gTT ^{3'}	5'-ACC ^{3'}	
3'-primer ³	272	201	544	282	3'-primer ³
	5'-TgT ^{3'}	5'-CTT ^{3'}	5'-ggT ^{3'}	5'-gCC ^{3'}	
Well No.	1	2	3	4	Well No.
HLA-B allele					HLA-B allele
*670101-670102	1		3		*670101-670102
*6702	2		3		*6702
*070201-070202, 070204-070207, 0704-0707, 0709, 0711, 0712, 0714, 0715, 0717-0726, 0728, 0730, 0731, 0733-0736, 0739-0746, 0748, 0749N, 0751-0755, 0757-0764, 0766-0768, 4201, 4204-4206, 4208, 8101-8104N	1			4	*070201-070202, 070204-070207, 0704-0707, 0709, 0711, 0712, 0714, 0715, 0717-0726, 0728, 0730, 0731, 0733-0736, 0739-0746, 0748, 0749N, 0751-0755, 0757-0764, 0766-0768, 4201, 4204-4206, 4208, 8101-8104N
*070203, 0747, 0756, 5510	1				*070203, 0747, 0756, 5510
*0710, 0713, 0729, 0765, 1576, 9501, 3576, 4202, 4209, 460101, 4603-4618, 5603, 5606, 8301				4	*0710, 0713, 0729, 0765, 1576, 9501, 3576, 4202, 4209, 460101, 4603-4618, 5603, 5606, 8301
3817, 39010101-39010102L, 390103-390602, 3908-391302, 3915, 3916, 3918-3920, 3922-3924, 3926-3928, 3931, 3935, 3937-3942, 3944, 3945, 5820	3				3817, 39010101-39010102L, 390103-390602, 3908-391302, 3915, 3916, 3918-3920, 3922-3924, 3926-3928, 3931, 3935, 3937-3942, 3944, 3945, 5820
*460102-4602, 7301				w	*460102-4602, 7301
HLA-B allele					HLA-B allele
Well No.	1	2	3	4	Well No.