

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-B allele ^{1,2,3}																
*78:01:01:01-78:01:02	1	2					7			10		12				
*78:02:01	1		3							10		12				
*78:02:02	1		3				7			10		12				
*78:03	1	2		4			7			10		12				
*78:04			3		5					10		12				
*78:05	1		3			6				10		12				
*78:06	1		3			6				10		12	13			
*78:07	1	2					7					12				15
*78:08	1	2			5		7					12				
*78:09	1		3				7			10		12	13			
*07:02:01:01-07:02:52, 07:04-07:07, 07:09-07:15, 07:17-07:26, 07:28-07:31, 07:33:01-07:36, 07:39-07:46, 07:48-07:49N, 07:51-07:64, 07:66-07:68:03, 07:70-07:84, 07:86-07:133, 07:135N-07:179, 07:181N-07:205, 07:208-07:218, 07:220-07:227, 07:229-07:231N, 07:233-07:298, 08:123, 08:132, 14:21, 14:57, 15:76, 15:101, 15:255, 27:01-27:11, 27:13-27:15, 27:17, 27:19-27:21, 27:24-27:28, 27:30-27:38, 27:40-27:51, 27:53-27:58, 27:60-27:76, 27:78-27:91, 27:93-27:101, 27:103-27:118, 27:120-27:128, 27:130-27:136, 27:138-27:152, 27:154-27:164, 38:26, 40:166, 42:01:01-42:02:01:02, 42:05:01-42:06, 42:08-42:10, 42:12-42:23, 44:90, 44:97, 45:06, 46:01:01-46:10, 46:12-46:17, 46:19-46:31, 46:33-46:39, 46:42-46:48, 46:51Q-46:72, 54:06, 54:09, 54:26, 54:38, 55:01:07, 55:04, 55:08-55:09, 55:14, 55:22, 55:24, 55:49, 55:51, 55:86, 56:02:01-56:04, 56:09-56:12, 56:18-56:19N, 56:31-56:32, 56:51, 67:01:01-67:07, 81:01-81:08, 82:01-82:03, 83:01									9							
*07:65		2					7		9			12				
*07:134, 35:76, 35:186, 35:229, 46:40:02									9			12				
*07:206		2					7		9							
*08:03, 08:52, 08:78, 15:361, 40:13, 40:19, 40:109, 40:117, 40:292, 40:340, 41:46, 48:18, 49:18:01-49:18:02, 57:02:01-57:03:03, 57:05, 57:07, 57:09, 57:17, 57:39, 57:42, 57:46, 57:57, 57:63, 57:66, 57:70, 57:80, 57:82, 57:84, A*24:146														14		
*08:09, 08:84, 13:46, 13:53, 39:06:01-39:06:05, 39:34, 39:57, 39:62, 39:64, 39:83, 39:90, 39:109, 39:127, 39:129, 55:78, 59:01:01:01-59:01:01:02, 59:03, 59:05-59:08, 59:10N, C*05:152, C*07:463, C*12:183															15	ontrol

*51:09:01, 51:09:03, 51:19, 51:73, 51:122	1				5		8	10	11		14					
*51:09:02	1				5			10	11		14					
*51:15								10	11							
*51:20, 51:208	1	w					8	10	11		14					
*51:21, 51:101, 51:153	1						8		11		14	15				
*51:22	1		4				8	10	11		14					
*51:31	1				5		8		11		14					
*51:36, 51:143							8		11		14	15				
*51:40, 51:108					5		8	10	11		14					
*51:54					w			10	11							
*51:93					5				11		14					
*51:104, 51:112, 51:147, 52:02:01	1						8	10			14					
*51:139, 51:194					5		8		11		14					
*51:144, 51:148, 53:19, 58:27-58:28:02, 58:65									11		14					
*51:157, 51:197							8		11			15				
*51:185	1	3					8	10	11		14					
*51:217	1				w		8	10	11		14					
*52:01:01:01-52:01:01:04, 52:01:03, 52:01:05-52:01:08, 52:01:10-52:01:18, 52:01:20, 52:01:22-52:01:28, 52:04-52:05, 52:06:03-52:07, 52:13, 52:15, 52:17, 52:19, 52:21:02-52:24, 52:27-52:31:01, 52:33-52:36, 52:38, 52:40-52:61, 52:64-52:66, 52:69-52:70	1					6		10	11		14					
*52:01:02:01-52:01:02:03, 52:01:04, 52:01:09, 52:06:01-52:06:02, 52:21:01, 52:31:02	1					6	8	10	11		14					
*52:01:19	1					6		10			14					
*52:01:21, 52:08	1					6			11		14					
*52:01:29, 52:68	1	3				6		10	11							
*52:02:02, 52:16, 52:32	1							10			14					
*52:03					5	6	8	10	11		14					
*52:09	1	3				6	8	10	11		14					
*52:10:01-52:10:02	1				5	6		10	11		14					
*52:10:03	1				5	6	8	10	11		14					
*52:11						6		10	11		14					
*52:12, 52:14, 52:20, 52:25:01-52:26, 52:39, 52:62-52:63	1					6		10	11							
*52:37	1					6			11		14	15				
*52:67						6			11		14					
*52:71	1				5	6			11		14	15				
*55:31			3					9				15				
*56:05:01	1	2					7	9	10							
*56:05:02, 56:34	1							9	10							
*56:06	1	2					7	9	10	12						
*56:21	1							9	10		14					
*56:45								9	10	12						
*57:01:20							8									
*59:02											14	15				
*73:01-73:02	1							9								
A*23:01:14		w					8									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control