

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-B allele^{1,2,3}																
*78:01:01:02-78:01:02	1	2					7			10		12				
*78:02:01	1		3							10		12				
*78:02:02	1		3				7			10		12				
*78:03	1	2		4			7			10		12				
*78:04			3		5					10		12				
*78:05	1		3			6				10		12				
*78:06	1		3			6				10		12	13			
*78:07	1	2					7					12			15	
*78:08	1	2			5		7					12				
*78:09	1		3				7			10		12	13			
*78:10	1		3		5		7			10		12	13			
*07:02:01:01-07:02:86, 07:04:01-07:07:03, 07:09:01-07:15, 07:17-07:26, 07:28-07:31:02, 07:33:01-07:36, 07:39-07:46, 07:48-07:49N, 07:51-07:64, 07:66-07:68:03, 07:70-07:84, 07:86-07:133, 07:135N-07:179, 07:181N-07:205, 07:208-07:218, 07:220-07:227, 07:229:01-07:231N, 07:233-07:305, 07:307-07:399, 07:401-07:435N, 08:123, 08:132, 14:21, 14:57, 15:76, 15:101, 15:255, 15:594, 27:01-27:11, 27:13:01-27:15, 27:17, 27:19:01:01-27:21:02, 27:24-27:28, 27:30-27:38, 27:40-27:51, 27:53-27:58, 27:60-27:76, 27:78-27:91, 27:93-27:101, 27:103-27:118, 27:120-27:128, 27:130-27:136, 27:138-27:152, 27:154-27:203, 27:205-27:214, 27:216-27:235, 27:237-27:241, 27:243N-27:248, 38:26, 40:166, 42:01:01:01-42:02:02, 42:05:01-42:06, 42:08-42:10, 42:12-42:23, 42:25, 42:27-42:31, 44:90, 44:97, 45:06, 46:01:01:01-46:10, 46:12-46:17, 46:19-46:31, 46:33-46:39, 46:42-46:48, 46:51Q-46:91, 54:06, 54:09, 54:26, 54:38, 55:01:07, 55:04, 55:08-55:09, 55:14, 55:22, 55:24, 55:49, 55:51, 55:86, 55:110, 56:01:13, 56:02:01:01-56:04:04, 56:09-56:12, 56:18-56:19N, 56:31-56:32, 56:51, 56:56-56:57, 56:60, 56:63, 56:66, 67:01:01-67:07, 81:01:01:01-81:10, 82:01:01:01-82:04, 83:01									9							
*07:65		2					7		9			12				
*07:134, 35:76, 35:186, 35:229, 46:40:02									9			12				
*07:206		2					7		9							
*07:400, 18:12:01-18:12:02, 35:69, 37:11, 37:37, 40:209, 40:350						6						12				
*08:03, 08:52, 08:78, 15:361, 38:73, 40:13, 40:19, 40:109, 40:117, 40:292, 40:340, 41:46, 48:18, 49:18:01-49:18:02, 57:02:01:01-57:03:05, 57:05, 57:07, 57:09, 57:17, 57:39, 57:42, 57:46, 57:57, 57:63, 57:66, 57:70, 57:80, 57:82, 57:84, 57:94, 57:96, 57:101, 57:124, 57:137, A*24:146														14		
*08:09, 08:84, 08:254, 13:46, 13:53, 13:131, 39:06:01:01-39:06:07, 39:34, 39:57, 39:62, 39:64, 39:83:01-39:83:03, 39:90, 39:109, 39:127, 39:129, 39:132, 39:140-39:141, 39:154, 39:178, 41:70, 55:78, 59:01:01:01-59:01:01:02, 59:03, 59:05-59:08, 59:10N, C*05:152, C*07:463, C*07:758, C*12:183, C*14:92															15	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*08:23, 08:54:01-08:54:02, 08:110, 13:09, 13:39, 13:67, 14:18, 14:73, 15:01:01:01-15:01:04, 15:01:06-15:01:29, 15:01:31-15:01:54, 15:01:56-15:01:64, 15:01:66-15:03:13, 15:05:01:01-15:12:03, 15:14-15:15:01:02, 15:18:01:01-15:21:01:02, 15:25:01:01-15:32:01, 15:33-15:35, 15:37-15:40:02, 15:44-15:56, 15:58, 15:60, 15:62-15:66, 15:69, 15:72-15:75, 15:77-15:82, 15:84-15:88, 15:90-15:93, 15:96-15:99, 15:102-15:114, 15:116-15:125:01, 15:126-15:129, 15:131-15:136, 15:138-15:147, 15:149N-15:154, 15:156, 15:158-15:161, 15:163, 15:165-15:167, 15:169-15:174, 15:176, 15:178-15:179:02, 15:181N-15:188, 15:190N-15:195, 15:197:01-15:201, 15:203-15:207, 15:209N-15:215, 15:217-15:221, 15:223-15:229, 15:231-15:235, 15:237-15:238, 15:240-15:253, 15:257-15:266, 15:269-15:272N, 15:274-15:290, 15:292-15:293, 15:295-15:302N, 15:304N-15:307, 15:309, 15:311-15:329, 15:331-15:338, 15:340-15:355, 15:357-15:360, 15:363:01-15:375N, 15:377Q-15:380N, 15:382-15:395, 15:398-15:402, 15:404-15:407, 15:409-15:410, 15:412-15:417, 15:419-15:422, 15:425-15:441, 15:443-15:445, 15:448-15:457, 15:459-15:461, 15:463N-15:484, 15:486-15:498, 15:501-15:504, 15:506-15:508, 15:510-15:515, 15:517-15:522, 15:524-15:530, 15:533-15:545, 15:547-15:549N, 15:551, 15:554-15:555, 15:557-15:562N, 15:564-15:569, 15:572-15:574, 15:576-15:578, 15:580-15:581, 15:583-15:585, 15:587-15:593, 15:595-15:601, 15:604N-15:605, 18:107, 35:01:39, 35:02:09, 35:08:02, 35:19, 35:25-35:26, 35:46-35:47, 35:49:01-35:49:02, 35:63, 35:154, 35:207, 35:210, 35:217, 35:222, 35:248, 35:271, 35:274, 35:293, 35:311, 35:317, 35:361, 35:424, 35:450, 35:522, 35:525, 40:01:01-40:01:05, 40:01:07-40:01:10, 40:01:12, 40:01:14-40:01:39, 40:01:41-40:01:75, 40:07, 40:10:01:01-40:10:02, 40:12:01-40:12:02, 40:14:01-40:16:01:02, 40:21-40:23, 40:25, 40:30-40:34, 40:36, 40:38, 40:42, 40:45-40:46, 40:48-40:49, 40:51-40:55, 40:59-40:63, 40:65, 40:67, 40:69, 40:72:01-40:73, 40:76-40:77, 40:79-40:81, 40:84, 40:87:01-40:88, 40:92, 40:100-40:102, 40:106, 40:108, 40:112-40:114:02, 40:116, 40:118N, 40:121, 40:123-40:126, 40:128-40:130:02, 40:132, 40:134-40:141, 40:146, 40:150-40:156, 40:158, 40:160:01-40:160:02, 40:163, 40:168, 40:170-40:172, 40:175, 40:178-40:179, 40:182-40:183, 40:186:01-40:187, 40:191-40:199, 40:204, 40:207-40:208, 40:212-40:213, 40:215-40:219, 40:221-40:223, 40:227-40:228, 40:231, 40:233-40:242, 40:245:01-40:245:02, 40:247, 40:249-40:253, 40:257-40:265N, 40:272-40:273, 40:277-40:282, 40:285-40:286N, 40:288, 40:298:01-40:301, 40:310, 40:312, 40:315-40:316, 40:319, 40:321, 40:323:01-40:329, 40:332-40:333, 40:336, 40:338N-40:339, 40:342, 40:344, 40:346, 40:352-40:353:02, 40:355, 40:357-40:358, 40:367-40:368, 40:370, 40:372N, 40:377, 40:379-40:383, 40:385-40:386, 40:391-40:392, 40:395-40:398, 40:400-40:404, 40:406, 40:416-40:422, 40:424, 40:426N, 40:428N, 40:430-40:431, 40:433, 40:436, 40:439, 40:443, 40:446-40:450, 40:452, 40:456-40:457, 40:461-40:463, 40:466, 40:469-40:470, 40:474-40:475:02, 40:477, 40:480, 40:483N, 40:485-40:486, 40:490-40:491, 40:495, 41:02:01:01-41:04, 41:08, 41:10-41:11, 41:13, 41:15, 41:17-41:19, 41:23-41:24, 41:31, 41:36, 41:38-41:43, 41:45N, 41:48-41:49, 41:51-41:52, 41:55-41:56Q, 41:58, 41:62-41:63, 41:65, 41:68-41:69, 41:71-41:72, 41:74, 44:02:45, 44:09, 44:46:01-44:46:02, 44:75, 44:129, 44:196, 44:254, 44:300, 44:345N, 44:438N, 44:454-44:455, 45:01:01:01-45:05, 45:07, 45:09-45:21, 45:23-45:25, 45:27, 47:02-47:03, 48:22, 50:01:01:01-50:01:04, 50:01:06-50:02:01:02, 50:04-50:20, 50:31-50:33, 50:35-50:42, 50:44, 50:46-50:48, 50:50-50:51, 50:53-50:57, 50:59-50:70, 50:72N-50:75, 57:01:38, 57:12, 57:38, 57:60, C*06:147, C*12:301			3													
*08:26:03, 08:56:03, 39:63		2					7									
*08:32, 08:180, 08:239, 18:01:01:01-18:03:01:03, 18:05:01:01-18:08, 18:10-18:11, 18:13-18:15, 18:17N-18:24, 18:26-18:36, 18:38-18:47, 18:50-18:53, 18:55-18:65, 18:69-18:78, 18:80-18:98:01, 18:99-18:106Q, 18:109-18:124, 18:126-18:131:01:02, 18:133-18:135:01:02, 18:137-18:141, 18:143-18:177, 18:180-18:193, 18:196-18:204, 18:206-18:212, 35:03:12, 35:03:18, 35:03:21, 35:08:12, 35:29:01-35:29:03, 35:50, 35:84, 35:136, 35:153, 35:162, 35:197, 35:231, 35:243, 35:249, 35:267, 35:269, 35:294, 35:297, 35:305, 35:316, 35:328, 35:393, 35:402-35:403, 35:411, 35:423, 35:472, 35:476, 35:489, 35:494, 35:497, 37:67, 39:19:01-39:19:02												12				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*08:133, 15:308, 18:04:02, 18:25, 18:68, 35:01:01:01-35:01:09, 35:01:11-35:01:38, 35:01:40-35:02:08, 35:02:10-35:03:11, 35:03:13-35:03:17, 35:03:19-35:03:20, 35:03:22-35:03:26, 35:03:28-35:03:31, 35:03:33-35:04:01:02, 35:04:03-35:08:01:10, 35:08:03-35:08:04, 35:08:06-35:08:11, 35:08:13-35:09:03, 35:11:01-35:12:04, 35:14:01-35:15:02, 35:17:01-35:18, 35:20:01-35:24:02, 35:30:01:01-35:36, 35:38-35:42:01, 35:43:01, 35:43:03-35:45, 35:48, 35:51:01-35:55, 35:57-35:59:02, 35:61:01-35:62, 35:64:01-35:68:02, 35:70-35:72, 35:74-35:75, 35:77-35:79, 35:81-35:83, 35:85-35:135, 35:137-35:152, 35:155-35:161, 35:163-35:184, 35:187-35:196, 35:199-35:202, 35:204-35:206, 35:208:01-35:209, 35:211-35:216N, 35:218-35:221, 35:223-35:225, 35:227-35:228, 35:230, 35:232-35:234, 35:236-35:242, 35:244-35:247, 35:250, 35:251:02-35:253, 35:255-35:264, 35:266, 35:268, 35:272, 35:275, 35:277-35:292, 35:295-35:296, 35:298-35:304, 35:306-35:310, 35:312-35:315, 35:319-35:327, 35:330-35:354, 35:356-35:360, 35:362-35:367, 35:369-35:392, 35:394-35:401, 35:404-35:410, 35:412-35:422, 35:425-35:449, 35:451-35:470, 35:473-35:475, 35:478-35:488, 35:490-35:493, 35:495-35:496, 35:498-35:521, 35:523-35:524, 35:526-35:529, 45:22, 58:64			3									12				
*13:03-13:04, 13:48, 13:120, 39:33, 44:20, 44:47, 44:339, 49:54, 56:17, 56:33, 59:04, 59:12, C*03:464										10						
*13:15, 49:09	1															
*13:58, 18:195, 37:01:01:01-37:01:05, 37:01:07-37:04:02, 37:07, 37:09-37:10, 37:12-37:13, 37:15-37:27, 37:29-37:36, 37:38-37:57, 37:59, 37:62-37:66, 37:68-37:81, 37:83-37:95, 37:97, 44:83, 44:134, 44:445, 52:88, 53:17:01-53:17:02, 53:28, 53:38						6					11					
*13:62, 27:153, 44:459	1									10						
*14:37	1														15	
*15:04:01:01-15:04:04, 15:137, 15:155, 15:303, 15:310, 15:499, 15:531, 15:552, 56:42, 56:44, 56:67			3							10						
*15:32:02, 18:125, 35:08:05, 40:01:11, 40:58			3				7									
*15:42, 15:83:01:01-15:83:01:02, 15:485, 41:01:01:01-41:01:07, 41:05-41:07, 41:09, 41:12, 41:14, 41:16, 41:20-41:22, 41:25-41:26, 41:28-41:30, 41:32-41:35, 41:37, 41:44, 41:50, 41:53-41:54, 41:57, 41:59, 41:61, 41:64, 41:66-41:67, 41:73, 45:08, 50:45, 50:49, 54:30, 55:18, 56:23, 56:54			3												15	
*15:202, 15:239, 15:376, 35:10, 35:13, 35:16, 35:28:01:01-35:28:01:03, 35:80, 35:226, 35:270, 35:318, 35:477, 40:185, 40:308, 40:365			3			6						12				
*15:509, 46:50			3						9							
*18:04:01, 35:03:27, 35:42:02			3				7					12				
*18:09, 18:54, 18:67, 18:136, 18:194:01:01-18:194:01:02, 35:471, 37:08, 38:06-38:07, 51:341, 53:01:01:01-53:01:12, 53:01:14-53:04, 53:07-53:15, 53:18, 53:20-53:27, 53:29, 53:31-53:32, 53:34-53:37, 53:40-53:66, 57:120, 58:01:01:01-58:02:02, 58:04-58:07, 58:09-58:26, 58:29, 58:31N-58:35, 58:37-58:43, 58:45:01-58:63, 58:66-58:95, 58:97-58:110, 58:112-58:121, 58:123-58:127, 58:129-58:133N											11					
*18:49, 18:79, 18:108, 18:132, 18:142, 35:203, 35:355												12	13			
*18:98:02							7					12				
*18:205					5							12				
*27:52, 27:236, 37:60									9	11						
*27:102									9						14	
*35:01:10, 35:03:32, 35:04:02, 35:251:01			3		5							12				
*35:27, 35:56			3									12	13			
*35:37, 35:235, 35:265, 35:273			3							10		12				
*35:60, 35:276			3									12			15	
*35:185			3				7					12	13			
*35:329, 53:05:01-53:05:02, 53:16, 53:33, 58:96			3								11					
*35:368			3						9			12				
*37:05						6						12	13			
*37:14						6					11	w				
*37:58, 37:96						6										
*40:01:06			3		5											
*40:394			3												14	
*41:60, 59:02, 59:11															14	15
*42:04, 46:18, 54:01:01:01-54:02, 54:04-54:05N, 54:07-54:08N, 54:10-54:25, 54:27-54:29, 54:31-54:32, 54:34-54:37, 54:39, 54:41-54:43, 55:01:01:01-55:01:06, 55:01:08-55:02:10, 55:02:12-55:03, 55:05, 55:07, 55:10-55:13, 55:15-55:17, 55:19-55:21, 55:23, 55:25-55:30, 55:32-55:48, 55:50, 55:52-55:77, 55:79-55:85, 55:87-55:109, 55:111-55:116										9						15
*44:03:40, 48:01:08, 57:01:32					5											
*44:06, 51:42, 51:45, 51:62, 51:135, 51:220, 51:227, 51:281, 51:319, 53:01:13, 53:06								8			11					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*46:11, 46:32:01-46:32:02, 54:03, 54:33, 54:40, 56:01:01:01-56:01:12, 56:01:14-56:01:17, 56:07-56:08, 56:13-56:16, 56:20:01-56:20:02, 56:22, 56:24-56:30, 56:35-56:41, 56:43, 56:46-56:50, 56:52-56:53, 56:55:01:02, 56:58-56:59, 56:61-56:62, 56:64-56:65, 56:68-56:84									9	10						
*46:40:01		w					7	8	9			12				
*49:03						6		8			11					
*50:58	1		3													
*51:01:01:01-51:01:01:78, 51:01:03-51:01:07, 51:01:10-51:01:14, 51:01:16-51:01:31, 51:01:33-51:01:38, 51:01:41-51:01:42, 51:01:44-51:01:58, 51:01:60-51:01:92, 51:03, 51:07:01, 51:08:01:01-51:08:05, 51:11N-51:12, 51:17-51:18, 51:24:01-51:24:05, 51:26-51:30, 51:32-51:33, 51:35, 51:38-51:39, 51:41N, 51:43-51:44N, 51:48-51:49, 51:51-51:53, 51:55, 51:57-51:58, 51:60-51:61:02, 51:65-51:72, 51:74-51:77, 51:79-51:80, 51:83-51:84, 51:86-51:88, 51:91, 51:94-51:96, 51:98N-51:100, 51:102, 51:105, 51:107, 51:109-51:111, 51:113-51:115, 51:117-51:118N, 51:120-51:121, 51:123-51:131, 51:133-51:134, 51:136-51:138, 51:140-51:142, 51:145, 51:149N-51:151, 51:154-51:156, 51:158:01-51:171, 51:173Q-51:174:02, 51:177-51:178N, 51:181, 51:183-51:184N, 51:186-51:191, 51:193, 51:195-51:196, 51:198, 51:200-51:207, 51:209-51:212, 51:215, 51:219, 51:221-51:222, 51:224, 51:228-51:237:01, 51:238-51:246, 51:248-51:250, 51:252-51:258, 51:260, 51:262-51:266, 51:268-51:273N, 51:275, 51:277-51:279, 51:282-51:283, 51:285-51:292, 51:294, 51:296-51:297, 51:299, 51:301-51:309, 51:311-51:318N, 51:320-51:336, 51:338-51:339, 51:342-51:344N, 58:08:01	1							8		10	11			14		
*51:01:02:01-51:01:02:02, 51:01:09, 51:07:02, 51:89, 51:103, 51:116, 51:119, 51:132, 51:176, 51:179, 51:218, 51:226, 51:237:02, 51:337, 52:18, 52:73, 58:08:02	1									10	11			14		
*51:01:08, 51:01:15, 51:01:43, 51:13:01-51:14, 51:37, 51:50, 51:63:01-51:63:02, 51:85, 51:90, 51:92:01-51:92:02, 51:97, 51:106:01-51:106:02, 51:172, 51:182, 51:192, 51:214, 51:216, 51:223, 51:225, 51:251, 51:261, 51:274, 51:298	1							8		10	11					
*51:01:32, 51:01:39, 51:01:59, 51:16, 51:82, 51:152	1							8			11			14		
*51:01:40, 51:78:01-51:78:02	1									10	11					
*51:02:01:01-51:02:01:03, 51:02:03-51:02:07, 51:23, 51:146, 51:175, 51:199, 51:213, 51:345								8		10	11			14		
*51:02:02, 51:180										10	11			14		
*51:04:01-51:04:02, 51:06:01:01-51:06:03, 51:10, 51:34, 51:46, 51:56:01-51:56:03, 51:59, 51:64, 51:81, 51:293								8			11			14		
*51:05					w					10	11			14		
*51:06:04, 51:144, 51:148, 53:19, 58:27-58:28:02, 58:65											11			14		
*51:09:01, 51:09:03, 51:19, 51:73, 51:122, 51:259, 51:276	1				5			8		10	11			14		
*51:09:02	1				5					10	11			14		
*51:15, 58:122										10	11					
*51:20, 51:208	1	w						8		10	11			14		
*51:21, 51:101, 51:153	1							8			11			14	15	
*51:22	1			4				8		10	11			14		
*51:31	1				5			8			11			14		
*51:36, 51:143, 51:295								8			11			14	15	
*51:40, 51:108					5			8		10	11			14		
*51:54					w					10	11					
*51:93					5						11			14		
*51:104:01-51:104:02, 51:112, 51:147, 51:267, 51:284, 52:02:01	1							8		10				14		
*51:139, 51:194					5			8			11			14		
*51:157, 51:197, 51:300								8			11				15	
*51:185	1		3					8		10	11			14		
*51:217	1				w			8		10	11			14		
*51:247	1		3					8		10	11					
*51:280, 51:340	1				5			8			11			14	15	
*51:310						6		8			11			14		
*52:01:01:01-52:01:01:17, 52:01:03, 52:01:05-52:01:08, 52:01:10-52:01:18, 52:01:20, 52:01:22-52:01:28, 52:01:30-52:01:39, 52:01:41-52:01:42, 52:01:44-52:01:46, 52:01:48, 52:01:50-52:01:51, 52:04-52:05, 52:06:03-52:07:02, 52:13, 52:15, 52:17, 52:19, 52:21:02-52:24, 52:27-52:31:01, 52:33-52:36, 52:38, 52:40-52:61, 52:64-52:66, 52:69-52:70, 52:72, 52:76-52:83, 52:86-52:87, 52:89N-52:96N, 52:98-52:101, 52:103N, 52:105-52:106	1							6			10	11			14	
*52:01:02:01-52:01:02:07, 52:01:04, 52:01:09, 52:01:40, 52:01:43, 52:01:47, 52:01:49, 52:06:01-52:06:02, 52:21:01, 52:31:02, 52:85, 52:102, 52:104	1							6	8		10	11			14	
*52:01:19	1							6			10				14	
*52:01:21, 52:08	1							6			11				14	
*52:01:29, 52:68	1		3					6			10	11				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*52:02:02, 52:16, 52:32	1									10				14		
*52:03					5	6		8		10	11			14		
*52:09:01	1		3			6		8		10	11			14		
*52:09:02	1		3			6				10	11			14		
*52:10:01-52:10:02, 52:74	1				5	6				10	11			14		
*52:10:03	1				5	6		8		10	11			14		
*52:11:01						6				10	11			14		
*52:11:02						6		8		10	11			14		
*52:12, 52:14:01-52:14:02, 52:20, 52:25:01-52:26, 52:39, 52:62-52:63, 52:75	1					6				10	11					
*52:37, 52:84	1					6					11			14	15	
*52:67						6					11			14		
*52:71	1				5	6					11			14	15	
*52:97	1					6		8		10	11					
*55:02:11					5				9						15	
*55:31			3						9						15	
*56:05:01	1	2					7		9	10						
*56:05:02, 56:34	1								9	10						
*56:06	1	2					7		9	10		12				
*56:21	1								9	10				14		
*56:45									9	10		12				
*57:01:20								8								
*73:01:01:01-73:03	1								9							
A*23:01:14		w						8								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control