

*07:13, 07:15, 67:02															18	19																																					
*07:28, 07:118, 27:32, 40:180, 40:414, 44:54, 44:135															17																																						
*07:29			3											15														39																									
*07:32, 07:207							8														24																41																
*07:55, 07:100, 07:113, 15:07:01:01-15:07:03, 15:45, 15:68, 15:126, 15:207, 15:331, 15:405, 15:431, 15:450, 46:12, C*02:60																																							43														
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48					
*07:56:01-07:56:02, 07:72, 14:54, 38:28, 39:63, 39:135, 49:03, 49:48, 58:108																																								40													
*07:86, 14:11, 15:03:01:01-15:03:08, 15:09:01-15:10:04, 15:18:01:01-15:18:06, 15:23, 15:37, 15:47:01-15:47:02, 15:49, 15:52, 15:54, 15:61-15:62, 15:64:01-15:64:02, 15:69, 15:72, 15:74, 15:80, 15:90-15:91, 15:93, 15:98-15:99, 15:103, 15:108, 15:114-15:115, 15:119, 15:123-15:124, 15:127, 15:131-15:134, 15:151, 15:153, 15:156, 15:158, 15:161, 15:173, 15:176, 15:186, 15:197:01-15:198, 15:200, 15:210, 15:220:01:01-15:221, 15:226N, 15:229, 15:235, 15:238, 15:242:01-15:243, 15:251:01-15:253, 15:255, 15:263, 15:266, 15:274-15:275:02, 15:281-15:282, 15:290, 15:292-15:294N, 15:306-15:307, 15:311-15:314, 15:323, 15:329, 15:335, 15:337-15:338, 15:351, 15:354, 15:369, 15:376, 15:379-15:380N, 15:388-15:389, 15:397, 15:414, 15:427, 15:433-15:436, 15:439, 15:444-15:445, 15:448, 15:451, 15:455, 15:475, 15:487N, 15:498, 15:502-15:503, 15:506, 15:509, 15:517, 18:25, 18:68, 38:33, 39:04, 40:12:01-40:12:02, 44:258, 48:22, 52:16, 52:27, 54:21, 56:16, 57:58, 57:80, 58:12															w																																						
*07:92-07:93, 14:63, 38:27, 39:27, 39:143, 40:73																		18																																			
*07:102, 38:79																				21																																	
*07:160																			19																																		
*07:186			3																																									39									
*07:195, 18:147, 40:281																													27																								
*07:206							8																																														
*07:253																		w																																			
*07:304							8											18																																			

Negative Control

*07:349, 41:10, 41:15, 41:52, 42:05:01-42:05:02, 42:08, 42:14	1																																																						
*07:352							8							w																																									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48							
*13:39, 13:67, 15:57, 15:135, 55:12, 56:41, 78:06, C*05:30, C*16:144																						23																																	
*13:46, 15:83, 15:179:01-15:179:02, 15:199, 15:218Q, 15:518, 18:56, 18:113, 27:47, 27:127, 37:01:01-37:01:16, 37:03N-37:06:02, 37:08, 37:10-37:11, 37:13-37:20, 37:22-37:38, 37:41-37:54, 37:56, 37:58-37:61, 37:63-37:65, 37:68-37:80, 37:82N, 37:84-37:87, 38:30, 44:02:01-44:02:59, 44:05:01-44:05:04, 44:08-44:09, 44:11-44:12, 44:14-44:25, 44:27:01-44:27:03, 44:33-44:34:02, 44:41:01-44:42, 44:44, 44:48, 44:51-44:52N, 44:55, 44:58N-44:59:02, 44:62-44:63, 44:66-44:68, 44:70-44:73, 44:78, 44:80, 44:83-44:84:02, 44:86-44:91, 44:93, 44:95, 44:97, 44:99-44:102, 44:104, 44:112-44:113, 44:116-44:119, 44:123, 44:126:01-44:127, 44:131-44:134, 44:136-44:140, 44:142, 44:145, 44:148-44:149N, 44:151-44:152, 44:162, 44:168-44:173, 44:176-44:177, 44:179, 44:185, 44:187, 44:190-44:191, 44:195N-44:196, 44:200-44:201, 44:206, 44:208, 44:211-44:212, 44:214, 44:216-44:220, 44:225-44:226, 44:229-44:230, 44:235-44:236, 44:238, 44:240-44:244, 44:249, 44:253, 44:255, 44:257, 44:260-44:265, 44:267N, 44:269-44:270:02, 44:273-44:274, 44:279, 44:287-44:289, 44:291-44:293, 44:296-44:297, 44:301, 44:303N, 44:306N, 44:308, 44:311-44:316, 44:320-44:322, 44:326-44:327, 44:329-44:331, 44:334N-44:335, 44:337-44:338, 44:341N, 44:344, 44:437, 45:01:01-45:25, 52:19, 55:56, 56:13, 57:09, 57:24, 57:119, 82:01:01-82:03, C*07:626								9																																															
*13:93				4	6																	20																																	
*14:28			3												15																																					44			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48							

Negative Control

*18:07:01-18:07:02, 18:29, 18:79, 18:136, 27:23, 35:50, 35:84, 35:162, 35:217, 35:231, 35:280, 35:323, 40:07-40:08, 40:13, 40:68, 40:106, 40:162, 40:220, 40:232, 51:103, 51:116, 51:151																					41				
*18:33, 27:12:01:01-27:12:01:02, 27:16, 27:92, 27:119, 27:153, 39:31:01:01-39:31:01:02, 39:128, 40:75, 40:322, 49:26																									44
*27:83					10		16						25										37		
*35:18, 35:44, 35:413, 51:05, 51:29, 51:54, 51:61:01-51:61:02, 51:82, 51:217, 51:279				8									24												
*35:27, 35:56, 35:185, 35:203, 35:355, 78:09-78:10													23	24											
*35:38, 35:115, 35:169, 35:260, 44:06, 51:08:01:01-51:08:04, 51:20, 51:36, 51:44N, 51:97, 51:153, 51:263, 51:274-51:275, 53:22				9									24												
*35:87	1												24										32		
*35:89													24			30							33		
*35:192, 35:240, 51:277													24											37	
*35:218, 35:256		2	3	5												29									
*37:21, 37:67																								41	
*37:57																									44
*37:62				6	9														20						29
*38:19, 39:03:01:01-39:03:01:02, 39:14:01:01-39:14:01:02, 39:24:01-39:24:02, 39:29, 39:37, 39:76, 39:120, 39:144																								33	
*39:12:01:01-39:12:01:02			3																						44
*39:78																									12
*40:21, 40:58, 40:137, 40:248, 40:282, 40:325, 40:351, 40:364																									
*40:228																									
*40:398																									32
*41:01:01:01-41:01:05, 41:05, 41:07, 41:14, 41:16, 41:20, 41:25-41:26, 41:28, 41:33-41:35, 41:37, 41:44, 41:53-41:54, 41:56-41:57, 41:59-41:61, 42:04, 55:20																									32
*41:02:01:01-41:02:08, 41:11, 41:13, 41:23, 41:27, 41:31, 41:38-41:43, 41:45N, 41:47, 41:55, 41:58, 42:01:01-42:02:01:02, 42:07, 42:11-42:12, 42:15, 42:22-42:25, 44:166																									32
*41:03:01-41:03:02																									32
*41:12																									40
*41:17																									32
*41:19	1																								32
*41:21					9	10	11																		37

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48
----------	---	---	---	---	---	---	---	---	---	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----	----