

INTERPRETATION TABLE				
HLA-B*82 SSP subtyping				
Amplification patterns of the B*82 alleles				
	Well			
	1	2	3	4
Length of spec.	195	135	230	195
PCR product				
Length of int.	800	1070	1070	1070
pos. control¹				
5'-primer²	420	557	368	105
	^{5'} -TTA ^{3'}	^{5'} -ggA ^{3'}	^{5'} -gTT ^{3'}	^{5'} -gCT ^{3'}
3'-primer³	572	3rd I	557	259
	^{5'} -gCg ^{3'}	^{5'} -TAT ^{3'}	^{5'} -ggC ^{3'}	^{5'} -gTT ^{3'}
Well No.	1	2	3	4
HLA-B allele				
*8201	1	2		4
*8202	1		3	4
*1506, 152701-152703, 1584, 9509			3	
*4410, 4415, 4418, 4501, 4505-4507, 5002	1			
HLA-B allele				
Well No.	1	2	3	4
Cw*0408, 0434			3	
Well No.	1	2	3	4