

INTERPRETATION TABLE			
HLA-B*82 SSP subtyping			
Amplification patterns of the B*82:01 to 82:03 alleles			
	Well		
	1	2	3
Length of spec.	195	140	230
PCR product			
Length of int.	800	800	1070
pos. control¹			
5'-primer²	420	557	368
	5'-TTA 3'	5'-ggA 3'	5'-gTT 3'
3'-primer³	572	3rd I	557
	5'-gCg 3'	5'-TAT 3'	5'-ggC 3'
Well No.	1	2	3
HLA-B allele⁴			
*82:01	1	2	
*82:02	1		3
*82:03	1	2	
*15:06, 15:27:01-15:27:03, 15:84, 15:109, 15:195, C*03:89, C*04:08, C*04:34, C*14:20			3
*44:10, 44:15, 44:18, 44:140, 45:01, 45:05-45:07, 45:11-45:14, 49:20, 50:02	1		
*51:01:28			
HLA-B allele			
Well No.	1	2	3

4	5
195	155
1070	1070
105	105
5'-gCT 3'	5'-gCT 3'
259	219
5'-gTT 3'	5'-ggg 3'
259	
5'-gTT 3'	
4	5
4	
4	
4	5
4	
4	5