

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
HLA-C allele ^{1,2}																																	
*01:02:01-01:02:17, 01:02:19-01:02:24, 01:06-01:08, 01:11-01:13, 01:15-01:20, 01:23, 01:25-01:33, 01:37N-01:42, 01:44-01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:80-01:89N	1																			20						26							
*01:02:18	1											12									20						26						
*01:03, 01:24, 01:34, 01:78	1																				20	21					26						
*01:04	1											12		14																			
*01:05, 01:22, 01:36, 01:79:01-01:79:02, B*39:76, B*54:18 ³	1																																
*01:09	1													14																			
*01:10	1	2																			20						26						
*01:14	1											11									20						26					31	
*01:21	1												12	13														27			30		
*01:35	1																													29			
*01:43	1	2										10									20						26						
*01:49-01:50, 01:55	1																					21											
*01:59	1											11									20						26						
*01:60	1																				19	20					26						
*01:67	1															16					20						26						
*01:90	1							8													17						20					26	
*02:02:01		2	3			6					11				w										24							31	
*02:02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:15-02:02:25, 02:07-02:09, 02:11, 02:15, 02:19-02:20, 02:24-02:25Q, 02:28-02:31, 02:34-02:35, 02:37-02:38N, 02:40, 02:42-02:44, 02:46, 02:48, 02:50, 02:52N-02:54, 02:56-02:57, 02:59-02:61, 02:63-02:64, 02:66-02:67Q, 02:69-02:71, 02:74		2	3			6					11				15										24								31
*02:02:05, 02:02:14, 02:10, 02:32, 02:51, 02:73		2	3			6					11				15																		31
*02:02:13		2	3								11				15																		31
*02:02:26, 02:21, 02:26:02-02:26:03, 02:39, 02:45		2	3								11				15									24									31
*02:03, 02:16:01-02:16:02, 02:72		2	3			6					11				15																		
*02:04		2	w			6					11				15									24									31
*02:05:01, 02:17		2	3			6					11			14	15																		31
*02:05:02-02:05:03		2	3			6					11			14	15									24									31
*02:06		2	3			6	8				11				15		17																31
*02:12		2	3			6							13	15														w		30	31		
*02:13		2	3			6					11				15		18							24									31
*02:14		2	w			6					11																						31
*02:18		2	3			6					11				15		18																
*02:22, 02:62		2	3			6					11				15						21			24									31
*02:23, 02:36, 02:68		2	3			6					11				15		17							24									31
*02:26:01		2									11				15									24									31
*02:27:01		2	3			6																		24				27					
*02:27:02		2	3			6																						27					
*02:33		2	3			6					11				15		18																31
*02:47		2	3			6	8				11				15									24									31
*02:49		2	3			6					11		13	15			18							24						30	31		
*02:55		2	3			6					11		13	15										24						30	31		
*02:58		2	3			6					11			15															29			31	
*02:65		2	3			6					11													24									
*03:02:01-03:02:09, 03:02:11-03:02:12, 03:05, 03:17, 03:25, 03:27, 03:33, 03:35-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226			3	4																			22										
*03:02:10				4																													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control

