

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
*02:55:02		2	3			6					11		13		15											26				30	31			
*02:58			3			6					11				15																	31		
*02:65		2	3			6					11													24										
*02:68		2	3			6					11				15		17							24								31		
*02:87		2								10														24			27							
*02:94		2	3			6	7				11				15																	31		
*02:101		2	3			w					11				15									24								31		
*02:107		2									11													24								31		
*02:115		2	3			6					11		13					18						24		26	27			30				
*02:126		2	3			6					11				15									24		24		w				31		
*03:02:01-03:02:09, 03:02:11-03:02:12, 03:02:14-03:02:17, 03:05, 03:17:01-03:17:02, 03:25, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:335, 03:338, 03:349-03:350, 03:365			3	4																			22											
*03:02:10, 03:264				4																														
*03:02:13, 03:16, 03:110, 03:332			3																				22											
*03:03:01:01-03:03:20, 03:03:22-03:03:37, 03:11:01-03:11:02, 03:18:02, 03:20N-03:22Q, 03:30, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:273, 03:275-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:351-03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367				4	5																	21	22											
*03:03:21, 03:04:26, 03:18:01, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182, 03:195, 03:344				4																			21											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*06:77								8			11			14									23								31	
*06:82, 06:160, 06:183-06:184								8			11			14									23									
*06:89								8			11			14			17						23								31	
*06:118								8			11			14			18	19					23									
*06:120								8			11			14									23		25						31	
*06:127:02								8			11												23							30	31	
*06:132:02			3								11			14	15								22			26					31	
*07:01:01:01-07:01:06, 07:01:08-07:01:40, 07:01:42-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:04:02:02, 07:04:04-07:06, 07:08, 07:11-07:13, 07:15-07:19, 07:21-07:25, 07:27:01, 07:29:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:138, 07:140-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:199:01, 07:200-07:209, 07:211-07:215, 07:217-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:262, 07:264N-07:266, 07:268-07:293, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:401, 07:404, 07:407-07:414, 07:417-07:448, 07:450-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518-07:525, 07:527-07:530, 07:532-07:535, 07:537-07:562, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599-07:600N									9													23										
*07:01:07, 07:04:03, 07:14, 07:27:02, 07:199:02, 07:526									9														23								30	
*07:01:41, 07:02:41, 07:415-07:416									9														w									
*07:02:08, 07:02:29, 07:33N, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:322, 07:405, 07:493, 07:517									9																							
*07:02:09						6			9				13				18						23				27					
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598									9		11										20		23									
*07:10, 07:28, 07:41, 07:43:01, 07:184, 07:196, 07:367, 07:536, 07:568									9													21	23									
*07:20, 07:263									9														23								29	
*07:26:01-07:26:02, 07:92, 07:314:01, 07:314:03, 07:317, 07:351									9																26							
*07:31:01, 07:514									9				14				18						23									
*07:31:02									9				14				18						23								30	
*07:39-07:40, 07:335									9													20										
*07:43:02									9														21	23							30	
*07:53, 07:216									9		11												23								30	
*07:64, 07:402									9						16								23								30	
*07:96:01-07:96:02				3					9														22			26				29		
*07:101, 07:148, 07:161			2						9	10													23			26						
*07:123, 07:173, 07:294									9								17						23									
*07:125, 07:356, 07:531						6			9		11		13	14										23								
*07:154									9								18						23									
*07:177									9				14				18															
*07:210, 07:238									9		11											20										
*07:242				4					9															23								
*07:314:02, 07:578				3					9														22			26						
*07:316, 07:338, 07:579		1							9														23									
*07:449									9																						27	
*07:563									9															23								
*07:583			2				6		9	10			13				18									26	27					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
*14:15, 14:87														14	16		18															
*14:16														14	16																	
*14:17															16		18	19	20							26		28				
*14:25		2													16				20		22				26		28					
*14:45															15	16			20						26							
*14:54															16						21						27					
*14:58	1														16					20					26		28					
*14:77										11					16																	
*15:02:01:01-15:02:01:07, 15:02:02-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:20, 15:02:22-15:02:24, 15:02:26-15:02:30, 15:06:01-15:06:03, 15:09, 15:13:01:01-15:13:01:02, 15:15, 15:18-15:19, 15:26, 15:28, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44:01-15:45, 15:47-15:53, 15:56-15:57, 15:60, 15:64, 15:67-15:68, 15:73, 15:78:01-15:84Q, 15:86-15:89, 15:91-15:101, 15:103, 15:106, 15:109, 15:112, 15:114, 15:119-15:123, 15:126-15:137, 15:139, 15:141-15:142, 15:145N-15:146			3					8		11					15	17						21									29	
*15:02:01:08N, 15:55, 15:58, 15:72			3					8		11					15	17															29	
*15:02:10, 15:02:17			3					8		11					15	17					21	22								29		
*15:02:21, 15:02:25, 15:61								8		11					15	17					21									29		
*15:03			3					8		11		13			15	17					21				26				29			
*15:04:01-15:04:03, 15:24, 15:30, 15:40, 15:66, 15:76-15:77			3							11					15	17						21								29		
*15:05:01-15:05:08, 15:05:10, 15:29, 15:46, 15:54, 15:59, 15:69-15:70, 15:90, 15:105Q, 15:108, 15:110, 15:115N, 15:117, 15:140			3					w		11					15	17						21								29		
*15:05:09								w		11					15	17						21								29		
*15:05:11-15:05:12			3					w							15	17						21								29		
*15:07			3					8								17						21						27		29		
*15:08, 15:74, 15:143			3					8		11					15	17						21							29	31		
*15:10:01								8		11					15	17						21										
*15:10:02-15:10:03, 15:17, 15:27, 15:113, 15:124			3					8		11					15	17						21										
*15:11			3		6			8		11					15	17						21										
*15:12				5				8		11					15	17	18					21								29		
*15:16			3					8		11			13		15	17						21			26							
*15:21			3					8							15	17						21					w		29			
*15:22, 15:125			3					w		11					15	17													29			
*15:23								w		11						17						21								29		
*15:25			3									13					18	19				21				26	27					
*15:31, 15:62			3					8		11					15	17	18					21							29			
*15:36					6			w		11						17						21								29		
*15:37			3					8		11					15	17												28	29			
*15:42		2	3					8		11					15	17						21							29			
*15:43			3					8								17						21	22									
*15:63								8		11						17						21										
*15:65			3							11					15	17	18	19											29			
*15:71			3					8							15	17						21							29			
*15:75			3					8		11			13		15	17	18					21							29	30		
*15:85			3					8		11						17						21							29			
*15:102			3							11					15	17						20				26		28	29	31		
*15:104	1		3					w		11					15	17						21							29			
*15:107								7	8		11				15	17						21										
*15:111			2	3				w		11					15	17						21							29			
*15:116			3					w							15	17						21				w		29				
*15:118			3			6		8		11					15	17						21						29				
*15:138								8		11						17						21						29				
*15:144					5			8								17	18					21					27	29				
*16:01:01:01-16:01:01:04, 16:01:03-16:01:15, 16:01:18-16:01:22, 16:06-16:08, 16:11, 16:13-16:14, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56, 16:58-16:59, 16:65, 16:72-16:73, 16:75-16:76, 16:79-16:81, 16:83, 16:86, 16:92-16:93, 16:95-16:98, 16:100, 16:105, 16:111-16:112, 16:114										11			13						18								27			30		
*16:01:02, 16:57, 16:85, 16:94, 16:106											11		13						18											30		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control

