



*02:55:02		2	3		6				11	13	15										26			30	31							
*02:58			3		6				11		15														31							
*02:65		2	3		6				11												24											
*02:68		2	3		6				11		15	17									24				31							
*02:87		2							10												24		27									
*02:94		2	3		6	7			11		15														31							
*02:101		2	3		w				11		15										24				31							
*02:107		2							11												24				31							
*02:115		2	3		6				11	13		18									24	26	27		30							
*02:126		2	3		6				11		15										24		w		31							
*02:164		2			6				11												24				31							
*02:180		2	3		6				11		15	18	19																			
*02:187		2	3		6				11												24				31							
*03:02:01-03:02:09, 03:02:11-03:02:23, 03:05, 03:17:01-03:17:02, 03:25, 03:33, 03:35:01-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:110, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258, 03:292, 03:296:01-03:296:02, 03:298-03:301, 03:315, 03:326, 03:329-03:330, 03:332, 03:335, 03:338, 03:349-03:350, 03:365, 03:371, 03:373, 03:386, 03:391N, 03:402-03:403, 03:406-03:407, 03:411-03:412, 03:425, 03:429, 03:449N, 03:452, 03:464, 03:468, 03:473, 03:475, 03:482, 03:489, 03:493-03:494, 03:497, 03:505			3	4																22												
*03:02:10, 03:264, B*15:96, B*44:274 <sup>3</sup>				4																												
*03:03:01:01-03:03:20, 03:03:22-03:03:50, 03:03:53-03:03:58, 03:11:01-03:11:02, 03:18:02-03:18:03, 03:20N-03:22Q, 03:30, 03:43:01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116:01-03:116:02, 03:119:01-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:273, 03:275-03:276, 03:284-03:285, 03:288-03:291, 03:304, 03:307, 03:312, 03:316N, 03:319-03:321, 03:325, 03:327, 03:336, 03:339, 03:341, 03:345-03:346, 03:352, 03:356-03:357, 03:360, 03:363N-03:364, 03:367, 03:370, 03:372:01-03:372:02, 03:374-03:375, 03:377N-03:378, 03:380N, 03:383, 03:389, 03:395, 03:398, 03:413-03:414, 03:418, 03:421N-03:422, 03:427-03:428, 03:430, 03:432N-03:433, 03:436, 03:438, 03:440, 03:444N, 03:447N, 03:457, 03:460, 03:471-03:472, 03:476-03:477, 03:481, 03:490, 03:495, 03:503, 03:506			4	5																21	22											
*03:03:21, 03:03:51, 03:04:26, 03:18:01, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182, 03:195, 03:344:01:01-03:344:01:02, 03:454, 03:459, 03:504				4																					21							
*03:03:52, 03:151, 03:165:01-03:165:02, 03:251, 03:267, 03:451					5																				21	22						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control







Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
*06:127:02, 06:264								8			11											23								30	31			
*06:132:02			3								11			14	15							22				26					31			
*06:204	1							8			11									20		23			26						31			
*07:01:01:01-07:01:02:08, 07:01:04-07:01:06, 07:01:08-07:01:40, 07:01:42-07:01:73, 07:01:75-07:01:92, 07:02:01:01-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:02:90, 07:02:92-07:04:02:02, 07:04:04-07:06:03, 07:08, 07:11-07:13:01:02, 07:15-07:19, 07:21-07:25, 07:27:01, 07:29:01:01-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93:01-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:132, 07:133:02-07:138, 07:140-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:199:01, 07:200-07:209, 07:211-07:215, 07:217-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:262, 07:264N-07:266, 07:268-07:292, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N, 07:330:02-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368:01-07:401, 07:404, 07:407-07:414, 07:417-07:448, 07:450-07:492, 07:494Q-07:513Q, 07:515-07:516, 07:518, 07:520-07:525, 07:526:02-07:530, 07:532-07:535, 07:537-07:562, 07:564-07:567, 07:569-07:577, 07:580-07:582Q, 07:584-07:597, 07:599-07:625, 07:627-07:680, 07:682:01:01-07:713Q, 07:715-07:722, 07:726N-07:747N, 07:749N-07:755, 07:757-07:774, 07:776N-07:820N, 07:822-07:844								9														23												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
*07:01:07, 07:04:03, 07:14, 07:27:02, 07:199:02, 07:526:01, 07:681									9														23									30		
*07:01:41, 07:02:41, 07:415-07:416									9													w												
*07:01:74						6			9				13				17	18					23				27							
*07:01:93, 07:02:08, 07:02:29, 07:02:91, 07:33N, 07:139, 07:152N, 07:156, 07:246:01-07:246:02, 07:267, 07:296, 07:322, 07:405, 07:493, 07:517, 07:714, 07:724									9				13				18						23					27						
*07:02:09						6			9				13				18						23					27						
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406, 07:598									9		11										20		23											
*07:10, 07:28, 07:41, 07:43:01, 07:184, 07:196, 07:367, 07:536, 07:568									9													21	23											
*07:20:01-07:20:02, 07:263									9														23								29			
*07:26:01:01-07:26:03, 07:92, 07:314:01, 07:314:03, 07:317, 07:351, 07:821N, 07:845									9																	26								
*07:31:01, 07:514									9					14			18						23											
*07:31:02									9					14			18						23									30		
*07:39-07:40, 07:335									9												20													
*07:43:02									9													21	23									30		
*07:53, 07:216									9		11												23										30	
*07:64, 07:402									9							16							22											
*07:96:01-07:96:02			3						9														22			26			29					
*07:101, 07:148, 07:161		2							9	10													23			26								
*07:123, 07:173, 07:294, 07:626, 07:725, 07:775									9								17						23											
*07:125, 07:356, 07:531							6		9		11		13	14									23											
*07:133:01, 07:242, 07:330:01, 07:748				4					9														23											
*07:154									9									18					23											
*07:177									9					14			18			20			23											
*07:210, 07:238									9		11										20													
*07:293, 07:519									9							16							23											
*07:314:02, 07:578, 07:756			3						9														22			26								
*07:316, 07:338, 07:579	1								9														23											
*07:449									9																							27		
*07:563									9												20		23											
*07:583		2				6			9	10			13				18										26	27						
*07:723						6			9				13														26	27						
*08:01:01:01-08:01:17, 08:01:20-08:01:29, 08:03:01-08:03:04, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22:01:01-08:22:01:02, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:105-08:106, 08:109, 08:117-08:119, 08:121N-08:122, 08:124, 08:127N-08:131, 08:133, 08:135-08:136, 08:141Q, 08:143-08:145, 08:147-08:149, 08:153-08:155, 08:157, 08:162, 08:164, 08:168, 08:173N-08:176, 08:178, 08:186-08:187, 08:189-08:190,										10							18				21						27							

Negative Contr

Negative Control









*16:25											11	13	15	18										26			30							
*16:29		2											14	18	20									25	27									
*16:34			3		6						11	13		18											27		30							
*16:35													14	18											27									
*16:37											11	13	15	18												w		30						
*16:40, 16:110					7																				27									
*16:43											11	13													27		30							
*16:45												13		18					21						27		30							
*16:48											11		14	15	18																			
*16:50											2		11	13	18											27		30						
*16:53															16										26	27								
*16:62												11	13		18									25	27		30							
*16:64												11	13		17	18									26	27		30						
*16:67													13		18										27									
*16:68															18	20									26	27								
*16:70												3		11	13	15	17	18									29	30						
*16:71														11	13		18									w		30						
*16:87													11	13		17	18								27		30							
*16:88													11		15	18	20								26				31					
*16:91													11		14	15	18	20											31					
*16:103													11		15	18																		
*16:109															14		17	18	20							27								
*16:113													7	11	13	16										27		30						
*16:121														3	11	13	15	18										30						
*16:139															10	11	13	18							24	27		30						
*16:140															11		15																	
*16:150																14	16	18	20						25	27								
*17:01:01:02-17:01:03, 17:01:05-17:01:09, 17:01:11-17:04, 17:06, 17:08-17:16:02, 17:18-17:19, 17:21, 17:24-17:26, 17:28-17:32, 17:34-17:39, 17:41-17:51:01:02																																		
*17:01:04																																		
*17:01:10																																		
*17:05																																		
*17:07, 17:23, 17:33																																		
*17:17																																		
*17:20																																		
*17:22																																		
*17:27N																																		
*17:40																																		
*18:01:01:01-18:02:04, 18:04-18:08, 18:10-18:13																																		
*18:03																																		
*18:09																																		
*18:14																																		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
A*02:42:01, B*15:539, B*39:76, B*54:18																																		
A*24:96, A*24:106																																		
A*24:174, B*57:129																																		
A*25:55, B*07:13, B*07:15, B*07:160																																		
B*07:02:77, B*08:01:36, B*15:01:50, B*15:436, B*18:03:02, B*35:08:02, B*35:08:05																																		
B*07:239, B*14:03																																		
B*08:145																																		
B*13:13:02, B*40:01:49-40:01:50																																		
B*15:125:02, B*44:02:44, B*56:01:09, B*58:74																																		
B*15:277, B*15:430, B*35:222, B*42:22																																		
B*15:510																																		
B*54:38																																		
B*57:01:24																																		
B*58:02:01:01																																		
B*73:01:01:01-73:02																																		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		

Negative Control

Negative Control