

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-C allele ^{1,2}																								
*01:02:01-01:02:17, 01:02:19-01:02:35, 01:06-01:07:01, 01:08, 01:11-01:13, 01:15-01:20, 01:23, 01:25-01:33, 01:37N-01:42, 01:44-01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:66, 01:68-01:77, 01:80-01:89N, 01:91-01:100, 01:103-01:107	1																			20				
*01:02:18	1											12								20				
*01:03, 01:24, 01:34, 01:78	1																			20	21			
*01:04	1											12		14										
*01:05, 01:07:02, 01:22, 01:35-01:36, 01:79:01-01:79:02, B*39:76, B*54:18 ³	1																							
*01:09	1													14										
*01:10	1	2																		20				
*01:14, 01:59	1										11									20				
*01:21	1											12	13											
*01:43	1	2								10										20				
*01:49-01:50, 01:55	1																				21			
*01:60	1																		19	20				
*01:67, 14:58 ³	1															16				20				
*01:90	1							8									17			20				
*01:101-01:102	1																17			20				
*02:02:01		2	3			6					11				w									
*02:02:02:01-02:02:03, 02:02:05-02:02:12, 02:02:14-02:02:25, 02:02:27-02:02:28, 02:03, 02:07-02:11, 02:15-02:16:02, 02:19-02:20, 02:24-02:25Q, 02:28-02:32, 02:34-02:35, 02:37-02:38N, 02:40-01-02:40:02, 02:42-02:44, 02:46, 02:48, 02:50-02:54, 02:56-02:57, 02:59-02:61, 02:63-02:64, 02:66-02:67Q, 02:69-02:74, 02:76-02:86, 02:88-02:93, 02:95		2	3			6					11				15									
*02:02:13, 02:02:26, 02:21, 02:26-02:02:26:03, 02:39, 02:45		2	3								11				15									
*02:02:29		2	3								11				15						22			
*02:02:30		2				6					11				15									
*02:04		2	w			6					11				15									
*02:05:01-02:05:03, 02:17		2	3			6					11			14	15									
*02:06		2	3			6		8			11				15		17							
*02:12		2	3			6							13		15									
*02:13, 02:18, 02:33, 02:75		2	3			6					11				15			18						
*02:14		2	w			6					11													
*02:22, 02:62		2	3			6					11				15					21				
*02:23, 02:36, 02:68		2	3			6					11				15		17							
*02:26:01		2									11				15									
*02:27:01-02:27:02		2	3			6																		
*02:47		2	3			6		8			11				15									
*02:49		2	3			6					11		13		15			18						
*02:55		2	3			6					11		13		15									
*02:58			3			6					11				15									
*02:65		2	3			6					11													
*02:87		2								10														
*02:94		2	3			6	7				11				15									
*03:02:01-03:02:09, 03:02:11-03:02:12, 03:05, 03:17, 03:25, 03:27, 03:33, 03:35-03:36, 03:60, 03:71, 03:84, 03:89, 03:95, 03:108, 03:121N, 03:135, 03:139, 03:146, 03:154, 03:169Q, 03:178, 03:190, 03:194, 03:197-03:198, 03:200-03:201N, 03:216, 03:221-03:222, 03:224N-03:226, 03:238, 03:249, 03:258			3	4																	22			
*03:02:10, 03:264				4																				
*03:02:13, 03:16, 03:110			3																		22			
*03:03:01-03:03:20, 03:03:22-03:03:28, 03:11-01-03:11:02, 03:18:02, 03:20N-03:22Q, 03:30-03:31, 03:43-01-03:43:02, 03:49-03:50, 03:52-03:53, 03:55-03:56, 03:59, 03:61-03:62, 03:66, 03:68-03:69, 03:75-03:76, 03:79, 03:83, 03:85, 03:88, 03:96, 03:102-03:103, 03:112, 03:116-01-03:116:02, 03:119-03:120, 03:122, 03:124, 03:126-03:127, 03:133, 03:141-03:142, 03:150, 03:152, 03:158, 03:160-03:161, 03:168, 03:171, 03:176, 03:185, 03:187-03:189N, 03:192, 03:196, 03:202-03:207, 03:214, 03:217, 03:220, 03:223, 03:227-03:230, 03:237, 03:241-03:242, 03:253-03:254, 03:272-03:276				4	5															21	22			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*03:03:21, 03:04:26, 03:18:01, 03:39, 03:64:02, 03:137, 03:144, 03:156, 03:166, 03:170, 03:182, 03:195				4																	21			
*03:04:01-01-03:04:06, 03:04:08-03:04:24, 03:04:27-03:04:46, 03:06:01-03:06:02, 03:08-03:09, 03:14, 03:19, 03:23-03:24, 03:26, 03:28, 03:32, 03:34, 03:37:01-03:38:02, 03:40:01-03:40:04, 03:42, 03:44, 03:46-03:48, 03:51, 03:54, 03:57, 03:63-03:64:01, 03:65, 03:70, 03:72-03:74, 03:77-03:78, 03:82, 03:87:01-03:87:02, 03:90-03:93, 03:98, 03:100-03:101, 03:104-03:107, 03:109, 03:111, 03:114-03:115, 03:117-03:118, 03:123, 03:125, 03:128-03:129, 03:131, 03:134, 03:136, 03:138, 03:143, 03:147-03:149, 03:153, 03:157, 03:159, 03:162, 03:164, 03:173-03:174, 03:179-03:181, 03:183-03:184, 03:186:01-03:186:02, 03:191, 03:193, 03:208N, 03:210-03:213, 03:215, 03:218-03:219, 03:232-03:236, 03:239-03:240, 03:244Q, 03:246-03:247, 03:250, 03:252, 03:255-03:257, 03:259, 03:261, 03:263, 03:265N-03:266, 03:269-03:270, 03:277N			3	4																21	22			
*03:04:07, 03:41:01-03:41:02, 03:80, 03:113, 03:260			3																		21	22		
*03:04:25, 03:145, 03:172			3	4																	21			
*03:07, 03:45, 03:163			3	4							11				15						21	22		
*03:10, 03:29, 03:209			3	4											15						21	22		
*03:13:01-03:13:02, 03:97, 03:132, 03:167, 03:177, 03:271				4	5																	22		
*03:15			3	4							11				15							22		
*03:58					5										15					20	21	22		
*03:67				4	w																21	w		
*03:81, 03:175				4	5												17				21	22		
*03:86					5										15					20		22		
*03:94			3												15					20	20	22		
*03:99															15					20				
*03:130, 03:140			3	4							11										21	22		
*03:151, 03:165, 03:251, 03:267					5																21	22		
*03:155			3	4												16					21	22		
*03:199, 03:245			3	4													17					22		
*03:231				4		6										16					21			
*03:243				4	5						11										21	22		
*03:248			3	4																				
*03:262				4	5																?	22		
*03:268				4	5						11				15						21	22		
*04:01:01-01-04:01:15, 04:01:17-04:01:22, 04:01:24-04:01:27, 04:01:29-04:01:60, 04:01:62, 04:01:64, 04:01:66, 04:04:01-04:05, 04:07-04:09N, 04:12-04:13, 04:15:01-04:15:03, 04:17-04:20, 04:24-04:28, 04:30-04:31, 04:34-04:35, 04:38-04:41, 04:43-04:53, 04:56-04:57, 04:59Q-04:67, 04:69-04:70, 04:73-04:76, 04:78-04:79, 04:81-04:102, 04:104-04:106, 04:109-04:111, 04:113, 04:115N, 04:117-04:128, 04:131-04:139, 04:141-04:146, 04:148-04:152, 04:154-04:159, 04:161-04:168, 04:170N, 04:173N-04:177, 04:179-04:189, 04:191N, 04:193					6						11		13					18			21			
*04:01:16, 04:23, 04:54, 04:72					6						11		13					18						
*04:01:23, 04:33, 04:116, 04:130					6						11		13								21			
*04:01:28, 04:01:61, 04:01:63, 04:01:65, 04:140											11		13					18			21			
*04:03:01-04:03:02, 04:06, 04:80, 04:147, 04:171			3		6						11		13		15			18			21			
*04:10, 04:153, 04:192, 04:194					6						11							18			21			
*04:11					6											16					21			
*04:14, 04:68					6						11		13					18			21		23	
*04:16			3								11		13		15			18			21			
*04:29					6								13			16		18			21			
*04:32, 04:77		2			6						11		13					18			21	22		
*04:36, 04:55					6											16		18			21			
*04:37					6						11		13		15			18		20	21			
*04:42:01-04:42:02											11		13	14				18			21		23	
*04:58					6						11		13					18	19		21			
*04:71	1				6						11		13					18			21			
*04:103											11		13		15			18			21			
*04:107			3								11		13		15						21			
*04:108					6						11		13				17	18						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*04:112						6					11		13				17	18			21			
*04:114						6					11		13			16		18			21			
*04:129						6	7				11		13					18			21			
*04:160			3			6					11		13		15			18	19		21			
*04:169, 15:36 ³						6					11						17				21			
*04:172						6							13			16					21			
*04:178						6		8			11		13				17	18			21			
*04:190						6					11		13		15			18			21			
*05:01:01:01-05:01:31, 05:03, 05:05-05:15, 05:18-01:05:19, 05:21, 05:24, 05:27-05:35, 05:37-05:41, 05:44-01:05:45, 05:47-05:54, 05:56-05:57, 05:58-02:05:61, 05:63-05:67, 05:69-05:75, 05:77, 05:80-05:84, 05:86-05:102, 05:104, 05:106, 05:108-05:111, 05:114							7				11				15						21			
*05:04:01-05:04:02, 05:22:01-05:22:02, 05:103							7				11				15									
*05:16, 05:85							7				11				15					20	21			
*05:17, 05:25, 05:55, 05:68, 05:76, 05:79							7				11				15			18			21			
*05:20, 05:43, 05:113N							7				11										21			
*05:23, 05:62							7				11				15				19		21			
*05:26						6	7				11				15						21			
*05:36							7				11				15		17				21			
*05:42							7				11		13		15			18			21			
*05:46							7				11		13		15						21			
*05:58:01			3				7				11				15						21			
*05:78						6	7				11										21			
*05:105		2					7				11				15						21			
*05:107							7				11				15					20				
*05:112							7				11		13		15						21		23	
*06:02:01:01-06:02:01:03, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:37, 06:02:40-06:02:42, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34, 06:36-06:39, 06:41-06:51, 06:53:01-06:66, 06:68-06:71, 06:73-06:75, 06:78, 06:80, 06:82-06:88, 06:90-06:93, 06:96-06:98, 06:100, 06:102-06:117, 06:119-06:121, 06:125-06:126, 06:128N-06:131, 06:133-06:135, 06:137-06:142, 06:145-06:146, 06:148-06:149								8			11			14									23	
*06:02:12, 06:11, 06:52, 06:122, 06:124, 06:147								8						14								23		
*06:02:16, 06:02:43, 06:04, 06:40, 06:79N, 06:81, 06:101, 06:123, 06:127:01-06:127:02, 06:136								8			11											23		
*06:02:38								8			11	12		14								23		
*06:02:39, 06:32, 06:94-06:95, 06:99											11			14								23		
*06:03:01, 06:132:01			3					8			11			14	15						22			
*06:03:02, 06:76:01								8			11			14	15									
*06:05							7	8			11						18		20			23		
*06:06								8			11								20			23		
*06:08		2						8			11			14								23		
*06:09								8			11			14							21	23		
*06:14, 06:35, 06:72, 06:143-06:144								8			11										21	23		
*06:17											11			14					20			23		
*06:31								8			11			14				18	20			23		
*06:67							7	8			11			14								23		
*06:76:02											11	12	13	14	15			18						
*06:77								8			11			14										
*06:89								8			11			14			17					23		
*06:118								8			11			14				18	19			23		
*06:132:02			3								11			14	15						22			
*07:01:01:01-07:01:40, 07:01:42-07:02:07, 07:02:10-07:02:28, 07:02:30-07:02:40, 07:02:42-07:06, 07:08, 07:11-07:25, 07:27-01:07:27:02, 07:29-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:63, 07:65-07:75, 07:77-07:91, 07:93-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:138, 07:140-07:147, 07:149-07:151, 07:153, 07:155, 07:157-07:160, 07:162-07:172:02, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:209, 07:211-07:222, 07:224-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:266, 07:268-07:293, 07:297-07:313, 07:318-07:321, 07:323-07:327, 07:329N-07:334, 07:336-07:337, 07:339-07:350N, 07:352-07:355, 07:357-07:366, 07:368-07:401, 07:404, 07:407-07:409									9													23		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*07:01:41, 07:02:41, 07:223									9															w
*07:02:08, 07:02:29, 07:26, 07:33N, 07:39-07:40, 07:92, 07:139, 07:152N, 07:156, 07:246, 07:267, 07:296, 07:314, 07:317, 07:322, 07:335, 07:351, 07:405									9															
*07:02:09					6				9				13					18						23
*07:07, 07:09, 07:49, 07:76:01-07:76:02, 07:247, 07:315, 07:328, 07:403, 07:406									9	11										20				23
*07:10, 07:28, 07:41, 07:43, 07:184, 07:196, 07:367									9												21			23
*07:31:01-07:31:02									9				14					18						23
*07:64, 07:402									9							16								
*07:96:01-07:96:02			3						9													22		
*07:101, 07:148, 07:161		2							9	10														23
*07:123, 07:173, 07:294									9								17							23
*07:125, 07:356					6				9				13	14										23
*07:154									9									18						23
*07:177									9				14					18						
*07:210, 07:238									9	11										20				
*07:242				4					9															23
*07:316, 07:338	1								9															23
*08:01:01-08:01:17, 08:03:01-08:03:03, 08:06, 08:08:01-08:08:02, 08:14, 08:16:01-08:16:02, 08:20, 08:22, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40-08:42, 08:44, 08:46, 08:50, 08:56-08:61, 08:65, 08:72:01-08:72:02, 08:78-08:82, 08:84-08:89N, 08:91, 08:96-08:99, 08:101-08:102, 08:105-08:106, 08:109										10								18			21			
*08:01:18-08:02:12, 08:04:01-08:04:03, 08:13, 08:15:01-08:15:02, 08:17-08:19, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62-08:63, 08:66-08:71, 08:73-08:77, 08:90, 08:92-08:94, 08:100, 08:103, 08:107, 08:110-08:112, 08:114-08:115										10											21			
*08:05, 08:25										10			13								21			
*08:07, 08:104										10								19			21			
*08:09, 08:11, 08:83										10								18						
*08:10							7			10	11				15			18			21			
*08:12										10										20	21			
*08:21										10			13					18			21			
*08:31		2								10											21			
*08:47										10								18	19		21			
*08:48, 08:108										10														
*08:95										10						16		18			21			
*08:113										10							17							
*12:02:01-12:02:10, 12:02:12-12:02:13, 12:10:01-12:10:02, 12:17-12:18:01, 12:20, 12:22, 12:27, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N, 12:83-12:86, 12:96, 12:103-12:106, 12:112, 12:114, 12:117-12:118, 12:124, 12:126, 12:128, 12:130, 12:134, 12:136-12:137, 12:142, 12:145, 12:148N												12	13											
*12:02:11								8				12												
*12:03:01:01-12:03:01:02, 12:03:03, 12:03:05-12:03:07, 12:03:10-12:03:12, 12:03:14-12:03:15, 12:03:17-12:03:19, 12:03:21-12:03:22, 12:03:25, 12:03:27-12:03:29, 12:03:31-12:03:33, 12:06-12:07, 12:11-12:13, 12:19, 12:23, 12:25, 12:29, 12:32, 12:34-12:35, 12:37-12:39N, 12:42Q-12:43, 12:45-12:48, 12:50-12:53, 12:57, 12:61-12:63, 12:65-12:66, 12:70-12:71, 12:75-12:79, 12:82, 12:87-12:95, 12:97, 12:99-12:102, 12:107-12:111, 12:115-12:116, 12:120-12:122, 12:125, 12:129, 12:131, 12:133, 12:138-12:141, 12:143												12	13	14										
*12:03:02, 12:03:08, 12:03:16, 12:03:20, 12:03:30, 12:26, 12:59													13	14										
*12:03:04, 12:98												12		14										
*12:03:09								8				12		14										
*12:03:13, 12:03:24												12	w	14										
*12:03:23			3									12	13	14										
*12:03:26								8				12	13	14										
*12:03:34														14										
*12:04:01											11			14	15									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control

*12:04:02, 12:60										11	12	13	14	15												
*12:05										11	12		14	15												
*12:08											12	13			17											
*12:09										11				15				20								
*12:14:01												13					18	19								
*12:14:02											12	13					18	19								
*12:15								8			12	13	14		17											
*12:16											12	13												23		
*12:18:02, B*67:02 ³												13														
*12:21							7			11	12			15												
*12:24											12	13						20								
*12:28, 12:58											12	13	14				18									
*12:31											12	13	14							21						
*12:33							7			11			14	15												
*12:41										11	12	13		15												
*12:44											12	13								21						
*12:49, 12:123											12															
*12:54										11	12	13	14													
*12:55											12	13	14		16											
*12:72											12	13		15												
*12:81, 12:113											12	13	14			17										
*12:119	2										12	13	14													
*12:127									10		12	13														
*12:132											12	13					18									
*12:135											12	13	14	15			18									
*12:144												13	14							21						
*12:146										11	12	13		15		18										
*12:147											12	13	14											23		
*14:02:01-14:02:07, 14:02:09-14:03, 14:05, 14:07N, 14:11, 14:13-14:14, 14:18-14:24:02, 14:26-14:27, 14:29-14:44, 14:46-14:48, 14:50-14:52, 14:55-14:57, 14:59-14:62, 14:65-14:69															16				20							
*14:02:08											12				16				20							
*14:04, 14:64										11					16				20							
*14:06, 14:08-14:09, 14:28:01-14:28:02, 14:53, 14:63, A*30:62, B*35:252 ³															16											
*14:10															16		18		20							
*14:12, 14:49										11									20							
*14:15													14		16		18									
*14:16													14		16											
*14:17															16		18	19	20							
*14:25	2														16				20				22			
*14:45															15	16			20							
*14:54															16					21						
*15:02:01:01-15:02:09, 15:02:11-15:02:16, 15:02:18-15:02:20, 15:08, 15:10:02-15:10:03, 15:13, 15:15, 15:17-15:18, 15:26, 15:28, 15:32Q-15:35, 15:38-15:39, 15:41, 15:44-15:45, 15:47-15:53, 15:56-15:57, 15:60, 15:64, 15:67-15:68, 15:73-15:74, 15:78-15:84Q, 15:86-15:89, 15:91-15:101, 15:103, 15:106			3					8			11				15		17				21					
*15:02:10, 15:02:17			3					8			11				15		17				21	22				
*15:03, 15:16			3					8			11		13		15		17				21					
*15:04-15:05:08, 15:05:10, 15:06:01-15:06:03, 15:09, 15:19, 15:24, 15:27, 15:29-15:30, 15:40, 15:46, 15:54, 15:59, 15:66, 15:69-15:70, 15:76-15:77, 15:90, 15:105Q			3								11				15		17				21					
*15:05:09, 15:61											11				15		17				21					
*15:05:11			3												15		17				21					
*15:07			3					8								17					21					
*15:10:01								8			11				15		17				21					
*15:11			3			6		8			11				15		17				21					
*15:12					5			8			11				15		17	18			21					
*15:21, 15:71			3					8							15		17				21					
*15:22, 15:55			3								11				15		17									
*15:23											11					17					21					
*15:25			3									13					18	19			21					
*15:31, 15:62			3					8			11				15		17	18			21					
*15:37, 15:58, 15:72			3					8			11				15		17									
*15:42		2	3					8			11				15		17				21					
*15:43			3					8								17					21	22				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		

Negative Control

*15:63								8			11					17			21					
*15:65			3								11			15		17	18	19						
*15:75			3					8			11		13	15		17	18			21	23			
*15:85			3					8			11					17				21				
*15:102			3								11			15		17			20					
*15:104	1		3								11			15		17				21				
*16:01:01-16:01:15, 16:01:18-16:01:19, 16:06-16:08, 16:10-16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39:02, 16:41, 16:44, 16:49, 16:51-16:52, 16:54, 16:56, 16:58-16:59, 16:62, 16:65, 16:67, 16:71-16:73, 16:75-16:76, 16:79													13				18				23			
*16:01:16								8					13				18				23			
*16:01:17, 16:15:02											12	13					18				23			
*16:02:01-16:02:12, 16:09, 16:12, 16:19, 16:25, 16:46-16:47, 16:60, 16:63, 16:69, 16:74, 16:77N											11		13		15		18				23			
*16:02:13											11	12	13		15		18				23			
*16:04:01, 16:04:03, 16:33, 16:42, 16:55, 16:61, 16:66, 16:78														14			18		20					
*16:18, 16:37													13		15		18				23			
*16:20, 16:64													13				17	18			23			
*16:29		2												14			18		20					
*16:34			3			6							13				18				23			
*16:35														14			18							
*16:40							7									16								
*16:43													13								23			
*16:45													13				18			21	23			
*16:48											11			14	15		18							
*16:50		2											13				18				23			
*16:53							7									16			20					
*16:57											11		13				18				23			
*16:68																	18		20					
*16:70			3								11		13		15		17	18			23			
*17:01:01:01-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:16, 17:18-17:19, 17:21, 17:24-17:28		2									11				15			19		21				
*17:01:04		2									11		13		15			19		21				
*17:05		2					7				11				15			19		21				
*17:07, 17:17, 17:23											11				15			19		21				
*17:20		2									11							19		21	23			
*17:22		2																19		21				
*18:01-18:02, 18:04-18:08											11								20	21	23			
*18:03		2									11								20	21	22	23		
*18:09											11				15				20	21				
A*24:64, A*24:96, A*24:106, A*24:174, B*08:16, B*08:123																					23			
B*07:239, B*14:03																			20					
B*15:277, B*35:222															15									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control