

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
HLA-C allele^{1,2,3}																																
*01:02:01-01:02:31, 01:02:33-01:02:35, 01:46-01:48, 01:51, 01:53, 01:57, 01:60-01:65, 01:68, 01:71-01:72, 01:75-01:76, 01:87-01:88, 01:92, 01:94-01:96, 01:100, 01:103-01:106	1												13																			
*01:02:32, 01:90-01:91, 01:97, <i>B*54:18</i> ⁴													13																			
*01:03, 01:78		2											13																			
*01:04			3					8					13									23										
*01:05, 01:77				4									13																			
*01:06, 01:38 ⁵	1				5								13																			
*01:07:01-01:07:02	1					6							13																			
*01:08	1						7						13																			
*01:09	1							8					13																			
*01:10, 01:52 ⁵	1									9			13																			
*01:11	1										10		13																			
*01:12:01-01:12:02, 01:79:01-01:79:02, 01:101-01:102												11	13																			
*01:13	1											12	13																			
*01:14	1												13																			31
*01:15	1	2											13																			
*01:16	1			4									13																			
*01:17, 01:41 ⁵	1												13	14																		
*01:18, 01:42 ⁵	1												13		15																	
*01:19, 01:43, 01:58 ⁵	1												13			16																
*01:20	1				5		7						13																			
*01:21													13	14																		
*01:22	1							8		10			13																			
*01:23	1												13	14		16																
*01:24		2											13				17															
*01:25, 01:59 ⁵	1												13				17															
*01:26	1												13					18														
*01:27, 01:45 ⁵	1												13						19													
*01:28, 01:56N ⁵	1												13							20												
*01:29, 01:33 ⁵	1												13								21											
*01:30	1									10			13						19													
*01:31, 01:44, 01:107 ⁵	1												13										23									
*01:32, 01:40 ⁴	1												13									22										
*01:34											11		13				18															
*01:35	1							8		10			13										23									

negative Control

