

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43		
*04:06:01-04:06:02										10			13			16													29													42			
*04:08, 14:20																														30													42		
*04:16, 04:103, 05:12, 05:18:01, 05:103:02, 06:04:02, 06:06, 06:40, 06:79N, 12:09, 16:02:01-16:02:02, 16:02:04-16:02:07, 16:02:09-16:02:14, 16:12, 16:19, 16:25, 16:37, 16:46-16:47, 16:60, 16:63, 16:69-16:70, 16:74, 16:77N, 16:84, 16:88-16:89N, 16:99, 16:101-16:104, 16:107-16:108, 16:115, 17:17, 18:01-18:02, 18:04-18:07N, 18:09-18:10																		18																									42		
*04:42:01-04:42:02												12						18																										42	
*04:107										10																			29														42		
*04:108, 16:20																								24																			42		
*04:110																																											?		
*04:140, 04:166																																											42		
*04:147										10			13																29	30													42		
*04:160										10			13			16													29															42	
*04:178, 15:02:21, 15:07, 15:36, 15:43, 15:85, 15:144																16								24																			42		
*04:198, 14:41				4																																							42		
*04:220												12						18													31												42		
*05:08, 05:52, 05:89, 06:02:27, 06:16N, 06:41, 06:62, 06:92, 06:112, 06:206	1																	18																											
*05:11, 05:17, 05:68, 05:70, 05:79, 05:115, 05:134, 15:65, 17:01:11, 17:06, 17:34																16		18																											
*05:18:04											11		14				17	18																											
*05:26, 07:272, 07:326										10																																			
*05:27, 05:151																16		18												30															
*05:39																														30															
*05:43												12						18																											
*05:98																		18													31														
*05:105					5													18																											
*05:106:01	1	2												14			17	18																											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43		

Negative Control

*05:106:02, 06:02:01:01-06:02:01:10, 06:02:03-06:02:08, 06:02:10-06:02:20, 06:02:22-06:02:26, 06:02:28-06:03:02, 06:07, 06:09-06:10, 06:12-06:15, 06:17, 06:19- 06:32, 06:34:01-06:39, 06:42:01-06:61, 06:63- 06:78, 06:80-06:81, 06:83-06:88, 06:90- 06:91, 06:93-06:95, 06:97-06:103, 06:105- 06:106:02, 06:108-06:111, 06:113-06:117, 06:121-06:124, 06:126-06:127:01:01, 06:127:02-06:146, 06:148-06:152N, 06:154- 06:167, 06:169-06:178, 06:180-06:187, 06:189-06:196, 06:198-06:203, 06:205, 06:207-06:209, 06:212, 12:04:01-12:05, 12:21, 12:33, 12:41, 12:60, 12:72, 12:135, 12:146, 12:154, 16:91, 18:03	1													18																											42					
*05:132														18		24																														
*06:04:01, 06:118, 06:153, 08:10, 15:02:18, 15:02:22, 15:04:01-15:04:03, 15:30, 15:40, 15:64, 15:66, 15:76, 15:105Q, 15:110, 15:129, 16:48, 17:01:01:02-17:01:08, 17:01:10, 17:01:12-17:05, 17:08-17:16, 17:18- 17:20, 17:23-17:33, 17:35-17:38											16			18																												42				
*06:05								w						18																											31			42		
*06:08	1				5									18																												42				
*06:89	1													18			24																									42				
*06:96	1													18																										30			42			
*06:107, 06:179	1													18				25																								42				
*06:168	1					6								18																												42				
*06:197, 15:77, 17:07													16	18																									30			42				
*07:02:09									12																														31							
*07:02:75				3					11					14				17																												
*07:101, 07:148, 07:161, 07:583, A*01:203, A*11:166, A*30:56, A*31:85, A*80:01:01:01- 80:03, B*18:96, B*44:157, B*55:74																																											33			
*08:14, 08:80														16																												31			42	
*08:24														16																									29		31			42		
*08:29, 12:03:16, 12:03:27, 12:03:30, 12:32, 12:144, 12:185	1																																													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43			
*08:31	1				5																																									
*08:37, B*27:84												12																																		

Negative Control

