

INTERPRETATION TABLE				
HLA-C*04 Add-on SSP typing				
	Well			
	1	2	3	4
Length of spec.	260	220	180	105
PCR product	360			
Length of int.				
pos. control ¹	800	1070	1070	1070
5'-primer(s) ²	28	1018	1018	900
	5' -TCA 3'	5' -gTg 3'	5' -gTg 3'	5' -CCg 3'
3'-primer(s) ³	112	1092	1052	963
	5' -CCA 3'	5' -TTA 3'	5' -Tgg 3'	5' -gCT 3'
	118			
	5' -gCT 3'			
	218			
	5' -gCT 3'			
Well No.	1	2	3	4
HLA-C allele ⁴				
*04:01:01:01-04:01:51, 04:03-04:08, 04:10-04:20, 04:23-04:29, 04:31-04:81, 04:83-04:145, C*05:78, C*16:34	1			
*04:09N	1	2		
*04:30	1		3	
*04:82	1			4
HLA-C allele				
Well No.	1	2	3	4