

Length of spec.	250	220	145	210	250	95	145	270	220	180	110	110	195	180	125	225	155	180	130	320	125	120	165	85	130	85	170																	
PCR product(s)									220	180	110		220	195	180		165	125	270	225	190	155	180	130	320	220	125	155	120	545	250	165	240	145	85	130	215	145	85	170				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24																				
HLA-C allele ^{4,5}																																												
*04:01:01:01-04:01:09, 04:01:11-04:01:22, 04:01:24- 04:01:27, 04:01:29-04:01:37, 04:82, 04:85-04:86, 04:90- 04:92, 04:97-04:99, 04:101- 04:102, 04:104, 04:109, 04:112	1	2	3									11																																
*04:01:10, 04:01:23	1	2	3																																									
*04:01:28	1		3									11																																
*04:03, 02:49 ⁶	1			4								11																																
*04:04:01-04:04:02		2	3		5							11																																
*04:05	1	2	3			6						11																																
*04:06				4	5							11																																
*04:07, 04:89	1	2					7					11																																
*04:08	1	2	3					8				11																																
*04:09N	1	2	3						9			11																																
*04:10	1		3							10	11																																	
*04:11	1		3							10		12																																
*04:12	1	w	3								11		13																															
*04:13		2	3		5						11			14																														
*04:14		2	3								11				15																													
*04:15:01, 04:15:03	1	2	3								11					16																												
*04:15:02, 04:44, 04:47 ⁸	1	2	3								11													20																				
*04:16	1										11			14																														
*04:17, 04:100	1	2	3								11										17		20																					
*04:18	1	2	3								11			14																														
*04:19	w	2	3								11																																	
*04:20	1	2	3								11																																	
*04:23, 04:38, 04:39, 04:108 ⁹	1	2	3								11																																	
*04:24	1	2	3								11																																	
*04:25, 04:40, 04:41 ¹⁰	1	2	3								11																																	
*04:26	1	2	3								11																																	
*04:27	1	2					7																																					
*04:28	1	2	3								11				15																													
*04:29	1	2	3								11	12																																
*04:30	1	2	3						9		11																																	
*04:31	1	2	3							10	11																																	
*04:32	1	2	3				7				11																																	
*04:33	1	2	3									12																																
*04:34		2	3		5			8			11																																	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24																				
Length of spec.	250	220	145	210	250	95	145	270	220	180	110	110	195	180	125	225	155	180	130	320	125	120	165	85	130	85	170																	
PCR product(s)									220	180	110		220	195	180		165	125	270	225	190	155	180	130	320	220	125	155	120	545	250	165	240	145	85	130	215	145	85	170				

200	215	285	215	185	75	180	275	180	200	175	115	135	135	330	270	280	85	260	135	160	300	240	Length of spec.	PCR product(s)
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	Well No.	HLA-C allele ^{4,5}
																								*04:01:01:01-04:01:09, 04:01:11-04:01:22, 04:01:24- 04:01:27, 04:01:29-04:01:37, 04:82, 04:85-04:86, 04:90- 04:92, 04:97-04:99, 04:101- 04:102, 04:104, 04:109, 04:112
																								*04:01:10, 04:01:23
																								*04:01:28
																								*04:03, 02:49 ⁶
																								*04:04:01-04:04:02
																								*04:05
																								*04:06
																								*04:07, 04:89
																								*04:08
																								*04:09N
																								*04:10
																								*04:11
									34															*04:12
																								*04:13
																								*04:14
																								*04:15:01, 04:15:03
																				45				*04:15:02, 04:44, 04:47 ⁸
																								*04:16
																								*04:17, 04:100
																								*04:18
																								*04:19
																								*04:20
																								*04:23, 04:38, 04:39, 04:108 ⁹
																								*04:24
																								*04:25, 04:40, 04:41 ¹⁰
																								*04:26
																								*04:27
																								*04:28
																								*04:29
																								*04:30
																								*04:31
																								*04:32
																								*04:33
																								*04:34
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	Well No.	HLA-C allele ^{4,5}
200	215	285	215	185	75	180	275	180	200	175	115	135	135	330	270	280	85	260	135	160	300	240	Length of spec.	PCR product(s)

