

## INTERPRETATION TABLE

### HLA-C\*04 Add-on SSP typing

	Well			
	1	2	3	4
Length of spec.	255	220	180	105
PCR product	360			
Length of int.				
pos. control <sup>1</sup>	<b>800</b>	<b>1070</b>	<b>1070</b>	<b>1070</b>
5'-primer(s) <sup>2</sup>	<b>28</b>	<b>1018</b>	<b>1018</b>	<b>900</b>
	5' -TCA 3'	5' -gTg 3'	5' -gTg 3'	5' -CCg 3'
3'-primer(s) <sup>3</sup>	<b>112</b>	<b>1092</b>	<b>1052</b>	<b>963</b>
	5' -CCA 3'	5' -TTA 3'	5' -Tgg 3'	5' -gCT 3'
	<b>118</b>			
	5' -gCT 3'			
	<b>218</b>			
	5' -gCT 3'			
Well No.	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>
HLA-C allele <sup>4</sup>				
*04:01:01:01-04:01:41, 04:03-04:08, 04:10-04:20, 04:23-04:29, 04:31-04:81, 04:83-04:127, 05:78, 16:34	1			
*04:09N	1	2		
*04:30	1		3	
*04:82	1			4
HLA-C allele				
Well No.	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>