

INTERPRETATION TABLE

HLA-C*04:09N SSP typing

	Well	
	1	2
Length of spec.	220	255
PCR product		360
Length of int.	800	1070
pos. control ¹		
5'-primer(s) ²	1018	28
	5' -gTg 3'	5' -TCA 3'
3'-primer(s) ³	1092	112
	5' -TTA 3'	5' -CCA 3'
		118
		5' -gCT 3'
		218
		5' -gCT 3'
Well No.	1	2
HLA-C allele		
*04:09N	1	2
*04:01:01:01-04:01:48, 04:03-04:08, 04:10-04:20, 04:23-04:144, 05:78, 16:34		2
HLA-C allele		
Well No.	1	2