

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48			
HLA-C allele ^{1,2,3}																																																			
*06:02:01:01-06:02:01:03, 06:02:03-06:02:43, 06:48, 06:51-06:52, 06:56, 06:61-06:64, 06:68, 06:75, 06:78, 06:80, 06:82, 06:84-06:85, 06:88-06:92, 06:94-06:95, 06:97-06:99, 06:102, 06:105-06:106:02, 06:108-06:110, 06:112-06:115, 06:117, 06:119, 06:121, 06:125, 06:129-06:131, 06:135, 06:137, 06:139-06:142, 06:145, 06:149	1	2																																																	
*06:03:01		2	3																																																
*06:03:02		2																																																	
*06:04	1			4																																															
*06:05	1				5	6														w																															
*06:06	1					6																																													
*06:07, 06:33⁵	1	2					7																																												
*06:08	1	2						8																																											
*06:09	1	2							9																24																										
*06:10	1	2								10																																									
*06:11, 06:147⁵	1	2									11																																								
*06:12	1	2										12																																							
*06:13	1	2											13																																						
*06:14, 06:143	1													14																																					
*06:15, 06:116N⁵	1	2													15																																				
*06:16N, 06:21⁵	1	2														16																																			
*06:17		2							9																																										
*06:18	1	2	3		5																																														
*06:19	1	2					7																																												
*06:20, 06:74Q⁵	1	2																																																	
*06:22	1	2								10	11																																								
*06:23	1	2			5			9																																											
*06:24, 06:37⁵	1	2															17																																		
*06:25, 06:36⁵	1	2																			21																														
*06:26	1	2																18																																	
*06:27, 06:29⁵	1	2																			20																														
*06:28	1	2																				19																													
*06:30	1	2																					22																												
*06:31	1	2																																																	
*06:32	1	2																			19			23																											

Negative Control

