

HLA-C07 2019 2H6.xlsx

HLA-C07 2019 2H6.xlsx

HLA-C07 2019 2H6.xlsx

HLA-C07 2019 2H6.xlsx

HLA-C07 2019 2H6.xlsx

HLA-C07 2019 2H6.xlsx

HLA-C07 2019 2H6.xlsx

Well No.
HLA-C allele ^{1,2,3}
*07:01:01:01-07:01:01:24, 07:01:04, 07:01:06-07:01:12, 07:01:14-07:01:24, 07:01:26, 07:01:28-07:01:35, 07:01:37- 07:01:51 , 07:01:53-07:01:61, 07:103, 07:111, 07:113, 07:115, 07:156, 07:162- 07:163, 07:182, 07:188, 07:190, 07:197, 07:200- 07:201, 07:204:01-07:204:02, 07:206, 07:215, 07:219, 07:223, 07:225, 07:228, 07:236-07:237, 07:246:01- 07:246:02, 07:249, 07:253, 07:255, 07:257:01-07:257:04, 07:266, 07:268-07:269, 07:277- 07:278, 07:280-07:282, 07:296- 07:300, 07:303, 07:313, 07:331-07:332, 07:362, 07:366, 07:369-07:370, 07:377, 07:380, 07:385, 07:408, 07:417-07:418, 07:422, 07:424, 07:427, 07:430-07:431, 07:435, 07:440- 07:444, 07:448, 07:454:01- 07:455, 07:460, 07:468- 07:471, 07:475, 07:479, 07:481, 07:488, 07:490, 07:492-07:493, 07:499, 07:502, 07:506, 07:508, 07:515, 07:521-07:522, 07:527- 07:528, 07:537, 07:540, 07:545-07:546, 07:549-07:550, 07:554-07:555, 07:557, 07:561, 07:570-07:571, 07:581, 07:588, 07:591, 07:602, 07:604, 07:610, 07:614, 07:617-07:619, 07:623- 07:624, 07:634
*07:01:02 , 07:01:27, 07:01:52, 07:547
*07:01:03
*07:01:05, 07:01:36 , 07:110, 07:230, 07:353, 07:463, 07:526, 07:587
Well No.
HLA-C allele ^{1,2,3}
*07:01:13
*07:01:25
*07:02:01:01-07:02:01:16,

07:02:01:18-07:02:86, 07:105, 07:114, 07:135, 07:144, 07:146, 07:157-07:160, 07:167, 07:169, 07:174- 07:175, 07:178, 07:185, 07:187, 07:193, 07:195, 07:208-07:209, 07:211, 07:213, 07:217, 07:221, 07:233, 07:239-07:240, 07:251, 07:261-07:262, 07:270- 07:271, 07:273-07:275, 07:283- 07:284, 07:286-07:287, 07:290- 07:291, 07:305-07:306, 07:309, 07:312, 07:318- 07:321, 07:325, 07:327, 07:333-07:334, 07:341:01- 07:341:02, 07:344-07:346, 07:363, 07:372, 07:376, 07:379, 07:382-07:384, 07:388- 07:389, 07:391-07:392, 07:398- 07:401, 07:405, 07:413- 07:416, 07:432-07:434, 07:436, 07:438, 07:449, 07:464-07:465, 07:472-07:473, 07:477-07:478, 07:482, 07:489, 07:495-07:496, 07:504- 07:505, 07:509-07:511, 07:516, 07:518, 07:520, 07:525, 07:529-07:530, 07:532, 07:539, 07:541- 07:543, 07:548, 07:553, 07:560, 07:564-07:567, 07:572, 07:574-07:577, 07:580, 07:584, 07:590, 07:592, 07:594-07:597, 07:601, 07:605-07:606, 07:608, 07:612, 07:616, 07:625, 07:628-07:631
*07:02:01:17N
Well No.
HLA-C allele^{1,2,3}
*07:03
*07:04:01:01-07:04:13, 07:139, 07:142, 07:323, 07:336, 07:354, 07:357, 07:364-07:365, 07:378, 07:394- 07:395, 07:420, 07:428, 07:459, 07:466-07:467, 07:501, 07:534-07:535, 07:569, 07:585-07:586, 07:622, 07:626
07:05, 07:538

*07:06:01:01-07:06:01:02, 07:613, 07:620
*07:06:02
*07:07
*07:08
*07:09
*07:10
*07:11
*07:12
*07:13, 07:289
*07:14
*07:15
*07:16
*07:17:01-07:17:04
*07:18, 07:517, 07:621
*07:19
*07:20
*07:21
*07:22
*07:23
*07:24
*07:25
*07:26:01-07:26:03
*07:27:01-07:27:02, 07:556, 12:181 ⁴
*07:28
*07:29:01-07:29:02, 07:347N ⁵
*07:30
*07:31:01-07:31:02
*07:32N
*07:33N
*07:35
*07:36
*07:37
*07:38:01-07:38:02
*07:39, 07:335
*07:40
*07:41, 07:447 ⁵
*07:42
Well No.
HLA-C allele^{1,2,3}
*07:43:01-07:43:02
*07:44
*07:45
*07:46
*07:47
*07:48
*07:49
*07:50
*07:51, 07:314:01-07:314:03⁵
*07:52, 07:337 ⁵
*07:53
*07:54

*07:55N, 07:514
*07:56:01-07:56:02
*07:57
*07:58
*07:59
*07:60
*07:61N, 07:194, 07:368:01-
07:368:02 ⁵
*07:62
*07:63
*07:64
*07:65
*07:66
*07:67
*07:68
*07:69
*07:70
*07:71
*07:72
*07:73:01
*07:73:02
*07:74, 07:308 ⁵
*07:75
*07:76:01
*07:76:02
*07:77
*07:78:01-07:78:02, 07:235Q,
07:559 ⁵
*07:79
*07:80
*07:81
*07:82
*07:83
*07:84
*07:85
*07:86
*07:87
*07:88, 07:474⁵
*07:89
*07:90
*07:91
*07:92
Well No.
HLA-C allele^{1,2,3}
*07:93
*07:94
*07:95
*07:96:01-07:96:02
*07:97
*07:98N
*07:99, 07:183 ⁵
*07:100, 07:161⁵
*07:101
*07:102
*07:104N, 07:119 ⁵

*07:106
*07:107
*07:108,01-07:108:02
*07:109,01-07:109:02
*07:112
*07:116
*07:117, 07:126 ⁵
*07:118, 07:203 ⁵
*07:120, 07:462 ⁵
*07:121Q
*07:122
*07:123
*07:124
*07:125, 07:402 ⁵
*07:127
*07:128, 07:279
*07:129, 07:531
*07:130
*07:131,01-07:131:02
*07:132, 07:179 ⁵
*07:133
*07:134
*07:136
*07:137,01-07:137:02
*07:138, 07:484N ⁵
*07:140
*07:141,01-07:141:02
*07:143
*07:145
*07:147
*07:148
*07:149, 07:381 ⁵
*07:150Q
*07:151
*07:152N
*07:153
*07:154
*07:155
*07:164N, 07:212 ⁵
*07:165, 07:180 ⁵
*07:166, 07:419 ⁵
Well No.
HLA-C allele ^{1,2,3}
*07:168
*07:170, 07:412
*07:171
*07:172,01-07:172:02
*07:173
*07:176
*07:177
*07:181
*07:184
*07:186, 07:451N ⁵
*07:189
*07:191N

*07:192
*07:196
*07:198N
*07:199:01-07:199:02
*07:202
*07:205
*07:207
*07:210
*07:214
*07:216
*07:218
*07:220
*07:222
*07:224
*07:226, 07:573 ⁵
*07:227N, 07:254 ⁵
*07:229
*07:231
*07:232
*07:234
*07:238
*07:241
*07:242
*07:243
*07:244
*07:245
*07:247
*07:248
*07:250
*07:252
*07:256
*07:258, 07:609
*07:259
*07:260:01-07:260:02
*07:263
*07:264N
*07:265
*07:267
*07:272
*07:276
Well No.
HLA-C allele^{1,2,3}
*07:285
*07:288
*07:292
*07:293
*07:294
*07:301
*07:302
*07:304
*07:307
*07:310
*07:311
*07:315
*07:316

*07:317
*07:322
*07:324
*07:326
*07:328
*07:329N
*07:330:01
*07:330:02
*07:338
*07:339
*07:340
*07:342
*07:343
*07:348
*07:349
*07:350N, 07:497 ⁵
*07:351
*07:352
*07:355
*07:356
*07:358, 07:523, 07:552
*07:359
*07:360
*07:361
*07:367
*07:371
*07:373
*07:374
*07:375
*07:386
*07:387
*07:390
*07:393N
*07:396
*07:397
*07:403
*07:404
*07:406
*07:407
Well No.
HLA-C allele ^{1,2,3}
*07:409
*07:410
*07:411
*07:421
*07:423
*07:425
*07:426
*07:429
*07:437N
*07:439
*07:445
*07:446
*07:450
*07:452N

*07:453
*07:456
*07:457
*07:458
*07:461
*07:476N
*07:480
*07:483N
*07:485
*07:486
*07:487
*07:491:01N-07:491:02N
*07:494Q
*07:498
*07:500
*07:503
*07:507N
*07:512
*07:513Q
*07:519
*07:524
*07:533
*07:536
*07:544
*07:551N
*07:558
*07:562
*07:563
*07:568
*07:578
*07:579
*07:582Q
*07:583
*07:589
*07:593N
*07:598
*07:599
*07:600:01N-07:600:02N
*07:603N
*07:607
*07:611
*07:615
*07:627
*07:632Q
*07:633N
Well No.
HLA-C allele^{1,2,3}
*01:17, 01:21, 01:23, 01:128, 01:152, 12:215
*01:60, 17:01:01:02-17:04, 17:06-17:19, 17:21-17:24, 17:27N-17:31, 17:33-17:34, 17:36-17:37, 17:39-17:40
*01:111N, 14:47N, <i>B</i> *15:209N, <i>B</i> *46:55N

*02:02:13, 02:02:29, 12:02:11, <i>B*15:434</i>
*02:02:28
*02:14:01-02:14:02, 02:107, 03:39, 03:67, 04:42:01- 04:42:02, 04:220, 06:02:01:01- 06:02:01:14, 06:02:03- 06:02:05, 06:02:07-06:02:11, 06:02:14-06:02:22, 06:02:24- 06:02:31, 06:02:33-06:02:54, 06:04:01-06:06, 06:08-06:10, 06:12-06:16N, 06:19-06:21, 06:23-06:35, 06:37-06:43:02, 06:45-06:51, 06:54-06:66, 06:68-06:69, 06:71-06:75, 06:78-06:100, 06:102:01- 06:117, 06:119, 06:121, 06:125-06:126, 06:128N, 06:130-06:131, 06:133-06:135, 06:137-06:143, 06:145-06:146, 06:148-06:152N, 06:154- 06:174, 06:176-06:188, 06:190- 06:193, 06:195-06:216, 06:218- 06:219, 15:23:01-15:23:02, 15:63, 15:138, <i>B*08:17</i> , <i>B*08:39, B*08:47</i>
*02:49, 02:75, 02:115, 03:15, 03:32, 03:45, 03:60, 03:136, 03:297, 04:03:01-04:03:04, 04:06:01-04:06:02, 04:16, 04:80, 04:147, 04:171, 04:190, 04:256, 04:286, 06:132:01- 06:132:02, 15:25, <i>A*01:110</i> , <i>A*11:135</i>
*02:57, 04:40, 04:242, <i>A*02:84, B*38:23, B*39:73</i>
*02:77
*02:99, 04:127, <i>B*07:87</i> , <i>B*15:272N, B*27:55</i> , <i>B*27:149, B*39:111</i> , <i>B*40:286N</i>
*03:03:29, 03:04:19
*03:14, 03:121N, 03:261, 03:361, 04:203N, 05:143, 12:14:01-12:14:02, 12:176, 14:17, <i>B*46:06, B*46:43</i>
Well No.
HLA-C allele^{1,2,3}
*03:69, 03:274, 04:27, 04:52, 12:58
*03:113:01

*03:137, 03:182, B*07:13, B*07:247, B*07:270, B*39:121, B*54:02
*03:162
*03:178, 03:300
*03:211:02
*03:220
*03:251, 03:314, 04:120, 05:01:01:01-05:01:16, 05:01:18-05:01:26, 05:01:28- 05:01:39, 05:04:01-05:08, 05:10-05:16, 05:18:01- 05:22:02, 05:24-05:41, 05:44:01-05:45, 05:47-05:48N, 05:50-05:51Q, 05:53-05:61, 05:63-05:83, 05:85-05:94, 05:96-05:109, 05:112-05:124, 05:126-05:132, 05:135, 05:137- 05:142, 05:144-05:145, 05:147- 05:148, 05:150, 05:152- 05:156, 05:158-05:162, 05:164, 08:01:01-01-08:04:03, 08:06, 08:08:01-08:14, 08:16:01-08:20, 08:22-08:37, 08:39-08:46, 08:48-08:50, 08:52N-08:63, 08:65-08:71, 08:73-08:75, 08:77-08:83, 08:85-08:103, 08:105-08:123, 08:125-08:136, 08:138- 08:141Q, 08:143-08:156, 08:158-08:159, 08:161N- 08:167, 15:130, A*29:10:01- 29:10:02, B*14:32, B*15:337, B*18:83, B*58:76
*03:342, B*56:14
*04:01:43
*04:34, 04:122, 04:212, 04:265, A*26:112, B*18:115, B*35:178, B*35:282, B*35:316
*04:58
*04:66
*04:112, 16:64
*04:129, 12:127, 12:203
*04:160
Well No.
HLA-C allele^{1,2,3}
*04:198
*04:222, B*07:189
*04:258
*04:273
*05:01:17
*05:03

*05:09, 02, 08:15:02
*05:17
*05:23, 08:07, 08:104
*05:42, 05:46, 08:05, 08:21, 08:137
*05:43, 06:129
*05:62
*05:84, 05:157
*05:111, 08:38
*05:134, 08:47
*05:146, 05:163, 08:142, <i>B*44:148</i>
*05:149
*05:151
*06:02:13, 06:44
*06:02:23
*06:07
*06:17
*06:18
*06:22
*06:36
*06:53:01-06:53:02
*06:67
*06:70:01-06:70:02
*06:101, 06:127:01:01- 06:127:02, 06:136, 06:144
*06:118, 06:175N
*06:120
*06:123
*06:153
*06:189
*06:194
*08:72:02
*12:03:10
*12:42Q, <i>B*46:51Q</i>
*12:115
*12:131
*12:211
*15:62
*15:65
*16:02:11
*16:10
*16:16Q, <i>B*15:218Q,</i> <i>B*35:333Q</i>
*16:31, <i>B*15:193, B*35:283</i>
*16:76, <i>B*40:304</i>
Well No.
HLA-C allele^{1,2,3}
*16:111
*17:05
*17:20
*17:26
*17:38

*18:01-18:03, 18:06-18:08,
18:10
*18:04
*18:05
*18:09
A'02:01:126
A'03:01:46, A'03:154:02,
A'11:32:02, A'24:02:62,
A'24:347:01
A'11:179
A'11:232, A'36:06
A'25:27:02
A'25:32
A'26:94, A'26:101
A'30:01:03
A'33:03:09
A'68:01:12
B'07:02:43, B'07:33:03
B'07:06:03
B'07:77, B'07:193
B'08:01:25
B'08:109
B'15:17:04
B'27:36, B'48:21, B'48:26
B'35:01:38, B'35:15:02,
B'58:01:17
B'35:03:15
B'35:08:04
B'35:364, B'51:142
B'37:04:02
B'38:43, B'49:29, B'51:181,
B'52:74
B'39:09:03
B'39:136
B'40:60, B'40:245
B'40:100
B'44:247
B'46:07N
B'51:01:63
B'56:08
B'57:01:18
B'73:01-73:02
Well No.
HLA-C allele^{1,2,3}