

## INTERPRETATION TABLE

### HLA-C\*07 Add-on SSP typing

	Well				
	1	2	3	4	5
Length of spec.	245	505	225	80	80
PCR product	425				
Length of int.	<b>800</b>	1070	1070	1070	1070
pos. control <sup>1</sup>					
5'-primer(s) <sup>2</sup>	47	992	5 <sup>th</sup> I	1049	1049
	5' -Agg 3'	5' -TAA 3'	5' -gTC 3'	5' -g 3'	5' -C 3'
	<b>648</b>				
	5' -CAC 3'				
3'-primer(s) <sup>3</sup>	302	1016	1043	1087	1087
	5' -ggC 3'	5' -CAC 3'	5' -CAA 3'	5' -AgC 3'	5' -AgT 3'
	<b>853</b>				
	5' -CAT 3'				
Well No.	1	2	3	4	5
HLA-C allele <sup>4</sup>					
*07:01:01:01-07:05, 07:07-07:10, 07:12-07:17:02, 07:20-07:33N, 07:35-07:49, 07:51-07:294, 07:296-07:307	1				
*07:06	1	2	3		
*07:11	1			4	
*07:18-07:19	1		3		
*07:50	1				5
HLA-C allele					
Well No.	1	2	3	4	5