

HLA-C07 Add-on 2018 2H3

Well No.	1	2	3	4	5	6	7	8
HLA-C allele ^{1,2}								
*07:01:01:01-07:05, 07:07-07:10, 07:12-07:17:04, 07:20-07:33N, 07:35-07:49, 07:51-07:294, 07:296- 07:330:01, 07:331-07:342, 07:344-07:348, 07:350N- 07:418, 07:420-07:516, 07:518-07:606, 07:608-07:612, 07:614-07:619. 07:622-07:634. 12:181	1							Negative Control
*07:06:01:01-07:06:02, 07:330:02, 07:613, 07:620	1	2	3					
*07:11	1			4				
*07:18-07:19, 07:517, 07:607, 07:621	1		3					
*07:50	1				5			
*07:343	1	2				6		
*07:349	1						7	
*07:419	1			4			7	
*16:111, 17:38					5			
Well No.	1	2	3	4	5	6	7	