

HLA-C08 2019 1K5

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE	AF	AG	AH	AI	AJ	AK	AL	AM	AN	AO	AP	AQ											
157	*04:106, 14:39, 14:96																					21	23																															
158	*04:120, 05:18:05			3										13								21						28																										
159	*04:129					5																21									31																							
160	*04:140, 04:166:01						6															21					26																											
161	*04:144, 06:02:52, 06:73, 06:214																					21								30																								
162	*04:220, 06:05						6															21		24		26																												
163	*04:223:01					5																					27		30		33																							
164	*04:229																					21																										40						
165	*04:277, 16:133																16					21																																
166	*04:291	1								9												21																																
167	*04:347, 06:161						6															21																																
168	*05:01:01:01-05:01:19, 05:01:21-05:01:45, 05:03, 05:06-05:07N, 05:13-05:15, 05:18:04, 05:19-05:22:02, 05:24, 05:26, 05:28, 05:30-05:38, 05:40-05:41, 05:44:02-05:45, 05:47-05:51Q, 05:53-05:61, 05:63, 05:65-05:67, 05:69, 05:71-05:75, 05:77-05:78:02, 05:80-05:84, 05:86-05:87, 05:90-05:91N, 05:94-05:97, 05:99N-05:102, 05:104, 05:108-05:111, 05:113N-05:114, 05:116, 05:118-05:132, 05:136-05:142, 05:144-05:145, 05:149-05:150, 05:152-05:162, 05:164-05:183, 05:185-05:187, 05:189-05:196, 05:198-05:200, 05:202Q-05:205, 05:208N-05:211, 05:213N-05:214						3								13																																							
169	*05:01:20			3										13																	31																							
170	*05:04:01-05:04:02, 05:135, 05:188, 05:207			3							11			13													25																											
171	*05:05:01-05:05:02, 05:147			3										13													25																											
172	*05:08, 05:89, 05:106:01			3															19																																			
173	*05:10			3										13																32																								
174	*05:11	1		3																																																		
175	*05:12													13								21				25																												
176	*05:16, 05:85			3							11			13																																								
177	*05:17	1			5																																																	
178	*05:18:01			3										13								21																																
179	*05:18:02-05:18:03			3										13									22																															
180	*05:23			3			7							13																																								
181	*05:25			3										13					17											30																								
182	*05:27, 05:184	1		3			7							13									22																															
183	*05:29:01-05:29:02			3										13	15																																							
184	*05:39			3			7							13									22																															
185	*05:42			3		5								13					17																																			
186	*05:43			3										13																																								
187	*05:44:01			3	4									13									20																															

Negative Control

